Simultaneous Estimation of Total Homocysteine and Methylmalonic acid in **Clinical Plasma/Serum Samples by Using Acquity UPLC-XEVO TQD**

Taposh G, Bhaskar K, Veeranjaneyulu P, Tirupateswar rao B, Rajesh PMN, Anil Kurup Waters India Applications Laboratory, Bangalore, 560058

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

INTRODUCTION


Clinical background:

Vitamin B12 (Cobalamin) deficiency is a common clinical problem in the elderly subjects, neonates and infants, it's deficiency may leads to irreversible neurological damage.


In healthy population the coenzyme form of vitamin B12 participates in two key metabolic pathways. One is the conversion of methylmalonylcoenzyme- A to succinyl-CoA & second one is formation of methionine from homocysteine (Hcy).

In case of vitamin B12 deficiency, methylmalonyl-CoA accumulates and

UPLC WITH XEVO TQD



RESULTS AND DISCUSSION



subsequently increases Methylmalonic acid (MMA) and homocysteine concentration in blood. Monitoring of MMA and tHcy's can be an indicative early biomarkers to characterize intracellular, functional vitamin B12 deficiency.

OBJECTIVE:

Development of simultaneous estimation of MMA & total Homocystiene using UPLC-XEVO TQD mass spectrometer.

CHALLENGES:

Both the analytes are highly polar and exist endogenously, having low molecular weight and low pka values.

Simple and fast LC-MS/MS, MRM based assay was developed for the simultaneous determination of MMA & tHcy in plasma and serum samples.

Sample preparation demonstrated is simple, rapid, and analogous for biological matrices.



HO	ЪЮ		γ 0		18 Homocystiene MMA_120416_ PA 2_018 19 Homocystiene MMA_120416_ PA 2_019 20 Homocystiene MMA_120416_ PA 2_020	EXT LQC EXT LQC EXT LQC	QC QC QC	0.70 1 0.70 1 0.70 1	81306 0.70 82991 0.70 87513 0.70	524864 518841 538208	0.345 0.353 0.348	750.00 711.15 94.3 750.00 750.97 100.3 750.00 727.43 97.0	Figure 5: Representative chromatogram Figure 6: Calibration curve of MMA of MMA-d3. <i>of MMA and MMA-d3. standards (50 ppb to 500 ppb).</i>
			NH_2		21 Homocystiene MMA_120416_ PA 2_021 22 Homocystiene MMA_120416_ PA 2_022	EXT LQC EXT LQC	QC QC	0.70 1 0.70 1	86725 0.70 84749 0.70	531736 520217	0.351 0.355	750.00 742.57 99. 750.00 764.39 101.	
			1112		23 Homocystiene MMA_120416_ PA 2_023 24 Homocystiene MMA_120416_ PA 2_024	EXT MQC EXT MQC	QC QC	0.70 3	827290 0.70 823760 0.70	570503 557036	0.574	2000.00 1963.69 98.3 2000.00 1999.88 100.0	CONCLUSION
		f Figure 2. Chamigal structure of			25 Homocystiene MMA_120416_ PA 2_025	EXT MQC	QC	0.70 3	326166 0.70	563121	0.579	2000.00 1994.01 99.	1) Demonstrated simultaneous quantification method for MMA and
Figure 1: Chemical structure of Figure 2: Chemical structur			ai structure o	J	26 Homocystiene MMA_120416_ PA 2_026 27 Homocystiene MMA_120416_ PA 2_027	EXT MQC EXT MQC	QC QC	0.70 3	316217 0.70	545132 548591		2000.001998.7599.92000.001970.9898.9	
MMA tHcy				28 Homocystiene MMA_120416_PA 2_028	EXT MQC	QC	0.70 3	318248 0.70	554655	0.574	2000.00 1970.30 98.		
					29 Homocystiene MMA_120416_ PA 2_029 30 Homocystiene MMA_120416_ PA 2_030	EXT HQC EXT HQC	QC QC	0.70 5	521036 0.70	542369		4000.00 4087.28 102.1 4000.00 3969.17 99.1	
METHOD SUMMARY					31 Homocystiene MMA_120416_ PA 2_030 31 Homocystiene MMA_120416_ PA 2_031	EXT HQC	QC	0.70 5	510723 0.70	540120 540101		4000.00 3909.17 99	2) The method demonstrated showed less matrix effect and good
Sample preparation: To the plasma/serum sample, two Internal				hal	32 Homocystiene MMA_120416_ PA 2_032	EXT HQC	QC	0.70 5	502725 0.70	529924		4000.00 4021.46 100.4	sensitivity up to 50 ppb for MMA & 500 ppb for tHcy.
					33 Homocystiene MMA_120416_ PA 2_033	EXT HQC	QC	0.70 5	02793 0.70	540574	0.930	4000.00 3919.59 98.	
standard was added with reduction reagent and vortexed for 5					Compound name: MMA								\neg 3) Opportunity for this method to be validated and used or routine
minutes. Extraction solution was added and centrifuged, the					# Name	D	Туре	RT	Area IS RT	IS Area	Response	Std. Conc Conc. %R	analysis.
supernatant was taken for LC -MS/MS analysis.					1 Homocystiene MMA_120416_ PA 2_001	EXT BLK		0.74	102 0.75	1	87.647	55882	
					2 Homocystiene MMA_120416_ PA 2_002 3 Homocystiene MMA_120416_ PA 2_003	EXT BLK IS EXT STD1	Standard	0.74 0.74	100 0.74 177 0.74	945 876	0.106 0.202	50.00 49.03 9	REFERENCES
					4 Homocystiene MMA_120416_PA 2_004	EXT STD2	Standard	0.74	222 0.74	894	0.248	75.00 78.47 10	.6
MRM CONDITIONS					5 Homocystiene MMA_120416_ PA 2_005	EXT STD3	Standard	0.74	241 0.74	873	0.277	100.00 96.64 9	1. Anne Lise Bjorke Monsen, et al. (2003). Homocysteine and
					6 Homocystiene MMA_120416_ PA 2_006 7 Homocystiene MMA 120416 PA 2 007	EXT STD4 EXT STD5	Standard	0.74	412 0.74 521 0.74	907 902	0.454	200.00 209.87 104 300.00 288.87 90	
Compound	Parent m/z	Daughter m/z	Cone (V) C	E	8 Homocystiene MMA_120416_ PA 2_008	EXT STD6	Standard	0.74	654 0.74	874	0.748	400.00 397.43 9	adolescence. American Society for Clinical Nutrition.
					9 Homocystiene MMA_120416_ PA 2_009 10 Homocystiene MMA_120416_ PA 2_010	EXT STD7 EXT blk	Standard	0.74 0.74	820 0.74 105 0.73	901 1	0.909 124 430	500.00 500.52 10 79368	
Homocysteine	136.08	90.03	10 1	.0	11 Homocystiene MMA_120416_ PA 2_011	EXTLOQ	QC	0.74	169 0.74	857	0.197	50.00 45.64 9	2. Xiaowei., et al. (2012).Simple, Fast, and Simultaneous Detection of
Homocysteine- d4	140.08	94.03	10 1	0	12 Homocystiene MMA_120416_ PA 2_012	EXTLOQ	QC	0.74	172 0.74	864	0.199	50.00 46.74 93	
				2	13 Homocystiene MMA_120416_ PA 2_013 14 Homocystiene MMA_120416_ PA 2_014	EXT LOQ EXT LOQ	QC	0.75 0.74	160 0.74 163 0.74	877 855	0.182	50.00 36.41 72 50.00 41.70 83	A Menigician actu using Liquia cinomatography antu mass
Methylmalonic acid	117.00	73.00	15 1	. 2	15 Homocystiene MMA_120416_ PA 2_015	EXT LOQ	QC	0.74	163 0.74	882	0.185	50.00 38.17 7	Spectrometry (LC/MS/MS) , 10.1007/8904 (2012) 205.
Methylmalonic acid -d3	120.00	76.00	15 1	2	16 Homocystiene MMA_120416_ PA 2_016 17 Homocystiene MMA_120416_ PA 2_017	EXT LOQ EXT LQC	QC QC	0.74 0.74	169 0.74 197 0.74	822 879	0.205 0.225	50.00 51.06 102 75.00 63.38 84	3. http://www.sigmaaldrich.com/catalog/search?term=516-05-
					18 Homocystiene MMA_120416_ PA 2_018	EXTLQC	QC	0.74	191 0.74	872	0.219	75.00 59.86 7	
MS Conditions:					19 Homocystiene MMA_120416_ PA 2_019	EXT LQC EXT LQC	QC	0.74	226 0.74	905	0.250	75.00 79.38 10	
MS System : Xevo TQD					20 Homocystiene MMA_120416_ PA 2_020 21 Homocystiene MMA_120416_ PA 2_021	EXTLQC	QC	0.74	200 0.74	857	0.230	75.00 70.66 94 75.00 69.14 92	N&focu s=product.
		22 Homocystiene MMA_120416_ PA 2_022	EXT LQC	QC	0.74	203 0.74	867	0.234	75.00 69.20 92	4. https://pubchem.ncbi.nlm.nih.gov/compound/L-Homocysteine.			
Mode : ESI +Ve & -Ve					23 Homocystiene MMA_120416_ PA 2_023 24 Homocystiene MMA_120416_ PA 2_024	EXT MQC EXT MQC	QC QC	0.74 0.74	379 0.74 365 0.74	848 906	0.447 0.403	200.00 205.46 102 200.00 177.03 8	.7 I I I I I I I I I I I I I I I I I I I
Cone gas : 50 L/Hr					25 Homocystiene MMA_120416_ PA 2_025	EXT MQC	QC	0.74	377 0.74	876	0.430	200.00 194.39 9	.2
					26 Homocystiene MMA_120416_ PA 2_026	EXT MQC	QC	0.74	402 0.74	916 042	0.439	200.00 200.12 10 200.00 173.18 8	
Desolvation Temp : 300 °C					27 Homocystiene MMA_120416_ PA 2_027 28 Homocystiene MMA_120416_ PA 2_028	EXT MQC EXT MQC	QC	0.74	393 0.74	943 905	0.397	200.00 175.18 8	
Source Temp : 150 °C					29 Homocystiene MMA_120416_ PA 2_029	EXT HQC	QC	0.74	620 0.74	872	0.712	400.00 374.25 93	The authors would like to thank Dr. Sujay Prasad (Chairman of Anand
					30 Homocystiene MMA_120416_ PA 2_030 31 Homocystiene MMA_120416_ PA 2_031	EXT HQC EXT HQC	QC QC	0.74 0.75	617 0.74 602 0.75	860 845	0.717 0.712	400.00 377.90 94 400.00 374.68 93	Diagnostics) and Mr. Prem Pal (Manager, Analytical Lab) for providing
Desolvation gas flow : 600 L/hr					32 Homocystiene MMA_120416_ PA 2_032	EXT HQC	QC	0.74	652 0.74	882	0.739	400.00 392.00 9	
Data processed through TargetLynx application manager.					33 Homocystiene MMA_120416_ PA 2_033 34 Homocystiene MMA_120416_ PA 2_034	EXT HQC EXT HQC	QC	0.74	601 0.74	895	0.671	400.00 348.51 8	Reference standards and patient serum & plasma samples.
					34 Homocystiene MMA_120416_ PA 2_034	ENTINUC		0.74	640 0.74	0/1	0.755	400.00 389.34 9	

© 2016 Waters Corporation

Contact details: PMN_Rajesh@waters.com