# Monitoring the Critical Quality Attributes of Antibody Drug Conjugates (ADCs) as Part of Biosimilar **Development: Case Studies of ado-trastuzumab emtansine**

Liuxi Chen<sup>1</sup>; Min Du<sup>1</sup>; Henry Shion<sup>1</sup>; Ying-Qing Yu<sup>1</sup>; Lan Wang<sup>2</sup>; Kai Gao<sup>2</sup>; Weibin Chen<sup>1</sup> <sup>1</sup>Waters Corporation, Milford, USA; <sup>2</sup>National Institutes for food and drug control, Beijing, China.

## INTRODUCTION

The recent clinical success of antibody-drug conjugates (ADCs) has invigorated research in the biotherapeutic field.<sup>1</sup> A broad variety of methods to functionalize antibodies with various payloads are currently available. The conjugation methods can influence the sample heterogeneity, and therefore impact the pharmacokinetic, safety and therapeutic efficacy of the product. As a result, it has become evident that quality attributes such as the site of modification and the drug-to-antibody ratio (DAR) need to be controlled to meet more stringent requirements for medical applications. In this study, we investigated the utility of an integrated high resolution analytics platform, consisting of a new IMS Q-Tof mass spectrometer and a targeted informatics system, to understand the critical quality attributes of Lysine-conjugated ADCs.

Figure 1. Surface exposed lysine residues on IgG1.



Figure 2: ADC Tmab-DM1 Structure Illustration

# **METHODS**

#### Liquid Chromatography

System:	ACQUITY UPLC H-Class Bio System
Detector:	ACQUITY UPLC Tunable UV (TUV)
Column:	ACQUITY UPLC BEH300 C18, 300Å, 1.7 μm, 2.1 mm x 100 mm

#### Mass Spectrometry

Mass Spectrometer: Xevo G2-XS QTof

Acquisition mode: MS

#### Informatics

UNIFI Scientific Information System (Waters Corporation)

Data Analysis Type in UNIFI:

- **1.** Peptide Mapping (MS<sup>E</sup>) workflow
- 2. Accurate Mass Screening (MS<sup>E</sup>) workflow

#### Sample Information

Multiple Antibody drug-conjugated samples (Tmab-DM1, Figure 1, 2) were denatured, alkylated, and digested by Asp-N endoproteinase, Leucine enkephalin (LeuEnk) was added to each sample at a final concentration of 50 fmol/µl as an internal standard.



Figure 3. UPLC H-Class Bio with Xevo G2-XS system controlled by UNIFI Informatics software.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

# **RESULTS AND DISCUSSION**

### **DAR Measurement on Intact ADC**

The distribution of the drug load is determined by MS intact analysis. The deconvoluted mass spectra contain 8 major peaks with mass difference of 957 Da between adjacent peaks, which is in agreement to the mass of covalently linked DM1 drug with one MCC linker. In both the innovator and candidate biosimilar ADCs, 8 major peaks correspond to Tmab with 0-7 DM1 drug and linkers respectively (label as +0 drug, +1 drug, etc). The less abundant peaks right next to the major peaks with 219 Da, which attributes to the unreacted linkers that modified the antibody but do not react with DM1.



### **Conjugation Site Identification**

Figure 6. LC/MS<sup>E</sup> chromatogram (BPI) of tryptic peptide mapping analysis for Tmab vs Tmab-DM1 in comparison mode.



Figure 7. LC/MS<sup>E</sup> chromatogram (BPI) of Asp-N peptide mapping analysis for Tmab vs Tmab-DM1 in comparison mode.



#### Table 2. Numbers of conjugation sites identified in different regions of Tmab using DDA and MS<sup>E</sup> methods.

Region	Total # of Lys	Trypsin		Asp N	
		DDA	MSE	DDA	MSE
Variable Fab	12	9	10	6	8
Constant Fab	14	13	13	9	7
Fc	20	16	17	4	6
Total	46	38	40	19	21

Figure 8. MS/MS spectra to confirm conjugations sites for positional isomers for Asp-N peptide <sup>224</sup> DKTHTCPPCPAPELLGGPSVFLFPPKPK<sup>251</sup>







# **ADC Peptide Analysis Workflow**



Table 1. The enzyme of choice for different quantification purposes. Trypsin digest was used to calculate relative abundance of conjugated peptides, while Asp-N digest was used to determine the relative site occupancy of individual site.



Figure 5. ADC Peptide level analysis identification and the

quantification workflows. Lysine-

conjugated ADC and unconjugated

control mAb were digested by

trypsin and Asp N respectively,

followed by MS<sup>E</sup> and DDA modes.

UNIFI Peptide Mapping workflow

was used to identify the conjugated

peptides and pinpoint the

conjugation sites. The same set of

MS<sup>E</sup> data were further analyzed

using Accurate Mass Screening

workflow in UNIFI to quantify the

relative site occupancy and relative

abundance of conjugated peptides

across different samples.



# **Vaters** THE SCIENCE OF WHAT'S POSSIBLE.





The CID fragmentation of the Tmab-DM1 generates a signature fragment ion (m/z 547.2, charge +1), commonly for all conjugated peptides. The signature fragment ion corresponds to a partial drug fragment.

# **Relative Site Occupancy Quantification**

UNIFI Accurate Mass Screening Workflow enables the site occupancy quantitation of ADC peptides (from AspN digestion).<sup>2</sup>



Figure 8. Site occupancy quantitation and cross sample comparison (3) ADC sample with duplicated injections were compared).

-> Peptide 2:D13 with 0 and 1 conjugation was illustrated here. The site occupancy ratio was automatically generated during the data processing. The diastereomeric peaks are combined in the relative% calculation. Trending plot is displayed for 2:D13 peptides across multiple ADC samples. These ADC samples were prepared via different conjugation methods.

# CONCLUSIONS

- 1. For Tmab-DM1, 80 out of 92 conjugation sites were observed.
- 2. UNIFI provided automated workflow for:
  - In-depth primary structure characterization of lysine-conjugated ADC.
  - Site specific localization of ADC conjugation (Peptide Mapping Workflow).
  - Quantification of relative site occupancy (Accurate Mass Screening Workflow)
  - While this presentation has focused on lysine-conjugated ADCs, these UNIFI workflows are directly applicable to other classes of ADC biotherapeutics.

#### Reference:

Wang L, Amphlett G, Blattler WA, Lambert JM, Zhang W, Protein Sci. 2005 Sep; 14(9): 2436-46. Waters Application Note (PN = 720005603). "Automated Quantitative Analysis of Antibody Drug Conjugates Using an Accurate Mass Screening Workflow in the UNIFI Scientific Information