## **ULTRA LOW-LEVEL DETECTION OF PERFLUOROALKYL SUBSTANCES USING LC MS/MS**

PEEK tubes

Lauren MULLIN<sup>1</sup>, Jennifer BURGESS<sup>1</sup>, Gordon FUJIMOTO<sup>1</sup>, <u>Sarah DOWD<sup>1</sup></u> 1.Waters Corporation, Milford, MA 01757

### **INTRODUCTION**

Perfluoroalkyl substances (PFASs) encompass a range of fully fluorinated alkyl compounds, typically with an anionic end group. As a result of their widespread use and subsequent leaching from materials, they have been found in various environmental and biological samples. Analytical approaches require increasingly lower levels of detection for risk assessment and legislative purposes. Additionally, PFASs are ubiquitous in the laboratory environment, namely in Teflon components used in analytical equipment, thus requiring special considerations in the LC system configuration. Here we present detection results for various chain length carboxylic and sulfonic perfluorinated acids, using enhanced MS/MS configuration with the Xevo TQ-S for ultimate sensitivity and a ACQUITY I-Class UPLC modified with the Waters PFC Analysis Kit to remove laboratory contamination sources.

The LC analysis was performed using a modified system comprised of PEEK tubing to replace Teflon tubing for solvent lines. Also, a BEH C18 isolator column and pre-cut stainless steel tubings were placed following solvent mixing and ensured that PFAS contamination was retained prior to injection of the sample, These modifications are captured in Figure 1A and B, and described elsewhere<sup>1,2</sup> MRMs were established and collision energies optimized for 17 native and isotopically labeled PFASs (Table 1) typically analyzed in biological and environmental samples. Standards were prepared over three orders of magnitude, and used as calibration curves. Tap and river water samples were filtered and extracted according to a previously described method<sup>1</sup>.



**MS Conditions** MS System: Waters Xevo TO-S Acquisition Polarity: ESI-Capillary voltage (kV): 0.4 Source temperature (°C): 150 Desolvation temp. (°C): 500 Cone gas flow (L/hr): 150 Desolvation gas flow (L/hr): 1000

Figure 1: Commonly monitored PFASs featured in this study.

#### **METHODS**

55			
MeOH 2 mM	ammoniu	um acetate	
1 ammonium	acetate		
Flow Rate	%A		%В
0.65	90.0		10.0
0.65	90.0		10.0
0.65	0.0		100.0
0.65	0.0		100.0
0.65	90.0		10.0
0.65	90.0		10.0

Table 1: RTs and	optimized	MRMs,	collision	energies	(CE)	and	cone	voltages

Compound	Formula	MRM (Quant. in Bold)	RT (min.)	CE (V)	Cone Voltag
PFBuA	$C_4HF_7O_2$	212.9 > 168.9	0.76	15	30
		212.9 > 212.9		2	30
PFPA	$C_5HF_9O_2$	262.9 > 262.9	1.94	2	30
		262.9 > 218.9		10	30
PFBuS	$C_4HF_9O_3S$	298.9 > 79.9	2.25	35	20
		298.9 > 98.9		25	20
PFHxA	$C_6HF_{11}O_2$	312.9 > 268.9	2.84	5	30
		312.9 > 118.9		20	30
PFHpA	$C_7HF_{13}O_2$	362.9 > 318.9	3.42	10	30
		362.9 > 168.9		15	30
PFHxS	$C_6HF_{13}O_3S$	398.9 > 79.9	3.49	40	30
		398.9 > 98.9		40	30
		398.9 > 118.9		30	30
PFOA	$C_8HF_{15}O_2$	412.9 > 368.9	3.84	10	30
		412.9 > 168.9		20	30
		412.9 > 218.9		15	30
PFNA	$C_9HF_{17}O_2$	462.9 > 418.9	4.18	10	30
		462.9 > 218.9		10	30
PFOS	$C_8HF_{17}O_3S$	498.9 > 79.9	4.20	45	30
		498.9 > 98.9		45	30
		498.9 > 229.9		35	30
PFDA	$C_{10}HF_{19}O_2$	512.9 > 218.9	4.46	20	30
		512.9 > 268.9		15	30
		512.9 > 468.9		10	30
PFUnDA	$C_{11}HF_{21}O_2$	562.9 > 518.9	4.69	15	30
		562.9 > 268.9		20	30
		562.9 > 318.9		20	30
PFDoDA	$C_{12}HF_{23}O_2$	612.9 > 568.9	4.90	15	35
0.571.4		612.9 > 168.9		30	35
PFTIA	$C_{13}HF_{25}O_{2}$	662.9 > 618.9	5.07	15	35
		662.9 > 318.9		20	35
		662.9 > 368.9	F 22	20	35
PFIEIA	$C_{14}\Pi F_{27}O_2$	712.9 > 668.9	5.22	15	35
13 0 0 5 0 0		/12.9 > 218.9	4.00	25	35
C PFOS		502.9 > 79.9	4.20	45	30
13		502.9 > 98.9		35	30
C PFOA		416.9 > 168.9	3.84	18	30
13 -		416.9 > 1/1.9		16	30
<sup>13</sup> C PFDA		514.9 > 269.9	4.46	15	30
		514.9 > 219.9		20	30

#### **RESULTS AND DISCUSSION**

The Xevo TQ-S MRM method showed exceptional sensitivity for the analytes. Measured limits-of-detection (based on peak-to-peak signal to noise measurements of 1:3) for the analytes were below 0.07 ng/mL in solvent standards. In the case of the most commonly monitored PFAS analyte, perfluorooctanesulfonate (PFOS), the limit-of-detection was 0.0125 ng/mL (Table 2 and Figure 3). These low limits of detection highlight the notable sensitivity of this platform for PFASs analysis. Analysis of water samples was performed and PFOS was found in both samples (Figure 4A and B).

Table 2: LOD/Qs, Linearity (R<sup>2</sup>) and linear dynamic range for selected PFASs monitored. Results for PFOS are highlighted in bold.

Compound	LOD (ng/mL)	LOQ (ng/mL)	R <sup>2</sup>	Range (ng/mL)
PFBuS	<0.0125	<0.0125	0.9999	0.0125 to 6.25
PFHxA	<0.0125	<0.0125	0.999	0.0125 to 12.5
PFHpA	<0.0125	<0.0125	0.999	0.0125 to 12.5
PFHxS	<0.0125	<0.0125	0.999	0.0125 to 12.5
PFNA	<0.0125	<0.0125	0.997	0.0125 to 12.5
PFOS	0.0125	0.0625	0.999	0.0625 to 12.5
PFDA	<0.0125	<0.0125	0.996	0.0125 tp 12.5
PFUnDA	0.0125	0.0625	0.992	0.0625 to 12.5
PFDoDA	<0.0125	<0.0125	0.998	0.0125 to 12.5
PFTiA	0.0625	0.0625	0.997	0.0625 to 12.5
PFTeTA	0.0625	0.125	0.991	0.125 to 12.5

# **Naters** THE SCIENCE OF WHAT'S POSSIBLE.



Figure 3: LOD for PFOS in solvent standard, and calibration curve covering 2.5 orders of magnitude.



Figure 4: PFOS detected at an on-column concentration 0.1 ng/mL in tap water (A), and on-column 2.1 ng/mL in river water. Samples were subjected to a 500x preconcentration factor. Ion ratio deviations were calculated for the from the mean of the analytical standards

#### CONCLUSIONS

- Accurate PFAS analysis free of laboratory background contaminants can be achieved using the Waters PFC Analysis Kit, adaptable to any ACOUITY UPLC system
- Sub-ng/µL levels of detection for solvent standards are easily attained with the developed MRMs on the highly sensitive Waters Xevo TQ-S

#### REFERENCES

- 1. Lee P. et al. Waters APNT10090111
- 2. Mullin L. et al. Waters APNT134893168