

ANALYSIS WORKFLOW FOR QUANTITATIVE PROTEOMICS, EMPLOYING TRIPLEX DIMETHYL LABELLING AND ION MOBILITY ASSISTED DATA INDEPENDENT ACQUISITION

Andrew Collins^{1,2}, Antony McCabe^{1,2}, Christopher J Hughes³, Ian Morns⁴, Johannes PC Visser³, Andrew R Jones^{1,2}, Jan Claereboudt⁵,
¹ Institute of Integrative Biology, University of Liverpool, UK; ² OMIC ANALYTICS LTD, Liverpool, UK; ³ Waters Corporation, Wilmslow, UK; ⁴ Nonlinear Dynamics Limited, Newcastle upon Tyne, UK; ⁵ Waters Corporation, Zellik, Belgium

INTRODUCTION

Stable isotope labelling can be routinely applied in LC-MS proteomics, for accurate and reproducible quantitative profiling on a large scale. Labels are incorporated metabolically, enzymatically, chemically or by stable isotope labelling. Here we demonstrate a novel informatics processing pipeline for data sets generated using dimethyl chemical labelling, applied in triplex. This method has the advantage of being broadly applicable to any sample type, and has quantitative reproducibility close to that achievable with metabolic labelling. We demonstrate the benefits of using software able to incorporate retention time (t_r) alignment and profiling, as well as profiling of ion mobility (IM) drift times (t_d) to increase confidence in peptide quantification and sensitivity.




Figure 1. Progenesis QI for proteomics/Proteolabels workflow.

METHODS

Sample preparation

Triplex dimethyl samples were prepared as previously described¹ and the proteomes (HeLa, yeast and *E.coli*) pooled and labeled as shown in Figure 2.

LC-MS conditions

Nanoscale LC separation of tryptic peptides was conducted with a trap column configuration using a M-class system and a 90 or 120 min gradient from 5-40% ACN (0.1% FA) at 300 nL/min using a BEH 1.7 µm C18 reversed phase 75 µm x 20 cm nanoscale LC column. MS data were acquired in triplicate in data independent analysis mode in ion mobility enabled data independent analysis mode (LC-IM-DIA-MS) using a Synapt G2-Si instrument.

Informatics

The LC-MS peptide data were aligned, peak detected and searched with Progenesis QI for proteomics using a reviewed UniProt protein sequence databases. Quantitative analysis of the peptides and protein grouping was conducted with Proteolabels.




Figure 2. Experimental designs (top: proteome unique; bottom: HeLa 100:100:100, yeast 50:25:10, and *E. coli* 50:75:90, 'light' : 'intermediate' : 'heavy', respectively).

RESULTS

Multidimensional data acquisition and analysis

Example data are shown in Figure 3, illustrating a two-dimensional representation of the data (intensity vs. m/z) and inset two three dimensional counterparts of the same set of triplets, showing intensity as a function of t_r and m/z, and intensity as a function of t_d and m/z. Progenesis QI for proteomics utilizes drift time to confirm ion alignment across multiple runs and to correlate DIA precursor and product ions, whereas Proteolabels uses the same information for pair/triplet detection and scoring.

