HDX-MS FOR ASSESSING COMPARABILITY BETWEEN INNOVATOR AND BIOSIMILAR BIOTHERAPEUTICS

Jing Fang¹, Ying Qing Yu¹, Alain Beck², Weibin Chen^{1,} and <u>Asish Chakraborty¹</u> ¹Waters Corp., Milford, MA, US ; ²Center of Immunology Pierre Fabre, Saint-Julien-en-Genevois, France



Figure 1: Waters integrated HDX-MS solution addresses each step of HDX-MS workflow for global (Intact), local (peptide) and residue (amino acid) levels analysis. The Enzymate along with back pressure regulator (BPR) enable high pressure pepsin digestion of proteins. DynamX 3.0: Industry leading HDX-MS informatics for automated processing of global (Intact), local (peptide), and residue (AA, ETD) levels of HDX MS data.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

Waters THE SCIENCE OF WHAT'S POSSIBLE.



Figure 6. Comparison of deuterium incorporation of innovator and biosimilar samples. Representative deuterium incorporation profiles of regions (residues 244-255, 245-254, 245-255, 285-303) in CH2 domains show minute difference. The red line represents the data from biosimilar product, the green, cyan and blue lines represent the data from the three batches of innovator samples. The experiments have been repeated in triplicate runs. The location of the region displayed minor difference among biosimilar and innovator sample is colored in red in the model structure of IgG1 (PDB: 1HZH). Glycosylation is shown in blue. Met255 is circled and shown in stick.

CONCLUSION

- Higher order structure of a protein is an important factor for its biological function.
- HDX-MS is reproducible and reliable method, and is able to localize the conformational changes.
- HDX-MS study showed that the conformation for the biosimilar product is highly comparable with the reference product overall. No distinct difference was observed in the CDR regions.
- Minute differences were observed in CH2 domain close to glycosylation area and attributed to minor differences in glycosylation between biosimilar and reference product.
- Our study shows great promise in adapting HDX-MS into biosimilar drug development process.