ADVANCES IN TARGETED OMICS QUANTITATION USING A NOVEL SCANNING QUADRUPOLE DIA METHOD

CJ Hughes¹, *K* Richardson¹, *J* Wildgoose¹, *LA* Gethings, *JPC* Vissers, *M* Green¹, *R* Chapman¹, *A* Grzyb¹, *P* Harapanahalli¹, *K* Craven¹ ¹Waters Corporation, Wilmslow, United Kingdom; ²Waters Corporation, Milford, MA

INTRODUCTION

Targeted LC-MS based assays are increasingly applied in the post-discovery omics area with emphasis on validation, the first of many phases in translational analyses, or in studies that are aimed at gaining the understanding of biological systems, drug development and treatment. Context is driving current omics experiments, thereby driving the development of LC-MS acquisition methods that can provide both qualitative and quantitative information in a single experiment. A novel Data Independent Acquisition technique that satisfies both these requirements will be demonstrated.

In the method, a low-resolution quadrupole mass filter is scanned repetitively and both precursor and MS-MS data are acquired at spectral rates approaching 2000 spectra per second. This method produces a high duty-cycle, highly specific and unbiased two-dimensional data set which can be viewed and processed using readily available informatics. The mode of operation has been implemented on a benchtop quadrupole / ToF mass spectrometer and has been applied to targeted lipidomics and proteomics experiments involving the quantitation of both lipids and proteins in human plasma samples derived from normal, obese and diabetic donor cohorts.

METHODS

Sample preparation

Standards. A proteomics standard comprised a four protein digest mixture (Waters Corporation) spiked into 1 μ g/ μ L *E. Coli* tryptic digest (Waters Corporation) matrix at a maximum concentration of 20 fmol/ μ L and then serially diluted to a lowest concentration of 100 amol/ μ L. A pre-mixed synthetic lipid standard mixture (Avanti Polar Lipids, Inc.) served as a lipidomics standard that was serially diluted in a human plasma extract.

Biological samples. Proteins and lipids were extracted from human plasma (Innovative Research Inc), which originated from 6 control, 6 obese and 6 diabetic patients. The protein extracts were prepared with 1% RapiGest SF (Waters Corpora-tion) prior to reduction, alkylation and overnight digestion with trypsin¹. Lipid extractions were performed as previously described². The extracts were centrifuged for phase separation and the lower fraction collected for LC-MS analysis.



Bioinformatics

SONAR[™] DIA data were processed using Progeneseis QI, Progenesis QI for proteomics (Nonlinear Dynamics) and ProteinLynx Global Server (Waters Corporation) using optimized threshold and search parameters. Quantitative analysis was performed with Skyline (University of Washington) using libraries derived from PQI, PQIp and PLGS compound or protein database searches.

RESULTS

Specificity

Figure 3 represents typical SONARTM data and demonstrates that the data format is the same as other multi dimensional datasets, eg ion mobility; hence, exhibits improved specificity. When viewing these data, the drift time axis is replaced by quadrupole mass. The gain in acquisition specificity *vs.* a non resolving first mass analyser DIA (MS^E) based method is described in Figure 4.



Figure 3. SONAR[™] DIA ToF vs. quadrupole m/z data, showing product ions (vertical bands) from peptides eluting over a 1 min window and the quadrupole sweep (diagonal line).



Discovery omics

The ability to search and quantify (relatively) SONAR[™] DIA omics data is shown in Figures 6 and 7, highlighting disease specific compounds, either lipids, peptides or proteins that show similar regulation trends.





Figure 6. Discovery lipidomics DIA 2DMS analysis illustrating identification, compound distribution, and grouped lipid quantitation/hierarchical clustering.





Figure 1. Xevo G2-XS QTof mass spectrometer instrument schematic.

LC-MS conditions

Label-free LC-MS was used for qualitative and quantitative peptide analyses. Experiments were conducted using a 90 min gradient from 5 to 40% acetonitrile (0.1% formic acid) at 300 nL/min using a M-Class system (Waters Corporation) and a HSS 1.8 μ m C18 reversed phase 75 μ m x 20 cm nanoscale LC column.

For lipid identification, the LC-MS experiments consisted of a 20 min gradient from 3 to 40% isopropanol:methanol (10 mM ammonium formate) at 500 μ L/min using a ACQUITY UPLC system. Here, a BEH 1.7 μ m C8 reversed phase 2.1 x 10 cm LC column was used.

A Xevo G2-XS QToF (Waters Corporation), Figure 1, was operated in SONARTM 2DMS mode. The optimized quadrupole window and the other parameters employed for lipid and peptide analyses are described in Figure 2.



Expt	lon Mode	Quad Scan (Da)	Quad Window (Da)	ToF Scan (Da)	Function Integration Time (s)	Fn1 CE (V)	Fn2 CE (V)
Lipidomics	+	500 to 1200	10	50 to 1000	0.1	6	20 to 50
Lipidomics	-	500 to 1200	10	50 to 1000	0.1	6	25 to 55
Proteomics	+	400 to 900	23	50 to 2000	0.5	6	14 to 40

Figure 2. SONARTM acquisition method and DIA acquisition parameters used in the different experiments.



Concept

Shown in Figure 5 is the conceptual use of the SONARTM 2DMS DIA data whereby precursor and product ions are extracted in multiple dimensions, *i.e.* retention time and first (quadrupole) and second (ToF) mass analyser m/z. As illustrated, the method is amenable to both metabolomics/lipidomics (A) and proteomics applications (B). The following sections will demonstrate that this data can be utilized for both qualitative and quantitative analysis of (multi) omics data.





Figure 5. Multiple dimension ion extraction using retention time, quadrupole (precursor) m/z and ToF (product ion) m/z for a (A) synthetic lipid mixture spiked into human plasma extract and (B) 4 protein digest mixture spiked and serially diluted into a tryptic digest of cytosolic E.coli.





Figure 7. Discovery proteomics DIA SONAR[™] analysis strategy showing peak detection, identification, grouped quantitative peptide analysis, and protein quantitation/resolution for pathway associated proteins (apolipoprotein plasma/serum complement).

CONCLUSIONS

- SONAR[™] DIA acquisition provides multi dimensional data sets exhibiting improved specificity and over other DIA methods
- SONAR[™] data acquired in multi omics experiments shows excellent qualitative and quantitative characteristics and can be applied in discovery, targeted and translational omics studies
- Application of SONAR[™] detected and quantitatively confirmed known diabetes/obesity proteins and lipid indicators

References

- 1. M'Basu et al. Advances in Quadrupole and Time-of-Flight Mass Spectrometry for Peptide MRM based Translational Research Analysis. Proteomics. 2016 Aug;16(15-16):2206-20.
- 2. Sarafian et al. Objective Set of Criteria for Optimization of Sample Preparation Procedures for Ultra-High Throughput Untargeted Blood Plasma Lipid Profiling by Ultra Performance Liquid Chromatography-Mass Spectrometry. Anal Chem 2014; 86; 5766-5774.
- van der Ham et al. Plasma Apolipoprotein CI and CIII Levels Are Associated With Increased Plasma Triglyceride Levels and Decreased Fat Mass in Men With the Metabolic Syndrome. Diabetes Care 2009; 32; 184-186.