A Streamlined Analytical Platform with RFMS Glycan GU Scientific Library Searching for Glycan **Comparability Assessment of Innovator and Biosimilar Infliximab**

William R. Alley, Jr., Ying Qing Yu, and <u>Min Du</u> Scientific Operations, Waters Corporation, Milford, MA 01757

INTRODUCTION

To analyze N-linked glycans derived from biotherapeutics, we have developed a streamlined UPLC-FLR-MS analytical platform. This platform utilizes *Rapi*Fluor-MS[™] from Waters Corporation to rapidly label glycans to improve both FLR and MS sensitivities. This allows lower-abundance glycans to be successfully detected. Following data collection and processing, raw retention times are converted to Glucose Units (GUs), which are then searched against the RFMS Glycan GU Scientific Library for structural assignment, which simplifies that task. To confirm the GU structural assignments, accurate MS data (sub 2-ppm) is used.

We have applied this analytical approach to compare innovator and biosimilar infliximab and we pay particular attention to the possible immunogenic glycans, those with Nglycolylneuraminic acid and a-linked galactose-galactose pairs. In general, we have found that the major glycan structures on both mAbs are very similar, but there appear to be some slight differences in some lesser abundant glycans based on lower % relative fluorescence peak areas.

EXPERIMENTAL WORKFLOW



Figure 1. Outline of the sample preparation workflow used in this study.



UNIFI WORKFLOW



INSTRUMENT CONDITIONS

LC Conditions

LC System: LC Column:

Column Temp.: Sample Temp.: Fluorescence Detection: Ex 265 nm/Em 425 nm

ACQUITY UPLC H-Class Bio ACQUITY BEH Glycan Amide, 2.1 x 150 mm, 130 Å pore size, 1.7 µm particles 60° C 10° C

Gradient Conditions

Mobile Phase A: 50 mM ammonium formate, pH = 4.5Mobile Phase B: Acetonitrile

Time (min)	Flow Rate (mL/min)	%A	%В
0.0	0.400	25.0	75.0
35.0	0.400	46.0	54.0
36	0.300	70.0	30.0
41	0.300	70.0	30.0
42	0.200	25.0	75.0
47	0.400	25.0	75.0
55.0	0.400	25.0	75.0

MS Conditions

MS System: Analyzer Mode:

Capillary Voltage: Cone Voltage: Source Temp.: Desolvation Temp.: *m/z* Range Scan Rate: Instrument Control: Scientific Library:

XEVO G2-XS QToF ESI+, sensitivity mode, ~30k resolution 3.0 kV 80 V 120° C 250° C 500-2000 2 Hz UNIFI 1.8.2 RFMS Glycan GU Scientific Library

CREATING THE RFMS GLYCAN GU SCIENTIFIC LIBRARY



This figure shows how the library was created. Glycans were first labelled with RFMS and analyzed by UPLC to determine their Glucose Unit (GU) value. The glycans were then subjected to a battery of exoglycosidase digestions to selectively remove specific monosaccharides to determine their linkages and locations. The structures were verified through accurate mass measurements. This work was done jointly with Mark Hilliard and Pauline Rudd of NIBRT.

LIBRARY

glycans sample

EXAMPLE LIBRARY ENTRY



CALIBRATING RETENTION TIMES



QUANTITATION REPRODUCIBILITIES



THE RFMS GLYCAN GU SCIENTIFIC

•177 glycans tailored to biopharma Created from: etanercept, infliximab, trastuzumab, erythropoietin, NIST mAb, human IgG, mouse IgG, bovine fetuin, yeast invertase •Contains high mannose, complex, and hybrid

•Contains glycans with N-glycolylneuraminic acid and N-acetylneuraminic acid, and a-linked galactoses •GU values were averaged from at least 10 injections per

The RFMS GU Glycan Scientific Library search is based on glucose units (GU) rather than raw retention times. We have found GU values to be more stable across different instruments. To convert retention times to GUs, a dextran ladder labeled with RFMS is first analyzed as a separation calibration. The resulting equation is then used to convert retention times of unknown samples.

The injection-to-injection and preparation-to-preparation reproducibilities were monitored through triplicate injections of one sample and triplicate preparations of all samples. For this particular glycan of interest, the %RSDs were less than 2%.

RESULTS

COMPARING INNOVATOR AND BIOSIMILAR INFLIXIMAB



In total, we identified 23 mass-confirmed glycans present on the innovator mAb and 21 on the biosimilar. All of the glycans on the biosimilar were present on the innovator molecule.

COMPARING GLYCAN CLASSES



While some slight differences were observed for different classes of glycans (sialylated and biantennary glycans being slightly elevated for the biosimilar), overall the different classes were present at comparable levels on each mAb.

IDENTIFYING DIFFERENCES IN ABUNDANCE LEVELS



This mirror plot shows that most of the higher abundance glycans are present at very comparable levels. However, some differences were noted for some lower abundance glycans, particularly those with N-glycolylneuraminic acid and a-linked galactoses.

GLYCANS WITH N-GLYCOLYLNEURAMINIC ACID

Example Library Search

Component name	Structure	Expected GU	ΔGU	GUStd Dev	Expected m/z	∆m/z	Mass Confirmed
F(6)A2(6)G(4)1Sg(6)1	<	8.1627	0.0011	0.0500	1122.4384	0.0005	True
M4A1	₽-{₽-₽-₽	8.1557	0.0081	0.0500	1122.9486	0.5107	False
F(6)A2[3]BG(4)1S(6)1		8.2239	0.0601	0.0500	1215.9806	93.5427	False
FA2G2F1	╃₋ <mark>҈╸<mark>╸</mark>╺ ┍─■╺</mark>	8.2365	0.0727	0.1000	1122.9486	0.5107	False
F(6)A3G(4)3		8.2393	0.0755	0.1000	1232.4857	110.0478	False
F(6)A2G(4)2Ga(3)1	╺╌ <mark>┊╸╸╸</mark>	8.2457	0.0819	0.1000	1130.9460	8.5081	False
F(6)A2G(4)2Ga1 iso	0 { 0 - ■ - 0 00 - ■ - 0	8.2995	0.1357	0.0500	1130.9460	8.5081	False

A2(6)G(4)1Sg(6)1.



N-glycolylneuraminic acid is not present on human glycans and can induce an immunogenic response. This is an important monosaccharide to monitor. These summary plots show that the complex glycans terminated with N-glycolylneuraminic acids are elevated in the biosimilar product when compared to multiple batches of the innovator drug. These types of glycans were present at low abundances in both mAbs.



aters/ THE SCIENCE OF WHAT'S POSSIBLE.

This figure shows the results of the library search for a GU value of ~ 8.16 \pm 0.2. Within this GU window, only one glycan could be mass confirmed and was measured at sub-2 ppm mass accuracy. This glycan was determined to be F(6)

Summary Plots for Glycans with N-Glycolylneuraminic Acid

Symbols: N-acetylglucosamine; Mannose; OGalactose; Fucose; N-glycolylneuraminic acid; N-acetylneuraminic acid;

Summary Plots for Glycans with a-linked Galactoses



Glycans with a-linked galactoses are nonhuman in nature and their presence in the human body has been shown to result in a severe immunogenic response. Therefore, this is another important group to monitor. This plot shows two glycans with a-linked galactose units were present at higher abundance levels in the innovator samples. Interestingly, we observed a glycan with a-linked galactoses and N-glycolylneuraminic acid that was elevated on the biosimilar mAb.

Summary Plots for the 4 Most Abundant Fucosylated Glycans



Fucose levels can influence the antibody-dependent cell-mediated cytotoxicity of a therapeutic mAb. The overall level of fucosylation was slightly elevated in the biosimilar product and was due its higher level of the F(6)A2[6]G(4)1 glycan, even though F(6)A2was slightly lower. The other more abundance fucosylated glycans were at approximately the same level in both mAbs.

CONCLUSIONS

- We have developed a streamlined platform for glycan analysis, from sample preparation to data collection and analysis using a scientific library for glycan structural determination.
- Overall, the glycan profiles were very similar, especially for the abundant structures.
- Some differences were observed for some lower-abundance glycans, including a complex glycan terminated with an Nglycolylneuraminic acid (elevated in the biosimilar sample) and a hybrid glycan possessing a-linked galactose units, which was observed only on the innovator samples.