

SIMPLE AND EFFECTIVE CLEANUP FOR UPLC-MS/MS
DETERMINATION OF VETERINARY DRUG RESIDUES IN EGG

Sujie Xia1, Kim Van Tran2, Dimple Shah2, Jeremy C. Shia2, Michael S. Young2 Jennifer Burgess2
1Shanghai Institute for Food and Drug Control, 1500 Zhangheng Road, Shanghai 201203, China; 2Waters Corporation, 34 Maple St, Milford, MA 01757, U.S.A
xiasujie@163.com., michael_s_young@waters.com

INTRODUCTION

Veterinary drugs are used in chicken farms to control diseases of laying hens. However, these compounds can be transferred to and accumulate in the eggs. The presence of veterinary drug residues in eggs is a potential health risk for the consumer because the residual drugs can provoke allergic reactions or induce pathogen resistance to antibiotics used in human medicine [1]. Sixteen representative veterinary drugs from twelve classes, most of which have MRPLs established in USA, EU and/or China, were chosen for this study [2][3]. Figure 1 presents the structures of a subset of these veterinary drugs.

Sample preparation is a challenging task for the multi-residue determination of veterinary drugs in eggs. The analyst must recover a wide variety of drug classes with different physico-chemical properties. Some of the target compounds may bind to proteins or other matrix components. Also, eggs contain high levels of lecithin (phospholipids) and fats; these co-extracted substances can lead to interference in the LC-MS analysis, contamination of the analytical column and other components of the UPLC system, and contamination of the mass spectrometer itself.

In this work, sample extraction, cleanup and analysis methods were developed for UPLC-MS/MS determination of a wide variety of veterinary drugs in eggs. Samples were treated with an acidified acetonitrile/water solvent to precipitate proteins, release bound residues and to extract the veterinary drugs of interest. Then, to remove fats and phospholipids, a simple pass-through cleanup was performed using a novel SPE device, the Oasis PRiME HLB cartridge.




Figure1. Structures of representative compounds from this study

Standard Compounds

Sixteen veterinary drugs from different classes were chosen for this study. Table 1 lists their formula, MWs, and MRPLs established in USA or China.

Table with 5 columns: Compounds, Formula, Monoisotopic MW, Market, MRL(ppb). It lists 16 veterinary drugs including Amprolium, Bacitracin A, Hygromycin B, Nystatin A1, Colistin B, Florfenicol, Flubendazole, Oxolinic acid, Tiamulin, Chlortetracycline, Erythromycin, Lincomycin, Oxytetracycline, Penicillin G, Tetracycline, and Tylosin.

Table 1. Veterinary drugs in this study (bacitracin, colistin and nystatin all contain a mixture of more than two components; one major component was chosen for analysis)

METHODS

Sample Preparation

Extraction: Two grams of homogenized whole chicken eggs were weighed into a 50 mL polypropylene centrifuge tube and 8mL of 0.2% formic acid in 80:20 acetonitrile/water were added. Following a 30s vortex and shaking for 30 mins, the mix was centrifuged at 4500rpm for 10min.

Pass-thru SPE Cleanup: An Oasis PRiME HLB cartridge (3cc,60mg) was mounted on a precleaned vacuum manifold. Cartridge conditioning is not required. 0.5 mL of the supernatant was passed-through the cartridge and collected using 1~2 psi vacuum. 0.2 mL of the collected extract was taken and diluted to 0.6 mL with aqueous 10 mM ammonium formate buffer (pH 4.5) prior to UPLC-MS/MS analysis.

LC and MS Conditions

Chromatography Conditions: System: ACQUITY® UPLC I-Class; Column: ACQUITY UPLC BEH C18, 2.1x100mm, 1.7µm; Column Temp.: 30 °C; Injection Volume: 10µL; Flow Rate: 0.4mL/min; Mobile Phase A: 0.1% formic acid in water; Mobile Phase B: 0.1% formic acid in acetonitrile; Gradient: The Initial composition was 85% A and 15% B. Phase B was increased linearly to 40% in the first 2.5min, and then linear ramp to 95%B in 1.4 min, maintained for 2.3 min, then returned to the initial composition and equilibrated for 2 min.

MS Conditions: System: Xevo® TQ-S; Ionization Mode: ES+ (ES-for Florfenicol); Capillary Voltage (kV): 3.00(2.50 for negative mode); Source Temp. (°C): 150; Desolvation Temp. (°C): 600; Cone Gas Flow (L/Hr): 150; Desolvation Gas Flow (L/Hr): 1000; Collision Gas Flow (mL/Min): 0.15; Nebuliser Gas Flow (Bar): 7.00

LC-MS/MS Parameters

MS/MS Transitions are presented in Table 2. Other LC-MS parameters are presented in Table 3. A typical chromatographic separation is presented in Figure 2.

Table 2. MRM Transition parameters for 16 veterinary drugs. Columns include Compounds, Precursor Ion (m/z), Product Ion (m/z), Cone voltage(V), Collision energy(eV) for two transitions.

Table 2. MRM Transition parameters for 16 veterinary drugs

Table 3. LC-MS Retention times and calibration data. Columns include Compounds, RT (min), LOD (ppb), Linear Range (ppb), and R2.

Table 3. LC-MS Retention times and calibration data




Figure 2. Overlay of quantitative MRM chromatograms of 16 veterinary drugs (at MRL level)

RESULTS AND DISUSSION

Oasis PRiME HLB Cartridge Pass-Through Cleanup

The Oasis PRiME HLB cartridge was evaluated with respect to analyte recovery and phospholipids removal from egg matrix. The total method recoveries ranged from 50-97%. However, the Oasis PRiME HLB cartridge cleanup contributes little to any method recovery losses. As shown in Figure 3, the measured recovery for the SPE cleanup step is better than 80 % for all compounds, with recovery for most compounds greater than 90%.




Figure 3. Recovery data for target veterinary drugs obtained using the Oasis PRiME HLB cartridge clean-up procedure (at 1MRL level)

Whole eggs contain significant amounts of fat and are among the highest sources of dietary lecithin (phospholipids). The total lipid content of chicken egg is about 11% by weight (excluding the shell) and the phospholipids content is about 0.35% [4]. Significant amounts of these potential interfering substances are extracted along with the target drugs in the initial sample preparation extraction step. Greater than 84% of total lipids were removed from the egg extract after pass-through cleanup with the Oasis PRiME HLB cartridge. The cleanup step was even more effective for removal of phospholipids. Figure 4 shows that the Oasis PRiME HLB cartridge cleanup removed greater than 95% of phospholipids from the egg extract.




Figure 4. Effective removal of phospholipids from egg extracts with Oasis PRiME HLB clean-up

Method Recovery and Precision

Recovery studies were carried out at three concentration levels (0.4MRL, 1MRL, 2MRL), six replicates per level. Matrix -matched standard calibration curves were used. Figure 5 shows the results. Recovery was greater than 70% for most target compounds (>70%) except for nystatin and hygromycin. Reproducibility was acceptable (RSD<20%) for all compounds except for hygromycin at 0.4MRL (RSD=34%).




Figure 5. Summary of recovery data (blank eggs samples spiked at 0.4MRL, 1MRL, 2MRL levels (Hygromycin B, Florfenicol, Penicillin G, Nystatin A1 have no corresponding MRLs, so they were studied at 40, 100, 200ppb levels)

CONCLUSIONS

- This analytical method meets the requirement for the simultaneous determination of several classes of veterinary drugs in eggs
- A simple pass-through cleanup procedure using Oasis PRiME HLB cartridge can remove more than 95% phospholipids from egg extracts
- The Oasis PRiME HLB cartridge cleanup procedure provided effective cleanup and good recoveries for the target veterinary drugs in egg
- The ACQUITY® UPLC I-Class coupled with Xevo® TQ-S MS offered good sensitivity for the veterinary drug residues in this study

References

[1]Antonia Garrido Frenich, et.al. Analytica Chimica Acta 661 (2010) 150-160.
[2] The Ministry of Agriculture Bulletin of PRC235, 2002.
[3] https://www.globalmrl.com/