Determination of Nanoparticle Sizes by Turbidity-Debye Method

Bei Niu, Xiangjin Song, Beatrice Muriithi, and Yuehong Xu* Waters Corporation, 34 Maple Street, Milford MA 01757, USA ACS 2017: Poster ANYL 173—2759475 *corresponding author: yuehong_xu@waters.com

THE SCIENCE OF WHAT'S POSSIBLE.

ers

INTRODUCTION

Interest in nanoparticles has risen consistently and rapidly in recent years, and the determination of the sizes of these materials has become crucial. There are various techniques for nanoparticle size analyses, and different techniques provide different information and accuracy. We have developed an analytical method to determine nanoparticle sizes based on Debye theory using turbidity measurements, and compared with other techniques such as dynamic light scattering (DLS), field emission scanning electron microscope (FE-SEM) and calculation from specific surface area (SSA).

THEORY

Attenuance and Transmittance	Turbidity from Dielectric Rayleigh S	catters Debye Approach
$D = -\log_{10} (T)^{[1]} Eq.1$ $T = I/I_0 = \exp(-NC_{ext}L)^{[2]} Eq.2$ $D = 0.434 \times NC_{ext}L Eq.3$	$C_{abs} = 0 \qquad \text{Eq.4}$ $C_{ext} = C_{sca} + C_{abs} = C_{sca} \qquad \text{Eq.5} \qquad \frac{\tau}{c} = \frac{D}{0.434Lc} = 4\pi^4 \left(\frac{1}{2}\right)^2$	The second seco
 D: attenuance, or absorbance for solution without scattering, measured by UV/ Visible spectrophotometer T: transmittance I: transmitted intensity I₀: incident intensity N: number of the same kind of particles per unit volume C_{ext}:extinction cross section of a particle L: path length 	$\begin{array}{c} C_{sca}^{n}: \mbox{ scattering cross section } \\ \tau: \mbox{ turbidity } \\ \lambda_{0}: \mbox{ wavelength } \\ \lambda_{0}: \mbox{ particle dian } \end{array}$	• $1/(Dk_d)_{c=0} = intercept$ • $1/(Dk_d)_{c=0} = intercept$ • $1/(Dk_d)_{c=0} = intercept$ • $c = \frac{\rho_{s,0} \times f_{w,0}}{k_d}$ Eq.9 Eq.9 Eq.9 Eq.9 Eq.9 Eq.10

EXPERIMENTAL

Samples were diluted in appropriate dilution factors in order to achieve linear correlation of $1/Dk_d vs. 1/K_d$ as well as sufficient signals for accurate measurements. All suspensions had a final buffer concentration of ~10mM and a pH of ~8.5 to keep the silica nanoparticles stable for measurements. UV/Visible spectrum was acquired from each suspension from 700nm to 300nm. The wavelength of 300nm is chosen because scattering cross section increases with the decrease of wavelength, and the silica nanoparticles behave as dielectric Rayleigh scatters at this wavelength. The measurements at 700nm are for reference.

RESULTS AND DISCUSSION




Figure 1. Example of UV/visible spectra of the diluted sample

Figure 2. Example Debye-like plot

Figure 3. Impact of particle concentration on measured diameter

Comparison of Nanoparticle Diameters Determined by Different Techniques:




Figure 4. Correlation of particle diameters by different techniques

Advantages and Disadvantages:

- Both turbidity-Debye method (T-D) and dynamic light scattering (DLS) are the measurements of the light scattered by nanoparticles, and thus, are more sensitive to larger agglomerates, if present. As shown in Figure 3, the particle concentration has great impact on the particle diameter measurement. The diameter determination at zero concentration is achievable by the T-D method with excellent correlation coefficient in the Debye-like plot (refer to Figure 2).
- The specific surface area (SSA) method has excellent reproducibility (refer to Table I). However, nanoparticle diameters can be underestimated if they are fused.
- The T-D method had good correlation with the SSA method (refer to Figure 4). The diameters by T-D were 14% greater than those by SSA in this study.
- The DLS method has poor correlation with the SSA method when the intercept is set to zero. The accuracy of the DLS measurement is in question.
- FE-SEM provides the capability of viewing nanoparticle morphology, but the accuracy and reproducibility in the measurements of such small particles is very low.

SUMMARY

An analytical method for the determination of nanoparticle sizes based on Debye theory using turbidity measurements has been developed. The method uses intercept at zero concentration to eliminate the affects of particle-particle interactions, and thus, improves the accuracy of the particle size measurement. In this study, we have also compared the turbidity-Debye method to dynamic light scattering, field emission scanning electron microscopy, and nitrogen sorption. Accurate and reproducible results have been achieved using the turbidity-Debye method.

REFERENCES

McNaught, A. D. and Wilkinson, A. (1997), "Compendium of Chemical Terminology – The Gold Book" (2nd Edition), Royal Society of Chemistry, Cambridge, UK.
 Kerker, M. (1969), "The Scattering of Light and Other Electromagnetic radiation," Academic Press, New York.

[3] http://www.sciner.com/Opticsland/FS.htm

[4] Hale, G.M. and Querry, M.R. (1973), *APPLIED OPTICS*, **12**, 555.

[5] Heller, W. (1965), J. Chem. Phys., 42, 1609.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS