# FAST ANALYSIS OF ISOFLAVONES IN DIETARY SUPLEMENTS

Column Tomp

Jinchuan Yang, Mark E. Benvenuti, Gareth Cleland, Ken Rosnack, Joe Romano Waters Corporation, Milford, MA

# INTRODUCTION


- Twelve major isoflavones were found in plants of soy (Glycine max), red clover (Trifolium pretense), and Kudzu (Pueraria lobata). These hormone-like compounds are often used in remedies to reduce menopausal and postmenopausal symptoms.
- Standard methods for isoflavones in dietary supplements use reversed-phase LC with C<sub>18</sub> columns and ultraviolet and visible light (UV-Vis) spectroscopy for separation and quantitation. The chromatographic run time is over 70 minutes long.
- This work demonstrates the transfer of the USP isoflavone method onto an ACQUITY Arc<sup>™</sup> UHPLC system and an ACQUITY UPLC H-Class system. The benefits of Mass Spectrometry in method transfer and sample analysis are highlighted.

GOAL


Transfer the USP isoflavone HPLC method to


faster UHPLC and UPLC methods to achieve

higher sample analysis throughputs



# **METHODS**







System:	ACQUITY Arc system with 2998 PDA and ACQUITY QDa Detector (Performance)
Software: Column:	Empower 3 CORTECS® C18 3.0 x 100 mm, 2.7µm
	(p/n 186007372)
Column Temp.:	30°C
Mobile phase A:	Water with 0.1% Formic acid
Mobile phase B:	Acetonitrile with 0.1% Formic acid
Inj Vol.:	2.0 µL
Run time:	18.0 min
UV detection:	260 nm
Elution Gradient:	

	Time (Min)	Flow Rate	%A	Curve
		(mL/Min)		
1	Initial	1.08	90	6
2	14.40	1.08	70	6
3	14.50	1.08	10	6
4	15.20	1.08	10	6
5	15.40	1.08	90	6
6	18.00	1.08	90	6

## **UPLC** conditions

System:	ACQUITY UPLC H-Class system with ACQUITY PDA and ACQUITY QDa Detector (Performance)		
Software:	Empower 3		
Column:	CORTECS <sup>®</sup> UPLC C18 2.1 x 75 mm, 1.6µm (p/n 186007094)		

## TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

Contact: Jinchuan\_yang@waters.com

	Colur	nn Temp.:	30°C					
	Mobil	e phase A:	Water with 0.1	% Formic	acid			
	Mobil	e phase B:	ACN with 0.1%	% Formic a	acid			
	Inj Vo	ol.:	0.7 µL					
Run time:		ime:	8.10 min	8.10 min				
	Elutio	on Gradient:						
		Time (Min)	Flow Rate	%A	Curve			
			(mL/Min)					
	1	Initial	0.55	90	6			
	2	6.48	0.55	70	6			
	3	6.53	0.55	10	6			
	4	6.84	0.55	10	6			
	5	6.93	0.55	90	6			
	6	8.10	0.55	90	6			

## **QDa Parameters:**

Polarity:	ES+
Capillary (kV):	0.8
Cone (V):	15
Probe Temperature (°C):	600

Table 1. Masses of isoflavone molecular ions

SIR Mass (Daltons)		
417.1		
447.0		
433.1		
254.9		
285.0		
270.9		
270.9		
503.4		
533.1		
459.1		
489.0		
519.0		
475.1		


# RESULTS

## 1) USP method transfer and optimization

The USP method (isoflavones powder extract) <sup>(1)</sup> was transferred to an ACQUITY Arc system with a CORTECS C18 column (2.7 µm 3 x 100 mm). The USP method's 74 minutes long mobile phase elution program was converted to a 18 minutes program using Waters UPLC Column Calculator. The mobile phase additive was changed from phosphoric acid (0.05%) to formic acid (0.1%), which is a mass spectrometry friendly additive. The column temperature was optimized to 30°C to meet the USP suitability criteria on peak resolution.

This USP method was also transferred to an ACQUITY UPLC H-Class system with a CORTECS UPLC C18 (1.6µm, 2.1 x 75 mm). The total run time in UPLC was shortened to 8.1 min, which is about one tenth of the HPLC run time.

Figure 1. The reference HPLC chromatogram for USP defatted powder soy RS. The run time is 74 min. This chromatogram was provided by






Figure 2. Chromatograms of USP defatted powder soy RS in A) UPLC and B) UHPLC conditions. The peak assignment were confirmed by mass spectrometry data that obtained from a Waters QDa mass detector.

## 2) Method performance

Table 2 shows the UV calibration results for the UPLC method. The square of the correlation coefficients  $(R^2)$  between the responses (peak area ratio) and the standard concentrations (ppm) for all compounds were better than 0.998. The relative standard deviation (RSD) in retention time for all compounds were less than 0.27%.

Table 3 shows the isoflavone results for the NIST 3238 SRM and a comparison to its certified and reference values. Relative difference of less than 14% was obtained for the genistin, glycitin, daidzin, genistein, and glycitein. The accuracy for the daidzein, genistein, and glycitein was also evaluated by a spiking experiment (Table 4). Recovery of 95% to 102% were obtained for these analytes.

Table 2. Isoflavones retention times and their relative standard deviation, calibration equations,  $R^2$ , and linear ranges for the UPLC method.

	Analyte RT Min	RT	Equation	B <sup>2</sup>	Range	
		Min	RSD (%)	Equation	n-	(ppm)
	Daidzin	2.03	0.21	$y = 2.12 \times 10^{-1} X + 3.93 \times 10^{-3}$	0.999	0.075-10
	Glycitin	2.27	0.27	$y = 2.45 \times 10^{-1} X + 4.54 \times 10^{-3}$	0.998	0.075-10
	Genistin	3.07	0.21	$y = 4.56 \times 10^{-1} X + 5.30 \times 10^{-3}$	0.998	0.075-10
	Daidzein	4.76	0.11	$y = 2.61 \times 10^{-1} X + 4.50 \times 10^{-3}$	0.998	0.075-10
	Glycitein	5.15	0.14	$y = 4.97 \times 10^{-1} X + 5.69 \times 10^{-3}$	0.998	0.05-10
	Genistein	6.32	0.09	$y = 4.88 \times 10^{-1} X + 4.87 \times 10^{-3}$	0.998	0.05-10

Daidzein Genistein

method were used in this label claim is unknown

## 3) Analysis of isoflavones in dietary supplements

## ample ng/kg)

Daidzin Glycitin Genistin Malonyl daio Acetyl daidz Acetyl glycit Malonyl gen Daidzein Glycitein Acetyl genis Genistein Total isoflav Label claim

 
 Table 5. UPLC-UV method determined Isoflavone content in dietary
supplements and their label claim values.

4) Mass detection helps to reduce co-eluting errors in UV detection It is important to note that due to the complex sample matrix, co-elution often occurs in isoflavones analysis. Figure 3 shows the UV peak fronting for the daidzein in NIST 3238 SRM. The quantitation error that is caused by

# THE SCIENCE OF WHAT'S POSSIBLE.

Table 3. Comparison of determined isoflavone values to the certified and reference values of NIST 3238 SRM (UPLC method).

	NIST value (mg/kg)	Determined value (mg/kg)	Rel. difference	
nistin	12700±530	11432	-10%	
/citin	3760±180	3669	-2%	
idzin	13400±2400	15102	13%	
idzein	241±5	262	9%	
nistein	108±10	93	14%	
/citein	211±5	190	10%	

## Table 4. Recovery results from a spiking experiment (UPLC method).

Original value (mg/kg)	Spiked level (mg/kg)	Determined value (mg/kg)	Recovery (%)	
317	499	827	102%	
159	499	634	95%	
165	499	665	100%	

## ANALYSIS OF ISOFLAVONES IN DIETARY SUPPLEMENTS

The isoflavone contents in three isoflavone dietary supplement samples were measured by this 8-minute UPLC-UV method. The sample

forms included tablets, capsules, and powder. The USP calibration and quantitation protocols were followed in the data processing The conversion factors for the acetyl and malonyl derivatives that were specified in the USP HPLC-UV

analysis. Table 6 shows the determined individual and total isoflavone contents, as well as the total

isoflavone contents on labels. For an

easy comparison, the label values were converted to the same concentration unit (mg/kg). Two of the three samples (C and D) showed

good agreement between the determined values and their label values, while one sample (B) had much ess measured total isoflavone

content than its label value. The reason for this low total isoflavone content in sample B compared to its

The isoflavone content in three isoflavone dietary supplement samples were measured by this 8-minute UPLC-UV method. The sample forms included tablets, capsules, and powder. The USP calibration and quantitation protocols <sup>(1)</sup> were followed in the data processing. The

	В		С		D	
	Mean	%RSD	Mean	%RSD	Mean	%RSD
	3095	1.0	12837	0.2	310	4.3
	617	2.0	4792	0.2	0	
	271	3.2	1766	0.6	664	1.1
idzin	97	7.8	655	2.0	0	
lzin	1634	1.7	7391	0.3	0	
itin	390	1.5	2202	0.04	0	
nistin	0		0		221	1.2
	6418	0.8	407	0.7	0	
	0		179	1.4	0	
istin	162	2.1	1042	0.4	0	
	0		0		0	
vones	12,684		31,270		1,195	
n value	>31	,250	>25,	,000	867	-2,600

conversion factors for the acetyl and malonyl derivatives that were specified in the USP HPLC-UV method were used in this analysis. Table 5 shows the determined individual and total isoflavone contents, as well as the total isoflavone contents on labels. For an easy comparison, the label values were converted to the same concentration unit (mg/kg). Two of the three samples (C and D) showed good agreement between the determined values and their label values, while one sample (B) had much less measured total isoflavone content than its label value. The reason for this low total isoflavone content in sample B compared to its label claim is unknown.




Figure 3. (A) Comparison of the UV (top) and the SIR (bottom) peaks for the daidzein in NIST 3238 SRM. A small peak fronting is observed in the UV, but not in the SIR peaks. (B) Comparison of the mass spectra (background subtracted) from the MS scan channel at daidzein peak's apex (top) and the peak front (bottom) for NIST 3238 SRM. The mass spectrum at peak front is different from the spectrum at peak apex. This confirms that the peak front is caused by a compound other than the daidzein.

# **CONCLUSION**

The USP isoflavone method was successfully transferred to UHPLC (on ACQUITY Arc system) and to UPLC (on ACQUITY H-Class system) svstems

The total run time is shortened from 74 min (HPLC) to 18 min (UHPLC) and to 8.1 min (UPLC)

Analysis results for the NIST reference material showed good agreement between the reference values and the experiment results

Mass detection is extremely useful in the analysis of complex samples. Its highly selective detection can reduce interferences, and the mass spectrum data is invaluable in troubleshooting

Dietary supplement samples have been analyzed by the UHPLC and the UPLC methods. One sample was found that its total isoflavone content was much lower than its label claim

## References

- 1. USP Monograph. Powdered Soy Isoflavones Extract, USP39–NF34 S1 [6841], The United States Pharmacopeial Convention.
- 2. L. X. Zhang, C. Q. Burdette, M. M. Phillips, C. A. Rimmer, and R. K Marcus. Determination of isofalvone content in SRM 3238 using liquid chromatographyparticle beam/electron ionization mass spectrometry, J AOA Int. 2015; 98(6): 1483

Please scan the QR code to download the relevant application notes

