

# Gel Permeation Chromatography of Phthalate in Walnut Oil Using Agilent EcoSpheres

## Application Note

Food Testing and Agriculture

### Author

Graham Cleaver  
Agilent Technologies, Inc.

### Introduction

Walnut oil from the common walnut, *Juglans nigra* L., Juglandaceae, is high-quality oil known for its delicate, nutty flavor and therefore used in cooking. The oil is produced mainly in France and is an expensive commodity that can become rancid if incorrectly stored. As part of a production process, walnut oil was to be packaged in plastic bottles manufactured from polyvinyl chloride (PVC). A stiff plastic, PVC is made flexible by the addition of plasticizers, typically phthalates, which help the polymer chains slide past each other. However, phthalates have been implicated as a potential health risk, raising concerns about phthalates leaching into foodstuffs.

This Application Note shows the separation of a phthalate additive from a sample of walnut oil as part of a quality control procedure, to isolate and measure the quantity and type of phthalate. To perform the separation, a 450 x 25 mm glass column was packed with Agilent EcoSpheres swollen in tetrahydrofuran.

EcoSpheres beads are designed for the economic, low pressure separations of pesticides, chlorinated hydrocarbons and other small molecules from high molecular weight organic matrices.




**Agilent Technologies**

## Method and Results

### Conditions

Column	EcoSpheres glass column, 450 x 25 mm (EcoSpheres loose media (100 g) p/n PL1460-4M03; complete glass column p/n PL1310-0054)
Sample	Walnut oil spiked with phthalate
System	Agilent 1260 Infinity Isocratic Pump and Manual Injector
Eluent	Tetrahydrofuran
Flow Rate	5 mL/min
Detector	Agilent Infinity Variable Wavelength Detector VL, 254 nm
Data collection/analysis	Cirrus GPC software and PL DataStream data capture unit



*Figure 1. Overlaid chromatogram of samples of spiked and un-spiked walnut oil produced by gel permeation chromatography with EcoSpheres microporous media. Clearly the phthalate is well resolved from the oil, allowing the small molecule to be isolated for further testing.*

## Conclusions

A GPC set up with EcoSpheres media is an easy and simple method for the analysis of food samples such as walnut oil. In this example, a phthalate was detected, suggesting such a system would be useful in the quality control of walnut oil packed in PVC containers.

For faster analysis times, the large particle size of EcoSpheres beads allows a 25 mm id glass column to be used at a flow rate of up to 10 mL/min, up to a maximum of 300 psi (21 bar). For maximum flexibility, EcoSpheres can be used with any liquid chromatography system capable of isocratic flow at the required flow rates. In addition, the pore size of EcoSpheres has been selected to maximize the separation of low molecular weight species such as phthalate from higher molecular weight interferences, which are excluded from the pores of the material.

## For More Information

These data represent typical results. For more information on our products and services, visit our Web site at [www.agilent.com/chem](http://www.agilent.com/chem).

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015  
Printed in the USA  
April 30, 2015  
5990-7623EN



**Agilent Technologies**