Waters

A COMPARISON OF 'H-NMR AND LC/MS(TOF) FOR A METABONOMICS EVALUATION OF RAT URINE FROM A TOXICOLOGICAL STUDY

Chris L. Stump^P, Maria Anthony², Robert S. Plumb¹, John N. Haselden², Jennifer H. Granger¹, Brian Sweatman², Jose Castro-Perez², and Hilary Major² ¹Waters Corporation, Milford, MA; ²GlaxoSmithKline, Ware, UK; ³Waters Corporation, Manchester UK

Presented at ASMS, Montreal, Canada • June 8-12, 2003

Introduction

The science of "Metabonomics" focuses on elucidation of markers of disease state, toxicological insult, etc. in biological fluids. These markers are important in drug discovery as they can streamline compound selection. Most drug candidates fail in Phase I or II clinical trials due to poor efficacy or toxicity. By utilizing a Metabonomics approach, poor drug candidates can be identified and eliminated much earlier in the drug discovery process. Metabonomics as practiced typically utilizes rapid screening analytical techniques such as ¹H-NMR or more recently LC/MS in combination with multivariate statistical analysis methodologies. In this study, rat urine from a toxicological study was analyzed by ¹H-NMR and LC/MS.

Methods

Samples:

Male and female rats were dosed with either vehicle alone, or a novel candidate pharmaceutical compound at 2 or 18 mg/kg. Urine samples were collected 0-8 and 8-24 hours post dose over the course of 90 days. Table 1 illustrates the dosing schedule, where day 33 represents the 0-8 hour collection point and day 34 represents the 8-24 hour collection. Samples were analyzed following a 1:4 dilution with water.

Table 1. Rat Dosing Schedule.

Rat Number	Dose Level	Gender	Time Points
1-8	Control	Μ	D33,D34,D84,D85
9-16	Low Dose	Μ	D33,D34,D84,D85
17-24	High Dose	Μ	D33,D34,D84,D85
25-32	Control	F	D33,D34,D84,D85
33-40	Low Dose	F	D33,D34,D84,D85
41-48	High Dose	F	D33,D34,D84,D85

LC Conditions:

LC System	Waters 2795 Separations Module
Column	Waters Symmetry $^{\scriptscriptstyle 8}$ C_{18} 3.5 μm
	2.1 x 100 mm
Flow Rate	600 µL/min. split to 150
	µL/min. into MS
Injection Volume	10 µL
Gradient	Linear, 0-20% B over 0.5-4
	min., 20-95% B over 4-8 min.,
	hold at 95% B for 1 min. then
	return to 0% B at 9 min.
Where $A = 0.1\%$ Fc	ormic Acid in Water,

and B = 0.1% Formic Acid in Acetonitrile

MS Conditions:

MC C	λ			
INIS System	vvaters /vlicromass [®] Q-lot			
	micro [™] with a LockSpray [™]			
	interface			
Ion Mode	Positive or negative			
	electrospray			
Cone Voltage	30 V			
Capillary Voltage	3.2 kV			
Desolvation Temp.	250 °C			
Source Temp.	120 °C			
Collision Energy	10 eV			
Detection Mode	Full scan 50-1400 m/z,			
	centroid mode			
Dwell	0.1 sec.			
Collision Gas	Argon, 5 x 10 ^{.3} mbar			
Lock Mass	Leucine Enkephalin at a			
	concentration of 50 fmol/µL, a			
	flow rate of 20 µL/min. was			
	employed with a frequency of			
	10 seconds			

¹H-NMR:

Single pulse ¹H-NMR spectra of all urine samples were obtained at 700 MHz , TSP in D₂O was added to the samples for referencing purposes.

Waters

Data Analysis:

PosterREPRIN

Two independent multivariate statistical analyses were conducted on the respective spectroscopic data.



Results

Typical LC/MS and ¹H-NMR data that were used in the Metabonomics approach is shown in Figures 1 and 2. A 12-minute reversed phase total ion chromatogram from a control rat urine sample is shown in Figure 1, while Figure 2 shows a typical rat urine ¹H-NMR spectrum.



Figure 1. LC/MS total ion chromatogram of control rat urine.





Detecting Dose Group Variability:



Figure 3. LC/MS data from male rats on day 85 from control, low dose, and high dose. Metabolites of GSK Compound A were located in the loadings plot at m/z 304 and m/z 319; these were independently confirmed by dosing of the radio-labeled compound and LC/MS/MS.



Figure 4. Extracted ion chromatograms of drug metabolites of m/z 304 and 319 from male rat urine where the rat was dosed with 18 mg/kg of GSK compound A and the urine was collected on day 85.

PosterREPRINT



Figure 5. The detected drug metabolites of m/z 304 and 319 were excluded from the analysis and PCA was re-run. The control and low dosage groups still cluster together away from the high dose group indicating that the GSK compound A is producing an effect on the endogenous metabolite profile of the rats in the high dosage group.



Figure 6. A comparison of the LC/MS (top) and ¹H-NMR (bottom) Scores and Loadings plots for the female rats on day 85, 8-24 hours for control, low dose, and high dose. The scores plot of both analytical techniques show that the control and low dose samples cluster together away from the high dose group indicating that the endogenous metabolite profile of the high dose group has changed.



Figure 7. LC/MS and ¹H-NMR Scores plots of males rats from control, low dose, and high dose on day 85. Again, high dose group clusters away from control and low dose groups in the scores plot obtained from both analytical techniques.

Detecting Gender and Diurnal Effects:



Figure 8. The metabolic profile of male and female control rats are sufficiently different to allow for differentiation by LC/MS and PCA (left). Subtle metabolic differences are discerned in the analysis of male rat urine collected in the morning and evening of day 33 (right). Similar group clustering is observed using ¹H-NMR.

Waters

Waters

Discussion

The LC/MS and 'H-NMR data shown here illustrate that GSK compound A induces a change in the endogenous metabolism of the rats when dosing at 18 mg/kg while the low dose group (2 mg/kg) remains similar to the control group. Because LC/MS and ¹H-NMR have different sensitivities with regard to analyte detection, it is likely that the group clusterings found in the PCA scores plots are not always the result of detecting the same chemical moieties. An advantage of the 3D nature of the LC/MS (TOF) data over ¹H-NMR is that after metabolite removal (Figures 3-5), more spectral data is retained for subsequent multivariate analysis. A further advantage of LC/MS (TOF) over ¹H-NMR is the ability to provide structurally relevant information that can be readily searched in a database. The elemental composition and MS/MS fragment ion information provided by LC/MS is amenable to providing a short list of potential marker compounds. Hence, to conclude this work, an example of metabolite detection and elucidation is illustrated using negative ion ESI LC/MS (TOF) data from the dose group variability study from day 85, 8-24 hours.



Figure 9. Scores and Loadings Plots of Negative ESI data from dose variability study on day 85, 8-24 hours. The structure of the identified marker at m/z 212 was elucidated.

			I ana I -	oth							
			885	<u>a</u>							
Single I	Mass Analy	sis									
Toleran	ce = 20.0 PF	M /	DBE:	min =	1.0, max = 10.0						
Isotope	cluster para	meters	Sena	aratio	n = 1.0 Abundar	$nce = 1.0^{\circ}$	%				
Manaica	tonic Mace	Odd and	d Even	Flech	ron lone	1.0					
000 /											
206 form	ula(e) evalua	ated wit	h 4 res	ults w	vithin limits (up to	50 closes	tres	ults	for e	ach	ma
206 form	ula(e) evalua	ated wit	h 4 res	DBE	vithin limits (up to	50 closes	t res	ults H	for e	ach 0	ma
206 form Mass 212.0016	Calc. Mass 212.0018	mDa -0.2	h 4 res	DBE	Formula C8 H6 N O4 S	50 closes Score	t res	ults H 6	for e	o 4	ma
206 form Mass 212.0016	Cak. Mass 212.0018 212.0004	mDa -0.2 1.2	h 4 res	DBE 6.5 7.0	Formula C8 H6 N O4 S C6 H4 N4 O3 S	50 closes	t res	H 6 4	for e	0 4 3	5 1 1
206 form Mass 212.0016	Calc. Mass 212.0018 212.0004 211.9977	mDa -0.2 1.2 3.9	PPM -0.7 5.6 18.2	DBE 6.5 7.0 2.5	Formula C8 H6 N 04 5 C6 H4 N4 03 5 C3 H6 N3 06 5	50 closes	t res	H 6 4 6	for e	ach 0 4 3 6	ma:

Figure 10. Elemental composition (C₈H₆NO₄S) of indoxyl sulphate determined from accurate mass LC/MS (TOF). LC/MS/MS and ¹H-NMR confirms structure.



Indoxyl sulphate, m/z 212.0018

Conclusions

Drug metabolites of GSK Compound A could be detected from the LC/MS and PCA loadings plot (from the dose group variability study).

The LC/MS and ¹H-NMR data produced similar group clustering in the dose group variability study.

Endogenous metabolite differences resulting from gender and diurnal effects could be discerned by using LC/MS and multivariate analysis (similar results have been reported with 'H-NMR).

Identification of potentially important biomarkers has been streamlined by using LC/MS-TOF exact mass and an integrated software browser. Hence, the structural elucidation stage of metabolite marker identification can be reached rapidly, allowing more effort to be focused on this challenging area.

WATERS CORPORATION 34 Maple St. Milford, MA 01757 U.S.A. T: 508 478 2000 F: 508 872 1990 www.waters.com





Waters, Micromass, Q-Tof micro, LockSpray, MarkerLynx, ChromaLynx and Symmetry are trademarks of Waters Corporation.

All other trademarks are the property of their respective owners. ©2003 Waters Corporation Produced in the U.S.A. June 2003 720000630EN LW-PDF

