#### THE USE OF PARALLEL LC/MS/MS WITH AUTOMATED OPTIMIZATION TO INCREASE THROUGHPUT FOR THE QUANTIFICATION OF INCUBATED DRUG CANDIDATE SAMPLES

Tabisam Khan', Kendon Graham', Paul Courtney<sup>2</sup> and Richard Weaver<sup>2</sup> Waters Corporation, Manchester, UK. <sup>2</sup>AstraZeneca R&D Charnwood, Physical and Metabolic Sciences, Loughborough, UK **Presented at ASMS, Montréal, Canada, June 8-12, 2003** 

#### OVERVIEW

To demonstrate the use of the Waters® Micromass® Quattro micro™ MUX-QuanOptimize™ together with multi-pump control (four Waters 1525µ Binary HPLC Pumps) pumps) in the analysis of incubated samples, in order to increase throughput for quantitative ADME studies in drug discovery.

Here we present a multiplexed electrospray ionization source interfaced to a tandem quadrupole mass spectrometer. The system is set up in such a way as to automatically generate compound optimization and integration parameters using QuanOptimize.

This resulted in the generation of intrinsic clearance plots for 76 drug candidates.

#### INTRODUCTION

High throughput quantification is an essential aspect of the drug discovery process. The demands on instrumentation and the chemist are growing at a rapid pace. There is a need to perform fast quantitative measurements for early determination of ADME properties. This information is then used to determine which of the drug candidates will progress to the development stage. The fail fast strategy is essential to prevent unnecessary work on compounds that do not fit the criteria. In this poster we meet the challenge of fast and efficient turnaround of multiple samples.

The main bottleneck was the ability to produce samples quickly enough for analysis but with the advent of multiple robotic sample production the bottleneck moves from preparation to the analytical stage. Since there are still many compounds to characterize at this stage of the development process, the optimization and quantification can be a time-consuming process. In order to ascertain the correct MRM transitions, looking at the precursor ion (or an adduct) and then obtaining a suitable product ion with the correct cone voltage and collision energy, can take an experienced mass spectroscopist a long time. With QuanOptimize, an automatic method development and quantification tool, samples can be efficiently analyzed in two stages. Firstly, each compound can be optimized to generate a multiple reaction monitoring experiment, where the product and precursor (with cone voltage and collision energies) ions are defined, in a run time of approximately 1 minute. Secondly, using generic HPLC fast gradients a quantification method can be produced. The whole procedure for a single compound can take as little as 2 minutes to turn around. With the addition of multiple parallel sprayers, MUX-technology<sup>™</sup>, interface and sample grouping the throughput can be boosted by a factor of four. With the introduction of 4 LC systems, under full MassLynx<sup>™</sup> 4.0 Software control, parallel chromatography using a 4 injection module (Waters 2777 Sample Manager) with 4 matched LC columns linked directly to a MUX source, we can easily adapt to the vast number of samples being generated.

In this poster we analyzed 76 individual compounds and 5 time points for each compound (a total of 380 samples) in two hours. The compounds were incubated with rat hepatacytes with sampling time points at 8,20,41,58 and 73 minutes that were extracted in methanol.

#### **EXPERIMENTAL SETUP**

Figure 1 shows independent LC pumps controlled under MassLynx 4.0 Software delivering independent gradients to four separate analytical columns into a 4-way MUXtechnology-equipped tandem mass spectrometer. The individual pumps are controlled under MassLynx 4.0 Software with serial connections via IEEE cables.

Figures 2, 3 and 4 show the configuration of MassLynx 4.0 Software for MUX-QuanOptimize and Four LC Control. Figure 2 shows the method editor for QuanOptimize. Under MassLynx 4.0 Software control, the automatic generation of an MS method and the





·	1		
Losses	Inlet Met	hods	MUX
Optimisation	Acquisition	Tune Files	Adducts
Cone Voltage min max 10 50	step 6	Ionisation M	fode +ve C -ve C both
Collision Energy min max 10 35	step 4	MS Method MS MSMS	l Creation
Do Optimisations			
Per sam	ple group	O Before any an	alyses
Optimisation Peak D	etection Parameters	s	agment Size –
MS Quan Method:	ms	Browse	min: 50
MSMS Quan Methor	t poerce	Browse	1

Figure 2. QuanOptimize Method Editor.

subsequent analytical method can be run. Figure 3 then highlights the configuration of the multiplexed electrospray interface (MUX) with randomized sample collection from the sample manager. Finally in Figure 5, the number of probes available in the MUX settings is defined. Also, random bottle locations in the sample list and the number of LC inlets are defined.

QuanOptimise Metho File	d Editor QuanOpti	imiseMUX.QLI	<b>م</b> اء ا	×
Optimisation	Acquisition	Tune File: thods	s Adducts MUX	
I⊽ Enab I⊽ Interla	le MUX Mode aced Injection Mode	S	prays 1 2 4	

Figure 3. QuanOptimize Enabled MUX.

×

Options	×
Multi-probe	
Multi-probe Source Multi-probe Capability Accurate Mass 4 Prob <u>e</u> s	s
- Dual Source	
Inlets Allow Random Bottle Locations in Sample List	
Multi-Inlet Capability	
4 Number of Parallel Inlets	
☐ <u>S</u> ample Prep	
ОК	Cancel

Figure 4. Configuration of Multiple LC Pumps.

124 2 1	11	E Shertest	30me (550					
			Queu	e Is Emp	ty			
Instrument of	Spe	ctrum Chromato	gram Map E	dt. San	ples +			
ment terment an		FieNate	File Text	Batte	Sarghe Type	Sarph Grog	OweA	QueiRel
	55	shameta_055	Cass11 20 may	SK1-C1-43	Avaite	*		
	- 56	athannera 0%	Canit2 20 mm	581-018	Analyte	5		
Inlet Method	\$7	ashammena_057	Casil 41 min	581-01:30	Analyte	8		
	58	athaoneca_058	Cass1241 min	561-01:90	Analyte	3		
and the second	58	advantation,059	Caniff 41 mm	581-01-55	Analyte	8.		
000	80	athairmate,000	Can1241 en	561-21.20	Andite	4		
ME Method	- 67	athaneneta_001	Cand Sil nin	581-01-42	Analyte	8		
	62	athatemeta (002	Cast1218 sin	561-017	Anales	4		
	62	amazoraca_063	Cau(11 58 mm	561-01-67	Andito	8		
	64	anhaberence, 064	Cass1250 ver	181-0132	Analys	1		
MS Tune	48	athainence_055	Cault 73 min	581-0154	Analyte	1		
	66	mitainneca_055	Case1373 view	383-0119	Analyte	4		
1241	87	athgoreta_057	Cau(11 73 min	583-0179	Analyte			
	64	athermeta_008	Can1273 min	381-01-44	Anales	6		
t Shutdown or Startup	100	athanemia 009	STD1 - Cault	581-022	Dordet		100	
-	20	athainmeia, 070	ST01 - Caselto	3816214	Stended	8	108	
	21	athamete 071	STD1 - Cault	581-0228	Standard	K	100	
and the second second	77	anhairemin, 072	STD9 - Caluft2	561-0238	Standard	4	105	
During the second second	73	athairmeia,073	Card 3 8 mm	581-01.56	Andite	*		
	24	athermore, 0%	Cault Ener	581-01-21	Analys			
	- 25	athamera 075	Cau158 mm	361-01.01	Andre	0		
Stirfue)	25	antaminia,0%	Castliller	561-01-46	Analite			
	31	amainmenta_077	Ear(13:20 mm	581-01.68	Analyte	14		
	28	athamete 070	Cau14 20 mm	581-01.22	Anales	*		
<b>E</b> 23	79	athaometa (079	Caroth 20 run	361-0133	Analys	c		
Options	00	atheoneca,000	Carolis 20 mm	5810158	Andre	*		
	44	athemere Off	Carril 241 mm	5814180	Andre	14		

Figure 5. QuanOptimize Analytical Sample List for Incubated Samples.

#### **EXPERIMENTAL CONDITIONS**

#### **Sample Preparation**

Standards for automated tuning by QuanOptimize were provided by AstraZeneca. These consisted of 76 compounds grouped into fours, diluted to approximately 5 mgm/mL in methanol, resulting in 19 different sample groups.

AstraZeneca provided samples from hepacyte incubations at 5 different time points, in sets of six per sample group. The samples were incubated at time intervals of 8,20,41,58 and 73 minutes and when the samples were taken they were quenched with methanol (2:1) to stop further metabolism.

#### LC CONDITIONS

Column:	Four Waters XTerra® MS C <sub>18</sub> , 3.5 µm, 3 x 20 mm		
Mobile Phase A:	(0.1% formic acid)	acetonitrile,	
Mobile Phase B:	(0.1% formic acid)	Water	
Gradient:	A%=100,	time=0 min	
	A%=0,	time=0.3 min	
	A%=0,	time=1.0 min	
	A%=100,	time=1.2 min	
Flow Rate:	1.5 mL/min, post- generate 60 µL/m	column split to in into each spray	
Injection Volume:	5 µL		
Solvent Delivery:	Four Waters 1525 gradient pumps	ōμ binary	
Sample Manager:	Waters 2777, four injection valves diluter option		

#### MS INSTRUMENT SETTINGS

MS	Quattro micro API
Polarity	ES+
Capillary (kV)	3.50
Cone (V)	35.00
Source Temperature	140 °C
Desolvation Temperature	400 °C
Cone Gas Flow	62 L/Hr
Desolvation Gas Flow	671 L/Hr
Peak width at half height	0.7 Da
LM 1 Resolution	14.0
HM 1 Resolution	14.0
Peak width at half height	0.7 Da
LM 2 Resolution	14.0
HM 2 Resolution	14.0
Gas Cell Pirani Pressure	5.74 <sup>e-3</sup> mbar

#### RESULTS

Sample group A			
MRM Transition	Dwel (secs)	Cone Volt/VI	Col Energy (eV)
272.2 > 147.1	0.03	35.0	30.0
285.1 > 193.0	0.03	35.0	30.0
415.2 ~ 179.0	0.03	25.0	30.0
415.2 > 170.0	0.03	35.0	27.0
500.0 > 100.0	0.03	35.0	28.0 (dummy injection
Sample aroup B			
MRM Transition	Dwell(secs)	Cone Volt (V)	Col.Energy (eV)
3911 > 914	0.03	35.0	30.0
422.9 > 213.9	0.03	35.0	23.0
446.9 > 92.4	0.02	25.0	22.0
482.1 > 273.1	0.03	35.0	23.0
Sample group C	D 11/ )	C 1/ 1/20	0.15 (10
MKM Transition	Dwell(secs)	Cone volt(v)	Col.Energy (eV)
425.0 > 125.2	0.03	35.0	20.0
462.9 > 108.3	0.03	35.0	23.0
333.9 > 255.9	0.03	35.0	19.0
446.9 > 92.3	0.03	35.0	20.0
Sample aroun D			
MRM Transition	Dwell(secs)	Cone Volt(V)	Col.Energy (eV)
437.0 > 169.9	0.03	35.0	10.0
469.0 > 92.4	0.03	35.0	10.0
449.0 > 310.8	0.03	35.0	10.0
334.0 > 255.8	0.03	35.0	10.0
Sample group E			
MRM Transition	Dwell(secs)	Cone Volt(V)	Col.Energy (eV)
366.1 > 121.3	0.03	35.0	35.0
392.2 > 204.2	0.03	35.0	30.0
376.2 > 323.1	0.0	35.0	32.0
443.1 > 218.1	0.03	35.0	35.0
Sample aroun F			
MRM Transition	Dwall(cocs)	Cone Velt0/	Col Enorgy (o)/)
202.2 - 240.0	Dweil(secs)	25.0	20.0
393.2 > 340.0	0.03	33.0	30.0
313.1 > 280.0	0.03	35.0	30.0
319.1 > 85.4	0.03	35.0	30.0
2/4.1 > 102.4	0.03	35.0	28.0
Sample group G			
MRM Transition	Dwell(secs)	Cone Volt(V)	Col.Energy (eV)
479 2 > 421 9	0.03	35.0	30.0
641.2 > 570.0	0.03	35.0	30.0
552 3 > 72 7	0.03	35.0	30.0
587.3 > 138.3	0.03	35.0	30.0
Sample group H			
MK M Transition	Dwell(secs)	Cone Volt(V)	Col.Energy (eV)
588.2 > 138.3	0.03	35.0	30.0
580.4 > 217.7	0.03	35.0	30.0
526.2 > 247.0	0.03	35.0	30.0
486.3 > 235.1	0.03	35.0	28.0
Sample aroun I			
MRM Transition	Dwell(secs)	Cone Vol t(V)	Col.Energy (eV)
589.3 > 140.3	0.03	35.0	30.0
524.2 ~ 157.2	0.03	25.0	20.0
J24.J ≥ 13/.Z	0.03	35.0	20.0
641.3 > 5/0.2 538.3 > 157.3	0.03	35.0	30.0
	0.00	05.0	2
Sample group J			
MRM Transition	Dwell(secs)	Cone Volt(V)	Col.Energy (eV)
371.3 > 135.2	0.03	35.0	37.0
402.2 > 153.0	0.03	35.0	37.0
377.3 > 112.3	0.03	35.0	30.0
277.2 > 112.2	0.03	35.0	28.0

Table 1. Automatic generation of MRM transitions by QuanOptimize.

#### 1: MS Data Acquisition Settings

Table 1 shows the fully optimized product and precursor ion with the correct cone voltage and collision energy. These transitions are then automatically placed into a MS Method editor and the quantification assay is then acquired. The results from the quantification are discussed in the result section.

The correlation between the MRM transition obtained by MUX-QuanOptimize and by single spray was excellent. Each compound was correctly assigned the precursor and product ion as well as the correct cone voltage and collision energy. The dwell time was automatically set from the software to take into account the amount of MRM transitions per sample group and provide sufficient data points (30 points per dalton) across a chromatographic peak to allow accurate peak integration.

#### 2: Analysis and Intrinsic clearance plots

The samples were run in accordance to the group numbers with each group containing four compounds. The final group consisted of a dummy run in order to make the sample list a divisible factor of four as the running of the MUX source with QuanOptimize makes this necessary.

The results for the optimization can be seen in the Results 1 section. Where sample was a low volume (insufficient for injections) the MS transitions obtained by AstraZeneca were used as the sample well plate had already been run at AstraZeneca and repeated at Waters Corporation. This allowed the running of all the samples in the 96-well plate. The samples were then ready for the analytical stage where the retention time would be calculated and the chromatograms integrated automatically using QuanOptimize.

The time taken for running all the samples was approximately 2 hours using the Waters Micromass Quattro micro MUX-QuanOptimize as opposed to a run time of 6.5 hours if the samples were run using the conventional single spray source.

A total of 76 compounds were analyzed and for each compound five time points were analyzed (380 samples analyzed in total with 1 dummy compound and 5 dummy time points to allow a factor of four for the software to recognize MUX is being used). Each compound contained five time points at 8, 20 41,58 and 73 minutes and the peak areas of each were plotted against time as shown in the analytical sample list (Figure 5).

The following clearance plot diagrams give an indication of the incubated samples clearance with time:



Figure 6. Example of the clearance of drug candidate A with time.



Figure 7. Example of the clearance of drug candidate B with time.

The clearance plots above (Figures 6 and 7) are some of the examples run and show a good correlation to the data obtained with a single sprayer instrument. The samples show the reducing peak areas with the set time points from the incubated samples. For the plots individual peak areas are taken from the MUX source. This data is then plotted onto logarithmic scale [log peak area] vs [time] and gave the following results (Figures 8 and 9).

All the data analyzed in this way resulted in R2 values of greater than 0.9, which shows excellent correlation between MUX and single spray.

# **Poster**REPRINT

# Waters



Figure 8. Intrinsic Clearance plot of drug candidate A.



Figure 9. Clearance plot of drug candidate B.

The combined data for all 76 compounds into 19 sample groups and the subsequent sampling time points was plotted for the MUX and directly compare to the results generated by a single sprayer, the graph below shows the excellent correlation (Figure 10).

The data in Figure 10 shows how the single sprayer compared with the MUX source on the Quattro micro. The correlation between the two was R2>0.8. This result shows that the single spray and MUX source are capable of close correlation for the analysis of compounds with a range of intrinsic clearance values. This means that the compounds that pass the initial stage of analysis can be moved onto further pharmacokinetic studies. It also highlights that the conclusion drawn at this stage of the drug discovery process can be arrived at much quicker with a MUX source than the single sprayer source whilst maintaining the integrity of the data.



Figure 10. Intrinsic Clearance Data: Single spray vs MUX.

#### CONCLUSIONS

The system that has been employed in this study was designed to meet the challenge of high throughput analysis, to come to a conclusion on whether or not to carry the drug candidates through to the next stages of the drug discovery process. Using the 4 LC pumps, the Waters 1525µ Binary HPLC Pumps under full MassLynx Software control, and the MUX source the data was reliable and accurate. The comparison with the data was excellent with good correlation, R2>0.8

The Quattro micro MUX-QuanOptimize demonstrated that the clearance plots could be generated in approximately two hours with fast and accurate information for the tuning conditions. The speed which samples for each group are analyzed allows fast interpretation of the clearance plots using QuanLynx<sup>™</sup> software.

The use of multi-pump control with the four Waters 1525µ Binary HPLC Pumps overcomes the issues associated with having potential blockage of columns. If one column does have a problem the data on the remaining columns is not affected. With the Waters XTerra MS C<sub>18</sub>, 3.5 µm, 3 x 20 mm columns a very fast generic gradient can be applied whilst being robust and reproducible enough to analyze many hundreds of samples.

There is also the possibility that multiple mobile phases can be used in order to set up different gradients on each pump, hence extending the scope of MUX by increasing the number of different compounds and associated chromatography that can be analyzed.

To whom all correspondence should be addressed: e-mail: Kendon\_Graham@waters.com

WATERS CORPORATION 34 Maple St. Milford, MA 01757 U.S.A. T: 508 478 2000 F: 508 872 1990 www.waters.com







Waters, Micromass, Quattro micro, MUX-QuanOptimize, MUX-technology, MassLynx, XTerra and QuanLynx are trademarks of Waters Corporation. All other trademarks are the property of their respective owners. ©2003 Waters Corporation Produced in the U.S.A. June 2003 720000645EN LW-PDF

