

Henry Y. Shion^{*1}; Stacy S. Powers²; ¹Waters Corporation, Beverly, MA²Lab Connections, a Mocon Company, Northborough, MA

Presented at ASMS, Montréal, Canada, 8th-12th June 2003

OVERVIEW

Poly-(DL-Lactide-co-glycolide) (PLG) polymers are biocompatible and biodegradable polyesters used in drug delivery. Molecular weight data for these polymers is of importance as it relates to performance characteristics.

PLG polymers are usually highly polydisperse ($P_d > 1.5$), thus presents a challenge for MALDI analysis.

GPC fractionation and MALDI-MS are coupled for the study of dispersed PLG.

INTRODUCTION

PLG polymers often have large polydispersities ($P_d > 1.5$).

Although (MALDI-MS) has been used extensively to provide molecular weight, structural and compositional information of synthetic polymers, it often fails to provide correct molecular weight values for polydisperse ($P_d > 1.2$) polymers.

The combination of GPC and MALDI-MS is used as a method for the characterization of polydisperse PLG by fractionation of the polymer into narrower molecular weight fractions before analysis by MALDI.

The polydispersity of the selected PLG sample is given as 5,000 to 15,000 Daltons by the manufacturer (Aldrich, Milwaukee, WI).

The GPC analytes (fractions) are deposited onto a matrix (Dithranol) pre-coated foil for subsequent MALDI analysis utilizing a LC-Transform[®] 600 instrument.




Figure 1. Poly-(DL-Lactide-co-glycolide) (PLG).

PLG is a copolymer of lactide and glycolide. As a result even low molecular weight oligomers exist as a complex mixture of molecular weights.

Lactide Glycolide Copolymer	MW: acid	MW: mer
Lactic acid: C ₃ O ₃ H ₆	90	72
Glycolic acid: C ₂ O ₃ H ₄	76	58

Table 1: PLG Copolymer repeating units information.




Figure 2. MALDI spectrum of PLG prior to GPC.

METHODS**Gel Permeation Chromatography (GPC)**

- Column: Waters Styragel[®] HR2 4.6 x 300 mm at 25°C
- HPLC Pump: Waters 616 pump with 600S controller
- Detector: Waters 2410 Refractive Index detector and 2996 Photodiode Array (PDA) detector
- Flow rate: Isocratic gradient of tetrahydrofuran (THF) at 0.3 mL/min for 60 min
- Injection volume: 20 mL

LC-Transform Sample Collection Conditions

Sample concentration:	10 mg/mL
Sheath gas temperature:	125°C
Sheath gas flow:	30
Nebulizer gas flow:	25
Foil:	Dithranol 1-Time®/with sodium salt
Moving rate of collection	
Stage:	5 mm/min

MALDI-MS

Instrument:	Waters MALDI LR
Source Voltage:	15 KV
Delay time:	500 ns
Laser:	337 nm Nitrogen Laser from LSI (Franklin, MA)
Laser Rate:	10 Hz
Mode:	Positive ions in Reflectron mode
Detector:	MCP

RESULTS**GPC PDA spectrum**

The GPC UV spectrum (Figure 3) shows that the polymer sample peak is between 10 min and 20 min.

The small peak close to the end of the run was not identified; no useful MALDI mass spectrum could be obtained from this fraction.




Figure 3. GPC PDA spectrum of the selected PLG co-polymer at UV wavelength of 254 nm.

MALDI-MS Spectra

Figure 4 shows a nice polymer distribution centered at about $m/z=2,000$, which indicates the GPC has separated the polydispersed PLG into narrow fractions. This is an example of one fraction.

Figure 5 is an expanded spectrum from Figure 4, which shows that the mass difference between the spectrum peak clusters is 14 amu, which is the mass difference between the monomers of Lactide $[\text{OCH}(\text{CH}_3)\text{C}(\text{O})]$ and Glycolide $[\text{OCH}_2\text{C}(\text{O})]$.





Figure 4. MALDI-MS spectrum from one of the GPC fraction spots on the MALDI target plate. The spectrum is a summation of 100 laser shots.



Sales Offices

AUSTRIA AND EXPORT (CENTRAL EUROPE, CIS, MIDDLE EAST,
INDIA AND INDIA SUBCONTINENT)

43 1 8771807

AUSTRALIA 61 2 9933 1777

BELGIUM AND LUXEMBOURG
02 726 1000

BRAZIL 55 11 5543 7788

CANADA 800 252 4752 X2205

CIS 7 095 931 9193

CZECH REPUBLIC 42 02 617 11384

DENMARK 46 59 8080

FINLAND 358 09 506 4140

FRANCE 1 3048 7200

GERMANY 06196 40 06 00

HONG KONG 852 29 64 1800

HUNGARY 36 1 350 5086

INDIA 91 80 837 1900

IRELAND 1 4481500

ITALY 02 274211

JAPAN 81 3 3471 7191

KOREA 82 2 3284 1300

MEXICO 52 55 5524 7636

THE NETHERLANDS 76 508 7200

NORWAY 6 384 6050

PEOPLES REPUBLIC OF CHINA

86 10 8451 8918

POLAND 48 22 833 4400

PUERTO RICO 787 747 8445

SINGAPORE 65 6278 7997

SPAIN 93 600 9300

SWEDEN 555 11 500

SWITZERLAND 62 889 2030

TAIWAN 886 2 2543 1898

UK 44 208 238 6100

U.S.A. AND ALL OTHER COUNTRIES:

WATERS CORPORATION

34 Maple St.

Milford, MA 01757 U.S.A.

T: 508 478 2000

F: 508 872 1990

www.waters.com

WATERS CORPORATION
34 Maple St.
Milford, MA 01757 U.S.A.
T: 508 478 2000
F: 508 872 1990
www.waters.com

Waters
RIGHT ON TIME.

Qpt-In! My Profile
www.waters.com

Waters and Styragel are registered trademarks of Waters Corporation.
All other trademarks are the property of their respective owners.
©2003 Waters Corporation Produced in the U.S.A. June 2003 720000656EN SD-PDF



001