

Michael Jackson, Michael P. Balogh, Warren Potts, Waters Corporation, Milford, MA 01757, USA

Presented at ASMS, Montréal, Canada, June 8-12, 2003

INTRODUCTION

- Full characterization of the physicochemical properties of new chemical entities (NCE's) typically takes place in pharmaceutical development laboratories. However, in an effort to minimize drug development costs, pharmaceutical companies are now pushing physicochemical screening back into the drug discovery laboratories so that unsuitable compounds fail earlier in the process.
- By generating retention time and quantitative data, LC/UV/MS provides a mechanism by which indicators of compounds' physicochemical properties—such as solubility, hydrophobicity, membrane permeability, stability, and protein binding—can be modeled and screened.
- The large number of compounds that require screening can be analyzed by high throughput parallel LC/UV/MS methods. This technology allows the analysis of more samples in the same time, or the analysis of multiple physicochemical properties at the same time. While many compounds can be tested, large volumes of data from these parallel analyses must still be processed and reviewed.
- The new software application presented here, ProfileLynx™ Application Manager, handles sample set management, data processing and results browsing, as well as provides exporting capabilities and interactivity with other software packages to simplify physicochemical property screening.
- The multiple streams of acquisition data from parallel LC/UV/MS systems are integrated into an easily managed and interpreted format by ProfileLynx.


PHYSICOCHEMICAL PROPERTY PROFILING

Compounds were screened for their physicochemical properties using three different indicators:

- Chromatographic hydrophobicity index using C₁₈ column retention (CHI) at various pH's
 - The CHI relates directly to the log D of a compound, which indicates the compound's solubility and permeability
- Chromatographic hydrophobicity index using immobilized artificial membrane (IAM) columns (CHI_{IAM})
 - The CHI_{IAM} is measured on IAM columns under gradient conditions and relates to the membrane permeability of a compound
- Solubility

EXPERIMENTS

Single Chromatographic System


*Waters ECP™ System for Early Candidate Profiling***Experiment 1: CHI Screening**

- Standards with known CHI values were run to determine the correlation between retention time and CHI
- Pharmaceutical compounds were then analyzed and their CHI values were calculated from the compounds' chromatographic retention time at pH 7.4

Experiment 2: Simulated Solubility Screening

- A wide set of pharmaceutical compounds were used to simulate a solubility study
- The analysis used a plate of standards of known concentration and then a plate of samples with varying concentrations
- Calibration curves were developed using the standards, and compounds were quantified to determine their "solubility"

Parallel Chromatographic System



Waters ECP™ System for Early Candidate Profiling

Experiment 1: CHI and CHI_{IAM} Screening

- The parallel system used four independent binary gradient pumps for simultaneous compound analysis under four different separation conditions
- Standards with known CHI and CHI_{IAM} values were run to determine the correlation between retention time and CHI or CHI_{IAM}
- Pharmaceutical compounds were then analyzed and their CHI or CHI_{IAM} values were calculated from the chromatographic retention time
- CHI values were calculated based on compound retention on C₁₈ columns at three different pH conditions (pH 2.0, pH 7.4, and pH 10.5)
- CHI_{IAM} values were calculated based on compound retention under simulated physiological conditions (pH 7.4)

Experiment 2: Simulated Solubility Screening

- The parallel system used four independent binary gradient pumps for the simultaneous analysis of four different compounds under the same conditions, resulting in a four-fold increase in sample throughput over the single chromatographic system.
- Calibration curves were developed using the standards, and compounds were quantified to determine their "solubility".

METHODS**Columns**

- CHI
 - Waters XTerra® MS C₁₈ 3.5mm 3.9x50mm
- CHI_{IAM}
 - Regis Technologies IAM Fast-Screen Mini Column 3x10mm with Ester IAM.PC.C10/C3 packing
- Simulated Solubility
 - Waters XTerra® MS C₁₈ 3.5mm 3.9x50mm

Mobile Phases

- CHI & CHI_{IAM}
 - A: CHI pH 2.0, 0.01% v/v phosphoric acid
 - A: CHI pH 7.4, 10 mM ammonium acetate
 - A: CHI pH 10.5, 10 mM ammonium acetate
 - A: CHI_{IAM} pH 7.4, 10 mM ammonium acetate
 - B: 100% ACN
- Simulated Solubility
 - A: H₂O
 - B: 100% ACN

HPLC GradientsCHI and CHI_{IAM}

Time	Flow	%A	%B	Curve
0	0.5 mL/min	100%	0%	-
10	0.5 mL/min	0%	100%	6
12	0.5 mL/min	0%	100%	6
14	0.5 mL/min	100%	0%	1

Injection Volume

- 2.0 µL

MS Conditions


- ESI+ and ESI-
- 100 to 600 amu
- Scan time 0.3 sec
- Interscan delay 0.1 sec
- Cone voltage 30V
- Capillary voltage 3.0 kV
- Source temp 120 °C
- Desolvation temp 300 °C

Simulated Stability

Time	Flow	%A	%B	Curve
0	0.5 mL/min	100%	0%	-
4	0.5 mL/min	0%	100%	6
6	0.5 mL/min	100%	0%	6

PHYSICOCHEMICAL PROPERTY PROFILING SOFTWARE

Step 1: Define the physicochemical property experiment



MassLynx - Default: ECP: Protein Data Set

ProfileLynx

Experiment Editor

Calibration Method

Process Calibration

View Calibration

Profile Method

Process Samples


View Results

File Name	File Type	File Size	Sample Name	Sample Time	Mass	Compound	Peak Method	Peak Value	Peak Name	Peak
13	Sample1	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
14	Sample2	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
15	Sample3	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
16	Sample4	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
17	Sample5	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
18	Sample6	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
19	Sample7	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
20	Sample8	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
21	Sample9	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
22	Sample10	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
23	Sample11	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
24	Sample12	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
25	Sample13	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
26	Sample14	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
27	Sample15	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
28	Sample16	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
29	Sample17	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
30	Sample18	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
31	Sample19	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
32	Sample20	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
33	Sample21	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
34	Sample22	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
35	Sample23	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
36	Sample24	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
37	Sample25	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
38	Sample26	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
39	Sample27	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
40	Sample28	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
41	Sample29	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
42	Sample30	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
43	Sample31	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
44	Sample32	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
45	Sample33	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
46	Sample34	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
47	Sample35	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
48	Sample36	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
49	Sample37	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
50	Sample38	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
51	Sample39	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
52	Sample40	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
53	Sample41	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
54	Sample42	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
55	Sample43	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
56	Sample44	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
57	Sample45	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
58	Sample46	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
59	Sample47	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
60	Sample48	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
61	Sample49	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
62	Sample50	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
63	Sample51	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
64	Sample52	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
65	Sample53	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
66	Sample54	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
67	Sample55	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
68	Sample56	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
69	Sample57	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
70	Sample58	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
71	Sample59	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
72	Sample60	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
73	Sample61	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
74	Sample62	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
75	Sample63	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
76	Sample64	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
77	Sample65	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
78	Sample66	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
79	Sample67	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
80	Sample68	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
81	Sample69	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
82	Sample70	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
83	Sample71	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
84	Sample72	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
85	Sample73	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
86	Sample74	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
87	Sample75	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
88	Sample76	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
89	Sample77	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
90	Sample78	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
91	Sample79	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
92	Sample80	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
93	Sample81	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
94	Sample82	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
95	Sample83	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
96	Sample84	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
97	Sample85	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
98	Sample86	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
99	Sample87	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
100	Sample88	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
101	Sample89	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
102	Sample90	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
103	Sample91	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
104	Sample92	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
105	Sample93	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
106	Sample94	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
107	Sample95	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
108	Sample96	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
109	Sample97	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
110	Sample98	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
111	Sample99	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		
112	Sample100	ECP: Run Run: ECP: Run	ECP: Run Run: ECP: Run	24	140	485	Unassigned	Not Calibrated		

Sample List

Using the MassLynx sample list, the physicochemical experiment parameters can be associated with appropriate samples. From the Sample List, all of the ProfileLynx Application Manager tools can be accessed.

Step 2: Define and set experimental criteria for physicochemical property calibration




EPCalib8.chp - EPCalibEditor

Compound List	Property	Value
Compound 1	Calibration Type	Protein Binding
Compound 2	Compound Values	
Compound 3	Compound Name	theophylline
Compound 4	Compound Value	18-400
Compound 5	Peak Location	
Compound 6	Function Number	Two
Compound 7	Trace	188
Compound 8	Chromatogram Mass Window (Da)	1,000
Compound 9	Start Retention Time	1,000
Compound 10	End Retention Time	9,000
Compound 11	Locate Peak Selection	Largest
Compound 12	Smoothing Enabled	<input checked="" type="checkbox"/> YES
Compound 13	Smooth Parameters	
Compound 14	Smoothing Method	Mean
Compound 15	Smoothing Iterations	2
Compound 16	Smoothing Width	3
Compound 17	Apex Track Enabled	<input checked="" type="checkbox"/> YES
Compound 18	Peak Detect Parameters	
Compound 19	Peak-to-Peak Baseline Noise	<input checked="" type="checkbox"/> 10.00
Compound 20	Peak Width at 5% Height	<input checked="" type="checkbox"/> 30.00
Compound 21	Baseline Start Threshold %	0.00
Compound 22	Baseline End Threshold %	0.50
Compound 23	Detect Shoulder Peaks	<input checked="" type="checkbox"/> NO
Compound 24	View Standard Peak Detection Parameters	
Compound 25	View Threshold Parameters	
Compound 26	Threshold Relative Height	<input checked="" type="checkbox"/> 1.50
Compound 27	Threshold Absolute Height	0.00
Compound 28	Threshold Relative Area	<input checked="" type="checkbox"/> 2.00
Compound 29	Threshold Absolute Area	0.00
Compound 30	Propagate Integration Parameters?	<input checked="" type="checkbox"/> YES

Calibration Editor

Developing the relationship between retention times and physicochemical property indicators is easily enabled through the Calibration Editor in ProfileLynx. In the Calibration Editor the user simply selects the type of retention time calibration to be used (CHI, CHI_{IAM}, IAM, or Protein Binding) and then provides a customizable name for the specific experiment. The user then sets criteria such as the equation to use for calculating the physicochemical property indicator, integration parameters, and retention time ranges for the retention time calibration.


Step 3: Review physicochemical property calibration



Data Review Tool

Retention time calibration review is facilitated through the Calibration Review Tool. Chromatograms can be manually integrated, peaks can be reassigned, and points can be ignored. The calibration is then stored and can be used during compound analysis.


Step 4: Select the type of experiment to be run and set pass/fail criteria




Experiment Editor

The type of physicochemical property screening experiment to be run is selected in the Experiment Editor. Options include CHI, IAM, solubility, permeability, stability, protein binding, partitioning, and metabolic inhibition.

Step 5: Review physicochemical property analysis results by compound or experiment




Data Review By Compound



Data Review By Experiment

Processed analysis results can be viewed by experiment or by compound using the data review tool. The sample table shows the calculated physicochemical property indicator values for the samples, as well as markings to indicate whether the sample passed or failed the test based on threshold criteria. The first column in the display shows a summary of all the pass/fail criteria. If the mouse indicator is passed over a failed field, a pop-up box indicates the reason for the failure. Manual adjustments to peak integration and peak assignments may also be made in this view. Physicochemical property analysis results may also be sorted and viewed by experiment.

Step 6: Export data*Data Export Manager*

Data export is facilitated through the Data Export Manager. The format of the exported data file is defined in the data review window. The export file itself can be either xml, comma delimited, tab delimited or html format.

REFERENCES

- Valko, K.; Du, C. M.; Bevan, C. D.; Reynolds, D. P.; Abraham, M. H. J. *Pharm. Sci.* 2000, 89 (8), 1085-1096.
- Valko, K.; Du, C. M.; Bevan, C. D.; Reynolds, D. P.; Abraham, M. H. *Current Medicinal Chemistry.* 2001, 8(7).
- Kangas, H.; Kotiaho, T.; Salminen, T.; Kostianinen, R. *Rapid Communications in Mass Spectrometry.* 2001, 15: 1501-1505.
- Liu, H.; Carter, G. T.; Tischler, M. *Rapic Communications in Mass Spectrometry.* 2001, 15: 1533-1538.
- Camurri, Giulio; Zaramella, Alessio. *Anal. Chem.* 2001, 73 (15), 3716-3722.

WATERS CORPORATION
34 Maple St.
Milford, MA 01757 U.S.A.
T: 508 478 2000
F: 508 872 1990
www.waters.com

Waters
RIGHT ON TIME.

CONCLUSIONS

- The ProfileLynx Application Manager is well suited to do multiple physicochemical property indicator measurements simultaneously.
- The ProfileLynx Application Manager easily interfaces with existing data systems and can report physicochemical data using a graphical interface.
- The ProfileLynx Application Manager allows easy data management and review and allows full automation of physicochemical indicator measurements with little or no operator intervention.

FUTURE WORK

- Open Access walk-up type analysis for fully automated analysis with automated reporting and e-mailing of results to chemists.
- Inclusion of time-based metabolic and stability calculations.
- Integration of a fully customizable sample import tool to allow connection with any type of sample list format.
- Development of database interactivity tools to allow warehousing and retrieval of the physicochemical data.

ACKNOWLEDGEMENTS

- Andy Organ, Klara Valko; GlaxoSmithKline, UK
- Mark Hemling, H. Y. Cheng; GlaxoSmithKline, US
- Larry Mallis, Chris Petucci; Wyeth Ayerst, US
- Graeme Wells, Paul Chilvers, Maqsood Adia, David Deans; Waters, UK

Opt-In! My Profile
www.waters.com

Waters, Micromass, ProfileLynx, ZQ, Masslynx, ECP, MUX-technology and XTerra are trademarks of Waters Corporation.
All other trademarks are the property of their respective owners.
©2003 Waters Corporation Produced in the U.S.A. June 2003 720000669EN LW-PDF

