A METABONOMIC NEPHROTOXICITY STUDY UTILIZING LC/MS

Major, H.¹; Castro-Perez, J.¹; Wilson, I.²; Lenz, E.²; Davies, D.²

¹Waters Corporation, Manchester, UK

²Dept. of Drug Metabolism and Pharmacokinetics, AstraZeneca, Mereside, Alderley Park, Macclesfield, Cheshire, SK10 4TG, UK

INTRODUCTION

Metabonomics involves the study of time-related metabolic profile changes that can be the result of exposure to a toxin or drug, to environmental effects or the onset of disease. These studies are concerned with the complement of endogenous compounds rather than with xenobiotics. The primary goal is to identify and quantitate small molecules in biological systems that are affected as the direct result of an external stimulus. Knowledge of these compounds, or biomarkers, can then be used for diagnosis, screening, or to direct further research. The metabonomics approach is now being widely investigated by large pharmaceutical companies in the area of drug discovery and development where an early indication of toxicity is of paramount importance in preventing the late attrition of a potential drug candidate. ¹H NMR has been the primary method of choice for analysis of changes to the profiles of complex biofluids, such as urine, as a result of toxic lesions induced by xenobiotics. We have investigated the use of exact mass LC/MS for the determination of potential biomarkers of nephrotoxicity. The gentamicin complex used in this study is an aminoglycoside antibiotic complex and a nephrotoxin.

EXPERIMENTAL

Animal Study

- Male Wistar-derived rats (n=5 per group), approx. 140 g in weight acclimatized in metabolism cages for 3 days prior to treatment
- Food and water provided ad libitum
- One group dosed with gentamicin at 60 mg/kg twice daily, subcutaneous
- Urine samples collected daily for 9 days pre- and post-dose from control and dosed animals and stored at -20 °C prior to analysis
- The excretion pattern of small organic molecules in the urine was studied using LC/MS



Figure 1. Waters® MS Metabonomics System

HPLC Conditions

HPLC:	Waters 2795XC Separations Module				
Column:	Waters Symmetry [®] C ₁₈ , 2.1x100 mm, 3.5 µm				
Mobile Phase:	A: water + 0.1% formic acid				
	B: acetonitrile + 0.1% formic acid				
Flow Rate:	600 µL/min split to 120 µL/min to MS				
Column Temp:	40 °C				
Injection Volume:	10 µL				
Gradient:					
Time (min) 0 0.5 4.0 8.0 9.0 9.1	%A %B 100 0 100 0 80 20 5 95 5 95 100 0				

MS Conditions

MS	Waters Micromass® Q-Tof micro™	
Ionization Mode	Positive Ion Electrospray	
Capillary	3200 V positive	
Sample Cone	30 V	
Source Temperature	120 °C	

Desolvation Temperature	250 °C				
Cone Gas Flow	50 L/hr				
Desolvation Gas Flow	500 L/hr				
Argon Collision Gas					
MS Acquisition Parameters - LockSpray [™] Enabled					
Acquisition Range	m/z 50-850				
Acquisition Rate	0.4 sec				
Inter-scan Time	0.1 sec				
Mode	centroid				
Lock mass frequency	5.0 sec				
Lock scans averaged	10				
Lock Reference	Leucine enkephalin				
	0.5ng/µL in 1:1 acetonitrile: water + 0.1% formic acid				
Lock Mass	556.2771				

Data Processing

- The analysis of large batches of complex biological samples such as these can generate a wealth of data that require multivariate statistical analysis and pattern recognition methods.
- The MarkerLynx[™] Application Manager incorporates a peak deconvolution package which allows detection and retention time alignment of the peaks eluting in each chromatogram.
- Data is collected into a single matrix by aligning peaks with the same mass/retention time along with their associated normalized intensities.
- Finally principal component analysis (PCA) is performed.

Waters

- The ability to exclude xenobiotics automatically has also been incorporated.
- Option to export the detected masses and ion intensities to third party multivariate software packages e.g. SIMCA (Umetrics AB, Sweden) and Pirouette (Infometrix, USA) for further analysis such as PLS-DA (partial least squares discriminant analysis).



Figure 2. MarkerLynx Report

RESULTS

The MarkerLynx results for the positive ion data are shown in Figure 2 and consists of:

- A sample bar which lists the samples analyzed
- A marker bar which lists the masses and retention times with associated intensities of all the components detected in each of the samples
- A TIC chromatogram view of the selected sample
- A trend view of the selected component across all samples (m/z 255.0868, RT 2.95 min.)

- A scores plot showing the separation and/or clustering between samples (annotation turned off for clearer visualization of the results)
- A loadings plot indicating the m/z and RT values of the ions responsible for the clustering/separation i.e. potential biomarkers of nephrotoxicity

Gentamicin is an antibiotic complex produced by fermentation of *Micromonospora purpurea* or *M. echinospora* and consists of 3 closely related components, gentamicin C1, C2 and C1a as shown in Figure 3.



Figure 3. Structure of gentamicin complex

Under the HPLC conditions employed the gentamicin complex and any metabolites formed were not retained on the column and eluted with the solvent front. A spectrum and elemental composition report for the gentamicin complex, eluting at 0.3 min, from a day 9 dosed rat urine sample are shown in Figure 4.



Figure 4. Spectrum and elemental composition report for gentamicin complex

The gentamicin and any metabolites are readily removed from the PCA plot leaving only endogenous metabolites.



Figure 5. PCA Scores Plot

The scores plot shows separation of the day 5 - 9 dosed samples from the controls and predose samples and the associated loadings plot (Figure 6) shows the ions responsible.



Figure 6. PCA Loadings Plot

The major species responsible for the separation of the day 5 - 9 dosed samples from the predose and controls are tabulated in Tables 1 and 2.

	Table 1. Positive lons Increased after Dosing							
Retention Measured		Calculated	Elemental composition	Postulated structure				
Time (min) Mass (Da)		Mass (Da)	of [M+H] ⁺ ion					
	3.74	3.74 105.0319		105.0340	C7H5O	fragment of hippuric acid		
	1.16 215.0181		215.0165	C ₄ H ₃ N ₆ O ₅	C ₄ H ₂ N ₆ O ₅			
	2.95 255.0868		255.0869	C12H15O6	C ₁₂ H ₁₅ O ₆			
į		_						
	Table 2. Positive Ions Decreased after Dosing							
	Retention		Measured	Calculated	Elemental composition	Postulated structure		
	Time		Mass	Mass	of [M+H] ⁺ ion			
	3.59		206.0434	206.0453	C10H8NO4	xanthurenic acid		
	3.95		190.0499	190.0504	C ₁₀ H ₈ NO ₃	kynurenic acid		
	1.02		267.1320	267.1345	C ₁₃ H ₁₉ N ₂ O ₄	C ₁₃ H ₁₈ N ₂ O ₄		
	1.08		245.1599	245.1614	C10H21N4O3	C ₁₀ H ₂₀ N ₄ O ₃		
	5.57		255.0653	255.0657	C ₁₅ H ₁₁ O ₄	C ₁₅ H ₁₀ O ₄		
	2.35		297.1447	297.1450	C ₁₄ H ₂₁ N ₂ O ₅	C14H20N2O5		
	2.73		220.1194	220.1185	C ₉ H ₁₈ NO ₅	pantothenic acid		
	4.70		149.0618	149.0603	C ₉ H ₉ O ₂	cinnamic acid		
	5.64		285.0746	285.0763	C ₁₆ H ₁₃ O ₅	C ₁₆ H ₁₂ O ₅		

For a few concentrated components where the ions are particularly intense the mass measurements are lower than expected. This is due to detector saturation, otherwise they are typically within 5 ppm of the calculated mass for the postulated elemental composition.

MS/MS exact mass analysis has been used for structural elucidation of the potential biomarkers. The MS/MS product ion spectrum from the ion at m/z 190, which was observed to decrease in intensity after dosing, is shown in Figure 6. This was shown to be kynurenic acid (4-hyroxyquinoline-2-carboxylic acid), a metabolic product of tryptophan.



Figure 7. MS/MS product ion spectrum from m/z 190

Xanthurenic acid, m/z 206, also identified as decreasing after dosing is also part of the tryptophan catabolism pathway (Figure 7). The tryptophan derived isoquinolines, kynurenic acid and xanthurenic acid, are not degraded further but are excreted in urine and are partly responsible for the yellow color of urine.



Figure 8. Part of Tryptophan Catabolism Pathway

Pantothenic acid, also seen to decrease after dosing, is a building block of coenzyme A. Work is ongoing to elucidate the structures of the other ions determined as being significant in the loadings plot and to evaluate their toxicological significance.

CONCLUSIONS

- LC/MS in conjunction with PCA analysis has been successfully used to screen rat urine after dosing with gentamicin.
- The MarkerLynx Application Manager simplifies the processing by incorporating peak deconvolution and data alignment with PCA in one software package.
- The control and pre-dose samples could easily be differentiated from the day 5 9 dosed samples.
- The m/z values of the ions responsible for the PCA separation were identified.
- MS/MS exact mass was used for structural elucidation of the potential biomarkers.

- Xanthurenic acid and kynurenic acid, part of the tryptophan catabolism pathway, were identified as decreasing after dosing.
- LC/MS data complementary to NMR which had identified glucose, which increased after dosing, as the main biomarker of nephrotoxicity.
- Further MS/MS structural elucidation is required for identification of the other potential biomarkers and their toxicological significance needs to be determined.

Vaters

Sales Offices:

AUSTRIA AND EXPORT (CENTRAL EUROPE, CIS, MIDDLE EAST, INDIA AND INDIA SUBCONTINENT) 43 1 877 18 07

AUSTRALIA 61 2 9933 1777

BELGIUM AND LUXEMBOURG 32 2 726 1000 BRAZIL 55 11 5543 7788

U.S.A. AND ALL OTHER COUNTRIES:

WATERS CORPORATION 34 Maple St. Milford, MA 01757 U.S.A. T: 508 478 2000 F: 508 872 1990 www.waters.com

WATERS CORPORATION 34 Maple St. Milford, MA 01757 U.S.A. T: 508 478 2000 F: 508 872 1990 www.waters.com







Waters, Micromass, Symmetry, Q-Tof micro, LockSpray and MarkerLynx are trademarks of Waters Corporation.

All other trademarks are the property of their respective owners. ©2003 Waters Corporation Produced in the U.K. November 2003 720000771EN CH-PDF

