## Waters

## A NOVEL METHOD FOR HTS SCREENING FOR PHASE II CONJUGATIONS FOR *IN VIVO* SAMPLES WITH AN ORTHOGONAL QUADRUPOLE TIME OF FLIGHT MASS SPECTROMETER USING NEUTRAL LOSS WITH EXACT MASS

G Harland', J Castro-Perez', I Beattie<sup>2</sup>, K Joncour<sup>2</sup> 'Waters Corporation, Atlas Park, Simonsway, Manchester, M22 5PP, UK. <sup>2</sup>AstraZeneca, R&D Charnwood, Bakewell Road,Loughborough, Leics. LE11 5RH UK.

## Introduction

During the drug discovery stages, screening and identification of *in-vitro* metabolites plays an important role, but remains very challenging due to the different complexity of structures being analyzed. LC/MS/MS is now a well-established technique and it plays a vital role in this early phase of Drug Discovery. With new advances in MS technology, and especially Q-Tof<sup>™</sup> technology and its exact mass capabilities, good accurate data will be ensured, which is very important for metabolite identification. By the use of exact mass data, and determination of metabolite structures it will allow medicinal chemists to make the necessary 'structural tuning' to achieve the desired compound activity. Moreover, as the demand to screen a large number of compounds in Drug Discovery increases, it is imperative to have information about Phase I and II reactions in one single experiment. In turn, this will facilitate compound selection for the next stage of drug development. Moreover, it will also provide a better predictive platform when the compound goes in-vivo. Incubation levels are also important since we want to have a reasonable concentration which may represent the levels seen in vivo. For this work we have developed a new approach to the detection of Phase II conjugations in which we can monitor Phase II conjugations such as glucuronides and so forth from a single LC run (Figure 1). To do this, conjugation moieties form a characteristic loss when analyzed by LC/MS/MS. This neutral loss can be monitored with great accuracy using exact mass. From this single injection the mass spectrometer applies two collision energies one at a lower level and one at a higher level (Figure 2, 3 and 4). This will also allow detection of Phase I metabolites in the low level collision energy. If a loss of a conjugation such as glucuronide with an exact mass of m/z 176.0321 ±20 mDa is detected then the instrument will automatically switch to carry out MS/MS and confirm the loss.

Otherwise, it continues switching between the two collision energy modes. In this paper, we will show the analytical strategy utilized to detect putative Phase II reactions based on the exact neutral loss data from a bile sample containing Midazolam and Fluanisone (Figure 5).

Test compound

### **Methods**

### Samples

Sample:

				Midazola Fluanison	ım and e		
Metabolizing System:				In vivo rat bile sample			
Dilution Factor:				Sample diluted in water			
				1/10			
LC con	dition	5					
Solvent Delivery System: Waters® 2795							
, ,				Separations Module			
Mobile Phase A:				Water +0.2 % Formic			
				Acid			
Mobile Phase B:			1	MeCN + 0.2 % Formic			
				Acid			
Gradie	nt:						
Time	A%	В%	C%	D%	Flow	Curve	
0.00	98.0	2.0	0.0	0.0	0.15	1	
1.50	98.0	2.0	0.0	0.0	0.15	1	
15.00	30.0	70.0	0.0	0.0	0.15	6	
17.00	5.0	95.0	0.0	0.0	0.15	6	
19.00	5.0	95.0	0.0	0.0	0.15	6	
19.10	98.0	2.0	0.0	0.0	0.15	6	
25.00	98.0	2.0	0.0	0.0	0.15	6	
Flow R	ate:			150 µL/n	nin		
Column:				Waters Atlantis <sup>™</sup> dC <sub>18</sub>			
				150 x 2.	1 mm id	3.5µ	
Injection Volume:				10 µL		·	
Mass Spec Conditions							

Waters Micromass®
Q-Tof micro <sup>™</sup> with
LockSpray™
ESI + ion mode

Centroid

20 and 45 eV

concentration

Leucine Enkephalin,

@ 556.2771 lng/µL

2.8 kV

100 °C

25 V

## 



High energy

## Results

M+H

Lock Mass:

Acquisition:

Cone Voltage:

Capillary Voltage:

Source Temperature:

Collision Energy Profile:

Desolvation Temperature: 250 °C



Figure 1. Neutral loss with a single injection.



Figure 2. Exact mass neutral loss on a Q-Tof.









• All major glucuronide metabolites were detected by using this approach and also other Phase I metabolites were detected. The overall mass accuracy for all detected xenobiotics was better than 2.1 ppm RMS (Table 1).



13.43 14.77 15.46

# Waters

Compound name	m/z observed	m/z theoretical	mDa +/-	ppm +/-
Midazolam	242.0904	242.0900	0.9	2.5
M+32	342.0001	342.0009	0.0	2.5
M+34	360.0915	360.0909	0.6	1.7
M+192	518.1130	518.1117	1.3	2.6
M+208	534.1079	534.1068	1.1	2.1
Fluanisone				
M+16	373.1927	373.1920	0.7	2.0
M-14	343.1822	343.1809	1.3	3.7
M+192	549.2248	549.2233	1.5	2.8
			Average mDa	RMS ppm
			0.9	2.1

**Poster**REPRIN1

Table 1. Mass accuracies for Midazolam and Fluanisone metabolites.

 In addition to the detection and confirmation of the glucuronide conjugations, we were also able to detect Phase I metabolites from the low CE energy scan from the two-dosed compounds (Figure 6).



Figure 6. XIC for Phase I metabolites using exact mass DDA<sup>™</sup> mode.

## Conclusions

- This is a highly selective and sensitive method for screening of specific conjugations by exact mass
- By this approach information from Phase I and other Phase II conjugations are also kept because the data is acquired in full scan exact mass
- With this approach there are no limitations to the number of exact mass neutral loss acquisitions which may be carried out on a single scan because the experiments are based upon DDA

WATERS CORPORATION 34 Maple St. Milford, MA 01757 U.S.A. T: 508 478 2000 F: 508 872 1990 www.waters.com







Waters, Micromass, Atlantis, Q-Tof, Q-Tof micro, LockSpray and DDA are trademarks of Waters Corporation. All other trademarks are the property of their respective owners. ©2003 Waters Corporation Produced in the U.S.A. October 2003 720000773EN SH-PDF

