ASMS

 \bigcirc

Protocols to Assure Reproducible Quantitative and **Qualitative Analysis of Tryptic Digests of Complex P**

Craig A. Waters Corporation, Milford, MA, USA Dorschel, Marc V. Gorenstein, Guo-Zhong Li, Jeffrey \mathbf{O} Silva, Scott

Overview

mixtures for the purpose of quantitative and qualitative proteomics. Purpose: Demonstrate reproducibility of tryptic digestion of complex protein

mass MS mode of acquisition. mended practice in a novel alternating low and elevated collision energy exact **Methods:** Replicate digestions of rat serum were carried out using a recom-mended practice. LC/MS analysis was performed, again following a recom-

digestion. and chromatography was obtained, demonstrating the reproducibility of the **Results**: A high degree of precision of signal intensity, mass measurement,

Introduction

tion of peptides should be consistent from run to run for valid quantitative comparisons. Maximizing the reproducibility of chromatographic retention times simplifies the overall task of correlating the massive data sets obtained in global proteomic experiments. protein in two samples. Accurate measurement of peptide masses allows more stringent searching of protein databases resulting in faster and more accurate qualitative identification of proteins. Ionizait is necessary that great care be taken in preparation of samples and in obtaining exact mass LC/ MS data. Digestions with proteases such as trypsin should be carried to completion, so that the ratio of concentrations of fragment peptides will accurately reflect the ratio of concentration of the original To obtain meaningful results in qualitative and quantitative proteomics, such as observation of differ-ent levels of expression of one or more proteins in a set of samples and identification of the proteins,

of rat serum using the Waters® Protein Expression System. We have developed protocols and software tools to enable us to reproducibly digest complex pro-tein mixtures, routinely obtain precise mass measurements on peptides (typically ≤ 5 ppm error), minimize run-to-run retention time drift at the capillary separation scale, and demonstrate the quality and reproducibility of our results. We illustrate this with digestion and analysis of replicate samples

Methods

Digests

Replicate digests of 5 µL (250 – 300 µg protein) aliquots of rat serum (Sigma, catalog num-ber R 9759) were prepared according to a recommended procedure employing a proprie-tary detergent and in which cystines are reduced and alkylated with iodoacetamide prior to digestion.

Samples

gests and 25 µL of a standard mixture of purified tryptic digests of five proteins. Seven samples were prepared by combining 75 μ L of each of seven replicate rat serum di-

Data Collection

LC/MS System: Waters Protein Expression System comprised of the Waters CapLC[®] System with the Waters Micromass[®] Q-Tof Ultima API Mass Spectrometer equipped with a NanoLockSpray[™] Source operated at 12,000 mass resolving power

•Column: Waters NanoEase[™] Atlantis[™] dC18 Column, 300 µm x 15 cm •Mobile Phase: A = 1% Acetonitrile in Water, 0.1% Formic Acid, B = 80% Acetonitrile in

and 20 min. re-equilibration at initial conditions Water, 0.1% Formic Acid •Gradient: 6% to 40% B over 100 min. at 4.4 µL/min, followed by 10 min. rinse (99% B)

Three replicate injections of 2.5 µL of each digestion sample solution were made, directly onto the column (no desalting trap was used). Mass spectral data were collected in alternating low and elevated collision energy using a1.8 second acquisition at each energy. The NanoLockSpray source was switched every 10 seconds to obtain a scan of the accurate mass standard ([Glu¹]-Fibrinopeptide B).

Expression System Informatics software tools developed in-house (**Poster TPR 354**, Li et al., "A Novel Algorithm to Track lons or Peptides Found by LC/MS in Multiple Injections of Multiple Complex Biological Samples"). Additional data analysis was performed using Microsoff[®] Excel and Spotfire[®] DecisionSite™ The raw data were processed using functions in ProteinLynx™ Global SERVER and Protein

Recommended Procedure for Tryptic Digestion of Ser 1. Dispense 5 µL of serum (250 – 300 vg protein) in a capped microcentrifu

Add 20 µL of 100 mM aqueous[1] ammonium bicarbonate. Add 25 µL of 0.2% solution of RapiGest[™] SF in water (add 500 µL of wc and vortex.

21. 22. 3. 5. 5. 7. 7. 7. 7. 8 Place tube in a block heater set at 80° C. Heat for 15 minutes, vortexing Remove from block and centrifuge to return condensate to bottom of tube Add 2.5 µL of 100 mM dithiothreitol (15.4 mg/mL in water), vortex. Place tube in a block heater set at 60° C. Heat for 30 minutes.

of tube. 9. Add Add 2.5 µL of 200 mM iodoacetamide (55.5 mg/mL in water), vortex. Remove tube from block, allow to cool to room temperature, and centrifuç

Place sample in dark at room temperature and allow 30 minutes reaction tin
Add 50 µL of a solution of Promega trypsin in 50 mM ammonium bicarbonc mM ammonium bicarbonate to one 20 µg vial of trypsin. Trypsin:protein ratio w
Place tube in a block heater set at 37° C. Incubate overnight.
Centrifuge to return condensate to bottom of tube.
Add 100 µL of water to the sample, vortex. Final sample concentration is eq tein/µL. Final RapiGest SF concentration will be 0.025%, and the RapiGest is n be refrigerated for short term storage, frozen for long term storage.
Water used in this procedure was purified by a Milli-Q A10 system (Millipor

Results

like peptides from multiple sample injections based on their unique exact r and compares their relative intensity. Peptides matched in this manner are Tryptic digest of a serum sample can be expected to produce a highly complex mixture of peptides. Waters Protein Expression System Informatics tools have been developed to extract the exact mass, tion time and intensity of each detected peptide or chemical entity in each data set. This exact mass retention time pair or EMRT. retention time information serves as a unique signature for each peptide. mass retention time signature Additional software matches to extract the exact mass, retended to a set. This exact mass and referred to as an exact mass

Representative base peak intensity chromatograms from 3 of the *7* samples (**Figure 1**) give a rough indication of the complexity of the samples and the reproducibility of the chromatography.



of the doubly charged precursor is within 10 ppm of the accurate mass (679.3265). The insert spectrum is the corresponding ele-vated energy spectrum in each case. Note that the precursor intensity and fragmentation pat-**Figure 2** represent the mass and intensity of the doubly charged precursor in each sample. virtually identical. The experimental exact mass tern, showing nearly the entire y'' series, The highlighted peak in the chromatograms represents peptide T73 of rat albumin, having sequence TVMGDFAQFVDK. The spectra in are

one injection of each sample gives another overview of sample similarity (Figure 3). Selected ion chromatograms for 6 masses in



Geromanos, and Timothy Riley

ipore Corporation).	s equivalent to 1.25–1.5 µg pro- is not hydrolyzed. Sample should	n time. sonate, vortex. (Add 400 µL of 50 'o will be approximately 1:100.)	ye to return condensate to bottom	periodically.	a ge tube, volume greater than 200
periodically. ge to return condensate to bottom n time. bonate, vortex. (Add 400 µL of 50 io will be approximately 1:100.) is equivalent to 1.25–1.5 µg pro- t is not hydrolyzed. Sample should	periodically. ge to return condensate to bottom n time. bonate, vortex. (Add 400 µL of 50 io will be approximately 1:100.)	periodically. ge to return condensate to bottom	periodically.		

Figure 2. Rat Albumin T73 in 3 Samples



Number of Occurrences

Mean Cv = 10%

3000

2000

1000

300

0

0

Cv of EMRT within

Sample (%)

54-

6-

10 20 30 40 50 Cv for EMRT in Experiment (%)

8-

100

200



EMRTs detected in 16 to 21 of the 21 injections account for 90% of the summed MS signal intensity for all the EMRTs. The low intensity of the less frequently detected EMRTs suggests these are near the limit of detection or are random events of chemical noise.

Off-diagonal plots correlating intensities of individual EMRTs in pairs of samples show a strong correla-tion in each instance. **Figure 7** is a representative example for sample 1 and sample 3.







The coefficients of variation for the intensities of each replicating EMRT were calculated, first for each sample (3 injections) and then for the entire experiment (21 injections). Within a sample (**Figure 8**), CV values are low, and reflect the fact that injection volume and ionization efficiency are very reproducible. CV values for the entire experiment (**Figure 9**) include such factors as sampling of the serum and diges-tion efficiency. These CV values also indicate excellent reproducibility. (Note that, for example, an EMRT found in all 21 injections will contribute 7 values to **Figure 8** but only one value to **Figure 9**.)































1000-

500

1500





Figure 8.

Intensity CV Frequency

Figure 9.

Intensity CV Frequency

for EMRTs in Samples































Figure

4500-

1

4000-





(M+H)⁺ Monoisotopic Mass

2500

8

00

2000

3000

8.

3500-

ŧ,

ų.

11

n,



dian Cv =

7%

The highest EMRT intensity CV values are associated with EMRTs of low intensity (Figure 10).



Mass measurement precision of replicating EMRTs is excellent. A great majority of such measurements have a precision within ± 5 ppm (**Figure 11**). The histogram of summed EMRT intensity for each mass precision bin (**Figure 12**) shows that the least precise mass measurement is associated with trace components.







While great care was taken to assure that the gradient composition was consistent from injection to injec-tion, some run-to-run variability in retention time was observed. The effect of even small compositional "ripples" in the gradient profile is a shift in retention time of all peptides eluting at that composition. The informatics tools used to correlate EMRTs from injection to injection are able to accurately account for this retention time shift, in effect greatly reducing the retention time variability in the data set and increasing the confidence in the correlated data. This is explored in much greater depth in **Poster TPR 354**.

Mass Precision (ppm)

Conclusions

we have demonstrated that: Through the use of the Waters Protein Expression System and Informatics tools,

- performed such that peptides produced vary in intensity with coefficients of Tryptic digestion of complex protein mixtures (sera) can be reproducibly variation typically well below 30%.
- Electrospray ionization efficiency of peptides is reproducible from run to run and sample to sample.
- Exact masses of peptides can be determined with precision typically less than 5 ppm throughout a series of samples.