

Waters Corporation, MS Technologies Center, Manchester, UK

### Introduction

All tandem quadrupole instruments are capable of monitoring both positively and negatively charged ions formed by analytes entering the ionization source. Occasionally a single assay will entail a mixture of compounds some of which produce positive charged ions and some of which produce negatively charged ions. However, the time taken to switch between these two modes of operation will generally prohibit the use of a single LC/MS/MS analytical run to monitor this mixture.

Tandem quadrupole mass spectrometers typically require a voltage settling period of 200 to 300 ms when switching from positive to negative ion modes (some tandem Quadrupole MS instruments require >700 ms). When performing an assay in multiple reaction monitoring (MRM) mode it is important to keep the MRM cycle time low enough to maintain a sufficient number of data points per chromatographic peak (typically 15-20 points per peak per MRM

channel monitored). If this is not done the quality of the resulting quantification data will decrease due to poor peak area definition. On most instruments polarity-switching time is too long to allow positive ion and negative ion monitoring to be used together effectively in a single analytical run while maintaining a low MRM cycle time.

The Waters Micromass Quattro Premier™ benchtop mass spectrometer has an unmatched capability to switch from monitoring positively charged ions to monitoring negatively charged ions in less than 20 ms without any degradation in data quality. Optimized high voltage supplies used in the Quattro Premier mean that positive to negative ion switching in full scan or MRM mode can now be achieved in 20 ms with no loss in performance. Mixed mode, short cycle time MRM analysis can now be performed without loss in signal intensity and without compromising the integrity of quantification data.

## Example Data

The data presented below was obtained using the Quattro Premier benchtop mass spectrometer. Direct loop injections of a solution containing sulfadimethoxine were monitored in the MRM acquisition mode. Separate acquisitions were performed in positive ion mode only (Figure 1) and negative ion mode only (Figure 2). A third acquisition was performed using a single method combining both positive and negative acquisition

methods using a polarity switching time of 20 ms (Figure 3). By comparison it can be seen that neither peak areas nor signal-to-noise ratios are affected by rapid polarity switching.

With the Waters Quattro Premier mass spectrometer you can confidently perform mixed ionization mode MRM analysis without compromising either sensitivity or MRM cycle time.




Figure 1. Replicate injections of sulfadimethoxine, monitored in positive ion mode only.




Figure 3. Replicate injections of sulfadimethoxine, monitored simultaneously in positive and negative ion modes.




Figure 2. Replicate injections of sulfadimethoxine, monitored in negative ion mode only.

## Sales Offices:

**AUSTRIA** 43 1 877 18 07

**AUSTRALIA** 61 2 9933 1777

**BELGIUM AND LUXEMBOURG** 32 2 726 1000

**BRAZIL** 55 11 5543 7788

**CANADA** 800 252 4752 X2205

**CZECH REPUBLIC** 420 2 617 11384

**DENMARK** 45 46 59 8080

**FINLAND** 358 9 506 4140

**FRANCE** 33 1 3048 7200

**GERMANY** 49 6196 400600

**HONG KONG** 852 29 64 1800

**HUNGARY** 36 1 350 5086

**INDIA** 91 80 2837 1900

**IRELAND** 353 1 448 1500

**ITALY** 39 02 27 4211

**JAPAN** 81 3 3471 7191

**KOREA** 82 2 820 2700

**MEXICO** 52 55 5524 7636

**THE NETHERLANDS** 31 76 508 7200

**NORWAY** 47 6 384 6050

**PEOPLE'S REPUBLIC OF CHINA** 86 10 8451 8918

**POLAND** 48 22 833 4400

**PUERTO RICO** 787 747 8445

**RUSSIA/CIS** 7 095 931 9193

**SINGAPORE** 65 6278 7997

**SPAIN** 34 93 600 9300

**SWEDEN** 46 8 555 11 500

**SWITZERLAND** 41 62 889 2030

**TAIWAN** 886 2 2543 1898

**UK** 44 208 238 6100

**US** 800 252 4752

# Waters

WATERS CORPORATION

34 Maple St.

Milford, MA 01757 U.S.A.

T: 508 478 2000

F: 508 872 1990

[www.waters.com](http://www.waters.com)



## For Complete Confidence

Waters, Micromass, Quatro Ultima, Quattro micro and Quattro Premier are trademarks of Waters Corporation.  
All other trademarks are the property of their respective owners.  
©2004 Waters Corporation Produced in the U.S.A. June 2004 720000914EN LL-PDF

