

Waters Corporation, MS Technologies Center, Manchester, UK

Introduction

Due to advances in HPLC and the arrival of UPLC™, chromatographic peaks are becoming extremely narrow. It is generally accepted that if reproducible peak area quantification is required chromatographic peaks should be defined by no less than 15 data points. To maintain this peak definition as peaks become narrower, tandem quadrupole mass spectrometers need to acquire faster. The Waters® Micromass® Quattro Premier™ benchtop mass spectrometer has been designed to operate in the multiple reaction monitoring (MRM) mode with the shortest of acquisition cycle times (10 milliseconds per data point, 100 data points per second).

An inter-channel delay period is required between successive MRM transitions to allow ions to be cleared from the collision cell. If two successive MRM transitions share a common daughter ion and the cell is not effectively cleared during the inter-channel delay a false signal contribution will occur between the transitions. This undesirable phenomenon is generally referred to as crosstalk.

At short MRM cycle times the ion transit time of ions in standard RF only collision cells (several milliseconds) is a significant problem in relation to crosstalk. The Quattro Premier mass spectrometer utilizes the patented T-Wave™ technology in the collision cell, enhancing the rapid acquisition MRM performance of the instrument. This new design of collision cell ensures that the undesirable phenomenon of MRM inter-channel crosstalk is virtually eliminated, even at the shortest of MRM cycle times.

Example Data

The data presented (Figure 1) is for the LC/MS/MS (MRM) analysis of verapamil. To measure inter-channel crosstalk a second dummy channel was simultaneously monitored (the two transitions share a common daughter mass). Rapid MRM acquisition was performed with a 10 ms dwell time and 10 ms MRM channel was considered to be a result of crosstalk occurring in the collision cell.

The data presented demonstrates that the travelling wave rapidly clears the collision cell between transitions, virtually eliminating crosstalk, even when using a very short inter-channel delay, enhancing the fast acquisition performance of the Quattro Premier mass spectrometer.




Figure 1. Rapid cycle time LC/MRM analysis of verapamil, demonstrating almost complete absence of inter-channel cross talk.

Sales Offices:

AUSTRIA 43 1 877 18 07

AUSTRALIA 61 2 9933 1777

BELGIUM AND LUXEMBOURG 32 2 726 1000

BRAZIL 55 11 5543 7788

CANADA 800 252 4752 X2205

CZECH REPUBLIC 420 2 617 11384

DENMARK 45 46 59 8080

FINLAND 358 9 506 4140

FRANCE 33 1 3048 7200

GERMANY 49 6196 400600

HONG KONG 852 29 64 1800

HUNGARY 36 1 350 5086

INDIA 91 80 2837 1900

IRELAND 353 1 448 1500

ITALY 39 02 27 4211

JAPAN 81 3 3471 7191

KOREA 82 2 820 2700

MEXICO 52 55 5524 7636

THE NETHERLANDS 31 76 508 7200

NORWAY 47 6 384 6050

PEOPLE'S REPUBLIC OF CHINA 86 10 8451 8918

POLAND 48 22 833 4400

PUERTO RICO 787 747 8445

RUSSIA/CIS 7 095 931 9193

SINGAPORE 65 6278 7997

SPAIN 34 93 600 9300

SWEDEN 46 8 555 11 500

SWITZERLAND 41 62 889 2030

TAIWAN 886 2 2543 1898

UK 44 208 238 6100

US 800 252 4752

Waters

WATERS CORPORATION

34 Maple St.

Milford, MA 01757 U.S.A.

T: 508 478 2000

F: 508 872 1990

www.waters.com



For Complete Confidence

Waters, Micromass, Quattro Premier, T-Wave and UPLC are trademarks of Waters Corporation.
All other trademarks are the property of their respective owners.
©2004 Waters Corporation Produced in the U.S.A. June 2004 720000917EN LL-PDF

