

Peter Alden, Steven A. Cohen, and Scott J. Berger Life Sciences R&D, Waters Corporation, Milford, MA

Overview

- Analysis of intact proteins by MS is a powerful approach for protein identification and assignment of protein modifications.
- Analyzing protein mixtures typically requires chromatographic processing to separate, concentrate, and desalt components prior to MS detection.
- Optimization of eluent composition plays an important role in LC/MS methods development, as sorbents possess differing selectivity and performance based on gradient shape and choice of acidic modifier.
- The positive chromatographic performance of an acidic modifier must be balanced with the potential for suppressing electrospray ionization.
- This poster describes a methodology where separation and ionization issues for protein LC/MS analysis are addressed in an efficient and automated fashion.







System: Waters[®] AllianceTM Bioseparations System composed of a 2796 Bioseparations Module, Column Heater Module, 2996 Photodiode Array Detector, ZQTM Mass Detector, and MassLynxTM Data System.

Chromatography:

BioSuiteTM pPhenyl Column (2.0 x 100mm,

1000A, 10µ polymeric particle) and Symmetry 300^{TM} C₁₈ Column (2.1 x 75mm, 300A, 5µ silica particle) operated at 150 µl/min. Separations began with a 5 min hold at 5% acetonitrile, followed by a 30 min linear gradient to 90% acetonitrile, maintaining a constant concentration of acidic modifier.

Sample: A mixture of seven proteins [bovine insulin (4 µg, 6 kD), horse cytochrome c (3 µg, 12 kD), bovine ribonuclease A (16 µg, 14 kD), horse myoglobin (2 µg, 17 kD), yeast enolase (3 µg, 47 kD), BSA (6 µg, 66 kD, apo-Transferrin (6 µg, 80 kD)) were prepared in 5% acetonitrile, and injected (10 µl) for each run.







Suppression Issues. A six protein mixture was analyzed by LC (polymeric RP)



Figure 4: Comparison of Acidic Modifiers for Protein LC/MS Separations-**Chromatographic Issues.** If the y-axis for each run is independently scaled, we can see from both the UV and MS traces that the choice of acidic modifier will have noticeable effects on peak retention time, column selectivity, and chromatographic resolution.

Automated Optimization of Intact Protein LC/MS Analyses

Figure 5: Effects of Acidic Modifier Concentration on LC/MS Separations. Varying acidic modifier concentration affects both the quality of the separation, and the extent of signal suppression. Proteins: (1) Rnase A, (2) Insulin, (3) CytC, (4) apotransferrin, (5) BSA, (6) Myoglobin, (7) Enolase.



Figure 6: Comparison of various acidic modifiers for LC/MS performance. The peak areas for the latest eluting component (Enolase) were determined for LC/MS analyses using a polymeric (LEFT), and silica based (RIGHT) reversed phase chemistry under a variety of acidic modifier concentrations. Similar performance is observed with both chemistries under the various conditions, except that MS response from the silica based column is more significantly improved with intermediate formic acid concentrations. Other acids show only increasing suppression with increasing concentration.

- An AutoBlendTM methodology permitted testing multiple LC/MS separation methods in an automated fashion.
- The type and concentration of an acidic modifier produce significant MS sensitivity and LC separation quality effects.
- The MS suppression effects are generally comparable for both polymeric and silica based RP separation chemistries.
- Formic acid demonstrates the least suppression of the acids tested, but resulted in poorer chromatographic performance.
- Higher formic acid concentrations increased MS response by producing narrower (more concentrated) peaks, without suppressing signal to a greater extent.
- This formic acid effect was more pronounced with a silica column than with a polymeric based RP chemistry.
- necessary.

Conclusions

• When building an LC/MS method for a protein mixture, tradeoffs between chromatographic and MS performance may be