

Asish B. Chakraborty, Steven A. Cohen, and Scott J. Berger Life Sciences R&D, Waters Corporation, Milford, MA

Overview

- Noncovalent protein complexes (NPC's) are the functional biological unit for many cellular activities.
- Thus, characterization of NPC's has become a major research effort in the post-genomic era.
- Ionization and detection of weakly associated NPC's is a challenging task, since a gentle transition of the complex from solution to the gas phase requires significant attention during methods development.
- Optimized electrospray conditions permit direct detection of intact NPC's, or its constituent subunits.
- Online separations can be coupled to ESI-MS to derive additional insights about complex structure and composition.

Introduction

The desire to describe the molecular machinery within the cell has led researchers to systematically characterize NPC's. Mass analysis of the intact complex is useful for discerning subunit stoichiometry and ligand interactions, while protein and peptide level analysis is more effective at characterizing and localizing subunit modifications. The gas phase ionization and detection of NPC's by MS methods is a challenging task since the transition of the complex between the solution and gas phase requires effective desolvation without disrupting complex structure (1,2). Furthermore, NPC analysis can be coupled with online separations to concentrate samples or separate complicated mixtures (3). We have developed a methodology to detect and characterize a tetrameric yeast alcohol dehydrogenase (ADH) complex using online SEC/MS analysis, while analyzing constituent subunits by a 2-D "heartcut" SEC/RP/MS methodology.

Yeast ADH forms a homotetrameric complex where each monomer (~37kD) is liganded to two divalent zinc ions. One zinc is essential for the activity and other is required to maintain structural integrity. Thus, the tetrameric complex can be expected to contain eight zinc ions for an intact complex mass of 147505 Da.





Mass Spectrometry Conditions (adapted from Ref. 4.5):

which be a set of the		
	MS System	Waters Micromass
	Ion Mode	Positive
	ESI source voltage	3 kV
	Cone Voltage	140 V (1D SEC/MS),
		45 V (2D SEC/RP/MS &
	Desolvation/Source Temp	200 °C /110 °C
	Post source pressure	7.1mbar (SEC-MS),
		1.7 mbar (SEC/RP
	Data Acquisition	500-10000 (m/z)
	Calibrant	2 mg/ml cesium iodide i
	Software:	Waters MassLynx TM 4.0

s[®] LCTTM ESI-Tof

& RP/MS)

P/MS)

in 50 % IPA & MaxEnt 1



Figure 4: RP/MS of the ADH complex. The TIC (Top) from the RP/MS analysis of yADH shows a high degree of complexity. If we sum spectra for the entire gradient (Bottom, inset), we produce a deconvoluted spectrum which reveals the presumptive 37kD yADH monomer, as well as a large number of lower mass species.

Figure 5: Characterization of ADH subunits by online "heartcut" 2D LC (SEC/RP)/ESI-Tof-MS. Combining SEC with the RP separation into an online 2D "heartcut" separation permits us to selectively transfer the ADH complex to the 2nd dimension for RP analysis. Thus, everything we observe in the 2D SEC/RP/MS run is directly associated with the complex observed during SEC/MS analysis.

Timing for the valve switches at the beginning and end of the SEC peak was derived from the initial 1D SEC/MS separation (Top), and produced a single RP/MS peak (Middle), that revealed only the expected monomeric ADH subunit (**Bottom**). The observed subunit mass is ~57 Da larger than predicted from the sequence, indicating a possible protein modification or sequence variation.

Conclusions

• The combination of SEC/MS and "heartcut" 2D SEC/RP/ MS experiments provided a fast and efficient approach for characterizing an ADH noncovalent protein complex, and it's constituent subunits.

• Online SEC-MS showed the tetrameric structure of the complex, while the heartcut 2D-MS analysis identified the subunits as homogeneous 37 kD components.

- Reversed phase analysis of the ADH complex was insufficient to describe the subunit structure, as low molecular weight contaminants were also present in the sample.
- The use of the SEC/RP heartcut separation permitted direct analysis of only those proteins present within the complex, while preventing these putative degradation products and/or non-ADH impurities from interfering with the subunit analysis.
- In future, this methodology should be directly applicable to study a wide range of nocovalent biomolecule interactions.

References

- 1. J.A. Loo, Int J Mass Spectrom 200, 175-186 (2000).
- 2. F. Sobott and C.V. Robinson, Curr Opin Struct Biol 12,729-734 (2002)
- 3. S. Martinovic, S.J. Berger, L. Pasa-Tolic, R. D. Smith, 72,5356-5360 (2000).
- N. Tahallah, M. Pinkse, C.S. Maier and A.J.R. Heck, Rapid Commun. Mass Spectrom. 15, 596-601(2001).
- B.N. Green, R.S. Bordoli, L.G. Hanin, F.H. Lallier, A. Toulmond and S.N. Vinogradov, J Biol Chem 274, 28206-28212 (1999).

Acknowledgements

We would like to thank Brian Green, Iain Campuzaano and Andy Jarrell from Waters for their valuable advice and assistance in obtaining optimal MS detection conditions for noncovalent protein complexes.