

Waters Corporation, MS Technologies Centre, Manchester, UK

## Introduction

The Waters® LCT Premier™ Mass Spectrometer is the highest performance bench-top orthogonal acceleration time-of-flight (oa-TOF) LC/MS system available today. Improvements in ion source, transfer optics, vacuum system, TOF optics and electronics have resulted in an oa-TOF instrument with unsurpassed ion transmission that gives an order of magnitude more sensitivity when compared to the previous Waters LCT™ Mass Spectrometer. With these improvements, Waters provides an instrument that has ability to detect analytes with full spectral acquisition that would not have normally been detected without a pre-concentration step.

Analytical application areas that benefit from sensitivity improvements include:

- Metabolite identification
- Impurity profiling in chemical synthesis
- Trace level component analysis in food or environmental samples
- Peptide/protein applications requiring full spectral sensitivity
- Natural product identification

## Experimental

### Infusion of leucine enkephalin

Leucine enkphalin (50 pg/µL in 50/50 MeCN/H<sub>2</sub>O + 0.1% formic acid) was infused at 5 µL/min and measured by electrospray positive ionisation on both the original LCT and the new LCT Premier mass spectrometers. Data was acquired in continuum acquisition mode at a rate of 1 spectrum/second. Figure 1 shows the comparative spectral data obtained from both instruments: the data from LCT Premier has typically 10 times more ion counts than the data from the LCT.





Figure 1. Infusion sensitivity of Leucine Enkephalin.

*Chromatographic signal-to-noise measurement on the LCT Premier*

Reserpine (1 pg/ $\mu$ L, 10  $\mu$ L injection) was repeatedly injected onto a HPLC column. Using electrospray positive ionisation, the ion at  $[M+H]^+ = 609.2812$  was monitored and the signal-to-noise measurement was calculated for each of the peaks. The analytical conditions were:

HPLC: Alliance<sup>®</sup> HT 2795 Separations Module  
 Mobile phase: 75% methanol/25% water + 5 mM ammonium acetate  
 Column: Waters Atlantis<sup>®</sup> dC18, 3  $\mu$ M, 2.1 x 30 mm  
 Injection vol: 10  $\mu$ L  
 Flow rate: 0.3 mL/min




*Figure 2. Signal-to-Noise Measurement of Reserpine (10 pg).*

Figure 2 shows the calculated signal-to-noise measurement for five repeat injections. The results show that the signal-to-noise measurement is typically in the region of 200:1 for 10 pg of reserpine on-column.

In another example of the LCT Premier's full spectral sensitivity, Glu-fibrinopeptide (GFP) has been analysed by nano-LC with detection by NanoLockSpray<sup>™</sup>. With positive ionisation, GFP produces a doubly charged ion at  $m/z$  785.8. Figure 3 shows that with the LCT Premier's sensitivity, 500 attomoles of GFP have been detected with a signal to noise of approx 30:1 (upper trace). The peak just after the GFP peak at 22.12 mins is an interference at the same nominal mass.

The LCT Premier: high sensitivity oa-TOF MS for LC/MS applications.



*Figure 3. 500 attomole of GFP on-column.*

## Sales Offices:

**AUSTRIA** 43 1 877 18 07

**AUSTRALIA** 61 2 9933 1777

**BELGIUM AND LUXEMBOURG** 32 2 726 1000

**BRAZIL** 55 11 5543 7788

**CANADA** 800 252 4752 X2205

**CZECH REPUBLIC** 420 2 617 11384

**DENMARK** 45 46 59 8080

**FINLAND** 358 9 506 4140

**FRANCE** 33 1 3048 7200

**GERMANY** 49 6196 400600

**HONG KONG** 852 29 64 1800

**HUNGARY** 36 1 350 5086

**INDIA** 91 80 2837 1900

**IRELAND** 353 1 448 1500

**ITALY** 39 02 27 4211

**JAPAN** 81 3 3471 7191

**KOREA** 82 2 820 2700

**MEXICO** 52 55 5524 7636

**THE NETHERLANDS** 31 76 508 7200

**NORWAY** 47 6 384 6050

**PEOPLE'S REPUBLIC OF CHINA** 86 10 8451 8918

**POLAND** 48 22 833 4400

**PUERTO RICO** 787 747 8445

**RUSSIA/CIS** 7 095 931 9193

**SINGAPORE** 65 6278 7997

**SPAIN** 34 93 600 9300

**SWEDEN** 46 8 555 11 500

**SWITZERLAND** 41 62 889 2030

**TAIWAN** 886 2 2543 1898

**UK** 44 208 238 6100

**US** 800 252 4752

# Waters

WATERS CORPORATION

34 Maple St.

Milford, MA 01757 U.S.A.

T: 508 478 2000

F: 508 872 1990

[www.waters.com](http://www.waters.com)

## For Complete Confidence

Waters, Alliance, Atlantis, LCT, LCT Premier and NanoLockSpray are trademarks of Waters Corporation.  
All other trademarks are the property of their respective owners.  
©2004 Waters Corporation. Produced in the U.S.A. September 2004 720000933EN LL-PDF

