

Tabisam Khan<sup>1</sup>, Lisa J Calton<sup>1</sup>, Cozette Cuppett<sup>2</sup>, Steve Preece<sup>1</sup><sup>1</sup>Waters Corporation, Floats Road, Wythenshawe, Manchester, M23 9LZ, UK<sup>2</sup>Waters Corporation, 34 Maple St., Milford, MA 01753

## INTRODUCTION

In this work we present a method of analysis for Amitriptyline, a common anti-depressant, in human plasma that has been validated on a HPLC/MS/MS system to the strict guidelines required by the FDA. This analysis was then compared to a method developed for a new UPLC™/MS/MS system that has made the analysis more sensitive and faster. Incorporating ultra performance liquid chromatography technology with the latest technology in tandem quadrupole mass spectrometry has resulted in a faster and more sensitive approach to quantitative bioanalysis.

## AIMS

The goal of these experiments was to develop a quantitative method for Amitriptyline in human plasma using Ultra Performance Liquid Chromatography (UPLC™) coupled to a tandem quadrupole mass spectrometer. UPLC is a novel technique that utilises columns packed with 1.7 µm particles, operated at high linear velocities resulting in increased chromatographic resolution and peak height compared to conventional HPLC. This technique has never previously been applied to a validated quantitative bioanalysis assay.

## EXPERIMENTAL

### SAMPLE EXTRACTION PROCEDURE

1. The method consisted of mixing human plasma (varying amounts depending on the standard used see Table 1) with internal standard (Doxepin at 100 pg/µL, 10 µL).
2. The resulting mixture was vortex mixed (ca 10 seconds) and centrifuged (ca 3000 rpm, ca 3 minutes).
3. The supernatant was placed into a glass vial and injected onto a Waters® Symmetry® C<sub>18</sub> (100 x 2.1 mm) HPLC column.
4. A 10 µL injection of the extract was run by gradient reversed phase chromatography with a mobile phase comprised of water (0.1% formic acid): acetonitrile (0.1% formic acid) at a flow rate of 0.3 mL/minute.

## PREPARATION OF STANDARD SOLUTIONS

### Amitriptyline primary stock solution (1mg/mL) and working standard solutions

A primary stock solution of Amitriptyline was prepared at a concentration of 1 mg/mL by accurately weighing out 5 mg of Amitriptyline into a flask and dissolving in 5 mL of mobile phase.

### Internal standard primary stock solution (100 pg/µL)

Internal standard solution was prepared by diluting a stock solution of Doxepin (1 mg/mL) 1 in 10 with acetonitrile/water (0.1% formic acid), 50/50, v/v.

## LC/MS/MS

LC system: Waters® 1525µ Binary HPLC Pump with a 2777 Sample Manager  
Column: Waters Symmetry® C<sub>18</sub>, 3.5 µm, 2.1 x 100 mm  
Isocratic: Mobile phase A: Water (30%)  
Mobile phase B: Acetonitrile (70%)  
Injection volume: 10 µL

## UPLC/MS/MS

LC system: Waters ACQUITY UPLC™  
Column: Waters ACQUITY BEH C<sub>18</sub>, 1.7 µm, 2.1 x 100 mm  
Gradient: Mobile phase A: Water  
Mobile phase B: Acetonitrile  
Time (min) Flow Rate [mL /min] % A % B  
Initial 0.3 95 5  
0.50 0.3 5 95  
1.50 0.3 5 95  
1.60 0.3 95 0.5  
Injection volume: 5 µL

## MASS SPECTROMETRY CONDITIONS

MS System: Waters Micromass® Quattro Premier™  
MRM Transition Dwell (secs) Cone Volt. Col.Energy  
278.30 > 91.00 0.05 25.0 26.0  
280.50 > 106.90 0.05 22.0 23.0

## RESULTS

### (A). LINEARITY

All calibration lines run over a five-day calibration gave correlation coefficients of greater than  $r^2 > 0.99$  with all standards within +/-15% deviation. Examples of two such calibration lines can be seen in Figure 1. This figure shows that the UPLC/MS/MS calibration line is equivalent to the standard HPLC/MS/MS method differs in the linearity as it covers the range of 0.01 to 100 pg/µL where as the standard HPLC/MS/MS is over the range of 0.5 to 500 pg/µL. This is due to the extra sensitivity that is obtained by running UPLC with the Quattro Premier. The top standard was too concentrated and saturated the detector. Therefore, the top standard was reduced to 100 pg/µL with the whole linear range going down to 0.01 pg/µL and maintaining the four orders of linear magnitude. Clearly, with the UPLC-Quattro Premier we have shown that increased sensitivity can be achieved while maintaining linearity over 4 orders of magnitude.




Figure 1. A comparison of the linearity of the UPLC-Quattro Premier with the standard HPLC-Quattro Premier.

### (B) QUANTIFICATION

The system approach to quantitative bioanalysis incorporates quantification software that fully supports UPLC and the Quattro Premier. The results shown in Figure 2 give excellent back calculated concentrations with deviations from the theoretical concentrations being +/-10%. This shows good correlation to the strict criteria as stipulated in GLP environments.

### (C). SENSITIVITY AND CHROMATOGRAPHIC RESOLUTION

The sensitivity differences have been highlighted on the QC 30pg/µL samples. In Figure 3 the signal to noise ratio is used as a guide to the relative differences in sensitivity that can be achieved. With the standard HPLC-Quattro Premier a relatively good signal to noise ratio of approximately 900:1 is achieved. This was compared to the UPLC-Quattro Premier and the result was 4500:1—an increase of five times in sensitivity with an increase of ten times in overall signal intensity.




Figure 3. A Comparison of QC 30 pg/µL on the UPLC-Quattro Premier with the HPLC-Quattro Premier.

### (D). ROBUSTNESS

For all good validated studies a robust method is required. The UPLC-Quattro Premier method was evaluated by running 200 QC's over a 24-hour period. This experiment was designed to see how robust the whole system and methodology was with complex matrices. The overall percentage variation expressed as a coefficient of variation (CV), was 2.9% as can be seen in Figure 4.




Figure 4. Evaluation of the Robustness for the UPLC-Quattro Premier Validated Method.

## CONCLUSIONS

The calibration lines all gave coefficients of determination of better than 0.995 over a linear range of 0.1 pg/µL to 5000 pg/µL. Back-calculated values for calibration and QC standards were all within 15%, including at the LOQ, thereby meeting the accepted criteria for bioanalytical method validation. No deterioration of column performance was observed for the duration of the analysis. Comparison with the HPLC/MS/MS data showed that, for the same analysis time, the chromatographic resolution was improved.

1.7 µm particle packed columns operated at elevated pressures have been used for the quantitative analysis of Amitriptyline in protein precipitated human plasma. Improved chromatographic resolution and sensitivity were observed compared to an equivalent HPLC/MS/MS method.




Figure 2. Quantification results window.

## PREPARATION OF QUALITY CONTROL SAMPLES

The QC samples were divided into suitable volumes and analysed on preparation. The remaining was stored at -22 °C until used for analysis.

Standard concentration used (ng/µL)	Volume of standard solution used (µL)	Volume of plasma (µL) + 500 µL of acetonitrile	Volume of Doxepin internal standard (µL)	Final concentration (pg/µL)
0.01	30	460	10	0.3 **
0.1	30	460	10	3 **
1	30	460	10	30 *
10	30	460	10	300 *
100	30	460	10	3000 *

Table 2. Detail for the preparation of quality control samples.

\* Used only for the LC/MS/MS calibration line.

\*\* Used only for the UPLC/MS/MS calibration line.