

Paul D. Rainville and Thomas E. Wheat
Waters Corporation, 34 Maple St. Milford, MA

Introduction:

LC-MS is a powerful tool used for the characterization of proteins. Reversed-phase HPLC is the most common mode of chromatography coupled with mass spectrometry. This separation technique can also be an effective method for the desalting of proteins. Pharmaceutical biotherapeutics and other proteins are often prepared or stored in salt containing buffer solutions such as PBS (phosphate buffered saline). These buffer solutions contain salts that form adducts with proteins, complicating data interpretation. These salts also suppress ionization during ESI-MS analysis. Therefore, separation of proteins from non-volatile salts is an important step prior to introduction into the mass spectrometer. Here we show a rapid, on-line desalting procedure with cycle times less than five minutes. The MassPREP™ On-Line Desalting Cartridge is a 2.1 x 10 mm device packed with polymer sorbent and was evaluated using acidic, (bovine serum albumin), basic, (cytochrome c), and large globular, (monoclonal antibody) protein samples. Parameters such as column carry over, loading, and lifetime were investigated. The results of this study show no sample carryover and excellent lifetime for repeated sample analysis under the conditions tested.




Figure 1A and 1B: Desalting of monoclonal IgG. The sample in PBS gives no recognizable protein mass spectrum. After desalting, the sample gives a typical protein spectrum that can be deconvoluted using the Waters® Micromass® MaxENT™ 1 algorithm to reveal isoforms of IgG.

Instrumentation:

Divterter valve:	Waters Selection Valve
HPLC System:	Alliance® 2796 Separations Module
Needle Wash Solution:	H ₂ O (35%) / Isopropanol (5%) / Acetonitrile (60%)
Mass Spectrometer:	Waters® Micromass® Q-tof micro
Ionization mode:	ES +

Divterter valve, HPLC system, and mass spectrometer were all controlled by Waters® MicroMass® Masslynx™ software



Sample preparation:

Bovine serum albumin (BSA), cytochrome c and monoclonal IgG, were obtained from Sigma Aldrich Inc. St. Louis Mo. All proteins were dissolved in phosphate buffered saline (PBS) at a concentration on 1ug/uL.



Separation method:

Eluent A: H₂O with 0.1% Formic Acid
 Eluent B: Acetonitrile with 0.1% Formic Acid
 Flow: 0.4 mL/minute

Time	Eluent A	Eluent B	Valve position
0	95	5	waste
0.5	95	5	waste
2.0	20	80	Mass spectrometer
3.0	95	5	waste


Chromatogram of desalting*Mass spectra of desalted proteins*

Figures 2A and 2B: The chromatographic elution profile is shown for three test proteins in figure 2A. Mass spectra summed across the chromatographic peaks are shown in figure 2B. These experiments show that the cartridge can be used over a wide range of protein types.

Chromatograms of carryover and lifetime study*Summed spectra of chromatograms in figure 3A.*

Figures 3A and 3B: A series of injections consisting of a blank, 1 µg BSA, 1 µg monoclonal IgG, and blank was repeated 25 times for a total of 100 injections. The performance of the cartridge is unaltered over this series. There is no evidence of carryover or background interference in the blank injection representing the 100th injection.

Chromatograms showing increasing injected mass of IgG. Summed spectra of indicated region in figure 4A.



Figures 4A and 4B: Total ion chromatograms and summed spectra for indicated region demonstrating loading capacity and carryover of monoclonal IgG. Blanks were injected between increasing loads of IgG1 to measure carryover (data not shown). The amount of sample can be increased to 10 µg in 10 µL of PBS without breakthrough and with no detectable carryover in the blank following a 10 µg injection.

Conclusion:

Waters MassPREP™ On-Line Desalting Cartridges are an effective tool for the successful desalting of acidic, basic, and large globular proteins. The desalted proteins yielded abundant ESI-MS signal suitable for deconvolution and protein characterization. Studies determining loading capacity and carryover showed that up to 10 µg of monoclonal IgG was successfully loaded with no breakthrough and no detected carryover in a subsequent analysis. The same results were obtained for BSA and cytochrome c for injected masses of up to 5 µg. Further, no loss of performance was observed after 100 consecutive injections were performed. The overall MassPREP™ desalting strategy provides an easy, fast (<5 minutes), and reproducible approach for the effective desalting of proteins stored in physiological salt buffers prior to high resolution structural analysis.

Sales Offices:

**Austria and European Export (Central South
Eastern Europe, CIS and Middle East)**

(43) 1 8771807

Australia (61) 2 9933 1777

Belgium (32) 2 7261000

Brazil (55) 11 5543 7788

Canada 800 252 4752

China (8610) 84518918

CIS/Russia (7) 095 336 7000

Czech Republic (42) 02 6171 1384

Denmark (45) 46 598080

Finland (358) 9 506 4140

France (33) 1 30 48 72 00

Germany (49) 6196 40 06 00

Hong Kong (852) 2964 1800

Hungary (36) 1 350 5086

India and India Subcontinent

(91) 80 28371900

Ireland (353) 1 448 1500

Italy (39) 02 27421.1

Japan (81) 3 3471 7191

Korea (82) 2 820 2700

Mexico (5255) 5524 76 36

The Netherlands 31 (0)76-50 87 200

Norway (47) 63 84 60 50

Poland (48) 22 833 4400

Puerto Rico (787) 747 8445

Singapore (65) 6278 7997

Spain (34) 93 600 93 00

Sweden (46) 8 555 11500

Switzerland (41) 62 889 2030

Taiwan (886) 2 2543 1898

United Kingdom (44) 208 238 6100

United States: 800 252 4752

Waters

WATERS CORPORATION

34 Maple St.

Milford, MA 01757 U.S.A.

T: 508 478 2000

F: 508 872 1990

www.waters.com

For Complete Confidence

Waters, Micromass, MassLynx, MassPREP, MassENT, and Alliance are trademarks of Waters Corporation.

720001077EN

