


Larry Zintek, Sr. Chemist, and Josh Neukom, Chemist, USEPA Region 5 Central Regional Laboratory, Chicago, IL 60605  
Jim Krol, Sr. Applications Chemist, Waters Corporation, Milford, MA 01757

## Homeland Security *Dealing with our Drinking Water*

"The focus is on catastrophic terrorism—threats to the security of our homeland that could result in large-scale loss of life and major economic impact."

[http://www.dhs.gov/dhspublic/theme\\_home5.jsp](http://www.dhs.gov/dhspublic/theme_home5.jsp)



"Homeland security is a current issue and deliberate contamination, terrorist attack, of the drinking water supply is a major concern."

- Homeland Security Presidential Directive / HSPD-9

## Initial EPA Region 5 60 Deleterious Analyte List

### Acidic Analytes:

2,4-D	Cacodylic Acid	Dalapon	Endothall
Fenamiphos	Glyphosate	Methidathion	Picloram
		Silvex	

### Neutral Analytes:

1-Naphthol	3-Hydroxycarbofuran	Alachlor	Aldicarb
Aldicarb Sulfone	Aldicarb Sulfoxide	Azinphos-Methyl	Bromadiolone
Buprofenzin	Carbaryl	Carbendazim	Carbofuran
Chlorpyrifos-methyl	Coumaphos	Cyclohexamide	Dichlorvos
Dimethoate	Diuron	EPN	Ethoprophos
Fenhexamid	Fensulfothion	Methamidophos	Methiocarb
Methomyl	Mevinphos	Oxamyl	Propoxur
Pyridaben	Pyridaphenthion	Scilliroside	Tebufoenpyrad
Tebufoos	Thiobencarb		

### Basic Analytes:


2-Aminobenzimidazole	Cyanizine	Cyprodinil	Dicrotophos
Imazalil	Imidacloprid	Nicotine	Pendimethalin
Pirimiphos-methyl	Propamocarb	Strychnine	Tebuconazol
Terbumeton	Terbutylazine	Thiabendazole	Toxaphene
Chloropicrin			

## Homeland Security Directive HSPD-9 EPA Initiative

- EPA Region 5 CRL is undertaking an initiative to develop a robust, comprehensive, and fully coordinated surveillance and monitoring system for drinking water quality that provides early detection and awareness of poisonous (deleterious) agents.
- EPA Region 5 CRL and Waters Corporation entered into a Cooperative Research and Development Agreement (CRADA) to develop a HPLC method utilizing Mass Spectrometry using MS/MS Spectral Library Confirmation.
  - No transferable LC/MS Libraries are available for wide-spread use
- LC/MS method to provide laboratories with a screening tool to quickly identify deleterious organic contaminants in drinking water, and final quantification using LC/MS/MS to confirm identification.
  - 280 analytes are listed into 4 phases based on toxicity;
    - Commonly used agrochemicals, pesticides and insecticides
  - Drugs of abuse, nicotine, LSD, methamphetamine, ecstasy

## MS Criteria for Library Spectra

- The cone voltage for a single quadrupole system was determined as giving a precursor ion  $[M+H]^+$  response at 10-25% response of the major product ion, and, if possible, yield 2 or more product ions.
- MS/CID Analyte full scan Spectra acquired during a chromatographic run of 1 ppm analyte standard in mobile phase.
- MS/CID Full Scan Spectra and relative retention time entered into library entry for 2 variable analyte confirmation.
- Results for single quadrupole yield equivalent results, and should be considered as a lower cost system alternative.
  - Single quadrupole full scan analysis is however less sensitive than product ion MS/MS



## The MS/CID and MS/MS Dual Strategy

- LC/MS/MS methods have been established for various multi-analyte subsets
  - MS/MS specificity and sensitivity in complex matrices
  - Detection of non-UV active analytes
- Subset Analytes into Acidic, Neutral, and Basic Groupings
  - Acidic pH for acidic analytes with -ESP
  - Neutral pH for neutral analytes with  $\pm$ ESP
  - Basic pH for basic analytes with +ESP
  - Minimize other ionization techniques, such as APCI, APPI, etc
- Analyte identification based on chromatographic retention time and MS/MS ESP Product Spectra, i.e., comparison of unknowns against a full scan analyte MS/MS product spectra for each analyte
  - Similar to Photodiode Array spectral processing
  - Analyte confirmation using conventional MS/MS transitions possible
- LC/MS Screening method
  - Create full scan Libraries with Single Quadrupole MS/CID
  - Establish cone voltage settings for analysis, and processing method parameters
- Analysis time less than 25 minutes with little sample preparation

## MS/MS Criteria for Library Spectra

- The cone voltage for MS1 was determined at greatest  $[M+H]^+$  or  $[M+NH_4]^+$  response using a 1 ppm pure standard prepared in mobile phase.
- The collision energy was determined giving a precursor ion  $[M+H]^+$  response at 10-25% response of the major product ion, and, if possible, yield 2 or more product ions.
- MS/MS Analyte Spectra acquired during a chromatographic run of 1 ppm analyte standard in mobile phase using daughter ion MS/MS functionality.
- MS/MS Product Ion Spectra and relative retention time entered into library entry for 2 variable analyte confirmation.

## Multi-Analyte LC/MS/MS Method Neutral pH "Universal Gradient"

System\*: Waters Alliance® 2695 System with Quattro micro™ API Triple Quadrupole, or with ZQ™ Single Quadrupole

Column: Waters Atlantis™ C<sub>18</sub>, 2.1 x 150 mm, 3.5  $\mu$ m

Buffer A: 5% AcCN / 5 mM NH<sub>4</sub>OAc, natural pH

Solvent B: 95% AcCN / 5 mM NH<sub>4</sub>OAc

Flow: 300  $\mu$ L / min

Col Temp: 30°C

Gradient: Hold at 100% A for 1 min  
Linear Grad to 100% B over 20 min  
Hold at 100% B for 5 min  
Re-equilibrate for 10 min  
100  $\mu$ L depending on LOQ requirements



Inj Vol: 10  $\mu$ L

Ionization: ESP


Data Proc: Waters MassLynx™, ver 4, with Chromalynx, QuanLynx, TargetLynx Options

\*Mention of Vendor Brands by US EPA Does Not Constitute Product Endorsement

## Examples of Library Product Ion Spectra



## Low Level 1 ppb Fensulfothion in Drinking Water



## Library Evaluation Summary

- The created ESP MS/MS and MS/CID Library have been evaluated by 2 outside laboratories using the same brand of triple quadrupole instrument.
- MS/MS (Product Ion) Spectra using Triple Quadrupole
  - All of the analytes can be successfully detected at 100 ppb
  - All of the analytes can be successfully detected at 25 ppb
  - 10 of 29 can be successfully detected at 1 ppb
- MS/CID Spectra Library
  - 25 out of 29 analytes can be successfully detected at 100 ppb
  - Detection limit is compound dependent and varies
  - 1 ppb is below detection limit for all MS/CID Library work
  - Analyte co-elution will bias the MS/CID Spectra because of additional co-elution ion fragments.
    - To be addressed using Chromalynx deconvolution spectral processing successfully used for GC/MS.
- Improvements and additions to the libraries are continually being made as more method research is conducted.
  - Adding additional analytes; acidic and basic analytes
    - Goal to have 200-300 analytes of homeland security concern
  - Adapt method protocol and libraries to single quadrupole instruments using in-source CID
    - Lower cost alternative for field labs
- Improvement in detection sensitivity and spectral quality
  - Analyte concentrations in drinking water are dependent upon the action limits defined as threat agent concentration.
  - Analyte enrichment using SPE may be necessary