

Weibin Chen<sup>1</sup>, Grace M. Credo<sup>1</sup>, Henry Shion<sup>2</sup>, John C. Gebler<sup>1</sup><sup>1</sup>Waters Corporation, Milford, MA 01757, United States, <sup>2</sup>Waters Corporation, Beverly, MA 01915, United States

## Overview

- A patterned self-assembled monolayer (SAM) that confers zones of differing wettability was constructed on a gold-coated silicon wafer to serve as a sample platform for MALDI TOF analysis.
- The usage of this substrate greatly increased the sensitivity of peptide analysis, resulting in 10 to 100-fold improvements compared to traditional MALDI TOF MS on a stainless steel plate.

## Introduction

The wide application of MALDI TOF MS in proteomics has become a driving force to improve the technology towards higher sensitivity and better reproducibility. Among many other factors, sample preparation is recognized as being a critical step in the overall workflow for high sensitivity analysis. In this poster, we present our approach to improve the sensitivity of MALDI TOF MS analysis via a novel MALDI substrate.

The substrate consists of a silicon wafer with a patterned self-assembled monolayer (SAM) surface that confers zones of differing wettability (see **Figure 1**). The molecularly-flat surface allowed for concentration and precise placement of samples for MALDI mass spectrometry. The variation of wettability permits a large volume of droplet to be applied. As the sample dries, it is focused into a relatively small analysis zone, thus dramatically improving the sensitivity of analysis. Evaluations of the substrate for reproducible, high sensitivity MALDI MS of peptides and protein digests are presented.

## Experimental

### 1. Substrate


The substrate was manufactured and provided by LumiCyte Inc. (Fremont, CA)

### 2. Sample Preparation Procedure


- Prepared sample solutions in aqueous 25% acetonitrile, 0.1% TFA (v/v)
- Deposited up to 5.0  $\mu$ L of sample (**Figure 2**) and waited for sample to dry completely
- Added 1  $\mu$ L matrix solution, CHCA (Waters) 0.07 mg/mL in solvent mixture ACN:EtOH:0.1%TFA(aq), 80:10:10 (v/v/v)
- Waited for matrix to dry and analyzed by MALDI TOF MS (Waters Micromass<sup>®</sup> MALDI micro MX<sup>™</sup>)

## Schematic Illustrations of the Substrate & Workflow in Sample Preparation

**Figure 1.** Sample well composed of differing wettability zones on the substrate surface. Order of wettability: Analysis Zone > Liquid Retention Zone > Boundary Zone




**Figure 2.** Sample Preparation Workflow using the Substrate




## Results


**Figure 3.** MALDI MS spectra of a peptide standard, [Glu<sup>1</sup>]-Fibrinopeptide, deposited on a high sensitivity SAM-based substrate exhibiting 0.01 fmol (10 amol) sensitivity. In a typical analysis, this sensitivity is achieved at least 10 of 12 wells.




**Figure 4.** MALDI MS spectra of a BSA tryptic digest deposited on a high sensitivity SAM-based substrate exhibiting 50 amol sensitivity. Spectra were subjected to peptide mass fingerprinting for positive protein ID.



**Figure 5.** Comparison of the sensitivity of a SAM-based substrate and a stainless steel target for analysis of diluted in-gel tryptic digests. In-gel tryptic digests from a SDS-PAGE band (1 pmol protein loaded) were recovered with 25  $\mu$ L of 25% acetonitrile. The digest extract was diluted sequentially with 25% acetonitrile/0.1% TFA solution at 10-, 100-, 1000-fold ratios. Only 1  $\mu$ L of each solution was applied to each target well.



**Figure 6.** Comparison of target plate sensitivity for the analysis of a blind protein digest sample, with a closer view of a selected peptide. Data collection was in automated mode.



## Conclusions

- A novel SAM-based MALDI substrate has been successfully constructed on a silicon wafer for high sensitivity peptide analysis.
- A sample preparation method has been developed for the usage of the target.
- An automated sample collection method has been developed for reproducible high sensitivity data collection on the target.
- The SAM-based substrate yields 10–100 fold improvements in sensitivity compared to a conventional MALDI stainless steel target.