

OVERVIEW

- Three off-line 2D-LC setups were compared including SCX-RP and two alternative systems: high/low pH RP-RP, and HILIC-RP.
- RP-RP (high/low pH) 2D-LC system provides the greatest number of identified peptides/proteins.
- Two methods of tandem mass spectrometry identification of peptides/proteins were compared: DDA, and a novel alternate scanning MS^E

INTRODUCTION

LC-MS methods are used for proteome research for disease markers and/or drug targets discovery. Due to the complexity of proteome samples and the wide protein dynamic range the analysis is difficult. An efficient LC separation and MS/MS analysis is required to enable the detection of lower-abundant proteins of interest in mixtures. One way to address the complexity of these samples is the application of multidimensional separation techniques like 2D-HPLC.

Currently the most common 2D-HPLC approach for separation of peptides is a combination of strong cation-exchange (SCX) and reversed-phase (RP) HPLC.

In this work we applied two alternative strategies for 2D-HPLC separation. First uses a combination of two RP-RP modes (using different pH), the second employs hydrophilic interaction chromatography (HILIC) and RP-LC. All 2D-HPLC methods were run in off-line mode using a 17 protein digests containing ~1500 peptides.

A novel alternate scanning mode for MS/MS using alternate scans at low and elevated energy has been used and compared to a traditional DDA methods.

EXPERIMENTAL

Sample Preparation

The standard proteins (Sigma) were reduced, alkylated and digested with modified porcine Trypsin (Promega) according to a standard protocol and then mixed in dynamic range 0.03-10 pmol/uL. ~8 Fractions (5 minutes) were collected in each mode (SCX, HILIC, RP high pH), volume was reduced by evaporation to the volume of original sample injected on 1stD LC (no concentration prior to 2ndD LC) and analyzed by LC-MS/MS.

HPLC Columns for fractionation

PolySULFOETHYL Aspartamide SCX 200Å, 5um 1x150mm (PolyLC Inc.)

Atlantis HILIC, 5um 1x150mm (Waters)

RP BEH-C18, 3.5um 1x150mm (Waters)

LC-MS Conditions

Tryptic peptides were analyzed by LC-MS/MS on Q-ToF Ultima (Waters). ProteinLynx Global Server V2.2 software was used for the data analysis.

RP column: NanoEase Atlantis dC18, 3um 0.3x150mm (Waters)

RP flow: 5 uL/min

Solvent A: 0.1% Formic Acid in water

Solvent B: 0.1% Formic Acid in 80% ACN

Temperature: 40°C

Gradient: 0-45 min 0-56% B; 1% ACN/min

RESULTS AND DISCUSSION

A. SCX mode, pH 2.6




Figure 2: LC-MS analysis of fractions

B. HILIC mode, pH 4.5




Figure 3: LC-MS analysis of fractions

C. RP high pH mode, pH 10.0




Figure 4: LC-MS analysis of fractions

Protein Expression Analysis vs. DDA - MS/MS acquisition




Figure 5: Data Directed Analysis (DDA). DDA employs an MS survey scan to identify abundant components (precursors) and then select these individual precursors for MS/MS fragmentation. The survey/analysis cycle is repeated throughout the run.




Figure 6: Protein Expression Analysis (MS^E). PE Analysis utilizes a fast cycle a low and elevated potential to collect parent and child ions. The child and parent ions are binned together via a accurate mass/intensity/retention time information. Data are then deconvoluted via software. This mode of data acquisition has an increased duty cycle over traditional DDA methods as well as a typical 5-50 fold increase in sensitivity.

Protein Expression Results

Proteins	pmol/uL	max. number of peptides *	1D LC-MS Number of peptides	SCX	HILIC	RP high pH
b-Galactosidase (E.coli)	0.03	31	-	-	-	4
b-2 Microglobulin (human)	0.05	7	-	-	-	-
a-2 Macroglobulin (human)	0.1	62	28	13	18	26
b-Lactoglobulin A	0.1	9	5	1	-	1
Ovalbumin	0.3	18	-	-	1	2
Alkoholdehydrogenase (yeast)	0.5	30	16	12	11	16
Enolase (yeast)	0.5	30	18	11	15	15
Hemoglobin (bovine)	0.5	24	8	4	4	10
Phosphorylase b (rabbit)	0.5	82	35	17	15	12
Peroxidase (horseradish)	0.5	11	-	-	-	-
Glyceraldehyd-3-phosphate dehydrogenase (rabbit)	0.5	24	3	3	4	7
Carboxypeptidase A (bovine)	0.75	12	6	-	4	6
Catalase (bovine liver)	1.5	36	27	18	23	25
Ubiquitin (bovine)	2	8	4	3	1	2
Transferrin (bovine)	2	54	38	24	27	38
a-1 Acid glycoprotein (bovine)	7.5	14	11	7	2	8
BSA	10	60	49	33	39	49
Total: 512 248 146 164 221						

Table 1: Off-line 2D LC-MS^E results for SCX, HILIC and RP high pH analysis of 17 protein digest (~1500 peptides).

* Maximum number of peptides represent peptides identified in a separate analysis of 2 pmol protein digest injection.

CONCLUSIONS

- A novel methods for 2D-HPLC separation of proteomic samples have been developed.
- A proposed high/low pH RP-RP 2D-HPLC approach allowed for a significantly increased number of identified peptides compared to HILIC-RP and SCX-RP modes.
- The advantages of RP-RP-HPLC are high peak capacity, high peptide recovery, and the use of salt free mobile phases.
- The performance of high/low pH RP-RP-HPLC compares well to current state-of-the-art SCX-RP-HPLC results and it could be used as an alternative to SCX-RP-HPLC.
- Novel MS/MS^E method doubled the number of identified peptides and provided a significantly better confidence for protein identification in complex sample compared to a traditional DDA mode.
- Novel MS^E method represents a promising tool for proteomic research.

DDA Results

Proteins	pmol/uL	1D LC-MS Number of peptides	SCX	HILIC	RP high pH
b-Galactosidase (E.coli)	0.03	-	1	-	-
b-2 Microglobulin (human)	0.05	-	-	-	-
a-2 Macroglobulin (human)	0.1	5	4	6	9
b-Lactoglobulin A	0.1	1	3	1	3
Ovalbumin	0.3	1	1	1	3
Alkoholdehydrogenase (yeast)	0.5	9	10	6	12
Enolase (yeast)	0.5	9	6	5	8
Hemoglobin (bovine)	0.5	2	4	1	4
Phosphorylase b (rabbit)	0.5	10	6	6	12
Peroxidase (horseradish)	0.5	1	1	1	2
Glyceraldehyd-3-phosphate dehydrogenase (rabbit)	0.5	-	2	1	2
Carboxypeptidase A (bovine)	0.75	-	1	-	2
Catalase (bovine liver)	1.5	15	16	17	18
Ubiquitin (bovine)	2	1	2	1	3
Transferrin (bovine)	2	21	13	13	24
a-1 Acid glycoprotein (bovine)	7.5	7	5	4	8
BSA	10	43	21	24	36
Total: 125 96 87 146					

Table 2: Off-line 2D LC-MS/MS (Data Directed Analysis) results for SCX, HILIC and RP high pH modes.