

Ignatius J. Kass¹, Hongji Liu², Kevin M. Millea², John C. Gebler², and Scott J. Berger²
¹Waters Corporation, Beverly, MA ²Life Sciences R&D, Waters Corporation, Milford, MA 01757

OVERVIEW

- Analysis of proteins by electrospray ionization mass spectrometry is complicated by the need to deconvolute multiply charged m/z spectra to produce neutral masses for sample components.
- This deconvolution process becomes significantly more complicated when analyzing data from online LC/MS protein separations.
- TIC information can be used to identify regions of MS data for summation and deconvolution with well resolved protein mixtures, but the approach fails with more complex protein separations.
- In this poster we describe a methodology for automating mass spectral deconvolution of complicated LC/MS data sets.
- The resulting data is output in a tab-delimited text file, and can be represented with intuitive and visually informative displays.
- Proper selection of processing parameters permit accurate intact protein mass determination, and retention of chromatographic profiles of each component.


METHODOLOGY

Yeast ribosomal proteins were analyzed by LC/ESI-TOF MS (Waters LCT) as described in work by Liu and coworkers (Ref. 1). LC/MS data was processed using a prototype software program called AutoME (Automated Maximum Entropy) within MassLynx 4.0 software (Waters) using the following conditions. 10 MS scan segments (1 scan/sec) were combined throughout the entire run. Segments containing a total ion intensity over 7000 counts were subjected to processing by the MaxEnt 1 (Ref. 2) spectral deconvolution algorithm (0.75 Da peak width, 3,000 to 45,000 output mass, 1 Da output resolution) until model convergence or 16 iterations was achieved. The resultant deconvoluted spectra was centroided, and components with intensities of greater than 20 counts were recorded to the output data file. Microsoft Excel 2003 was used to generate 2-D Gel bubble plot and Mass Fingerprint displays.







E.coli cytosolic proteins were analyzed by LC/ESI-TOF MS (Waters LCT Premier) under similar conditions, but with a gradient of only 60 min. Spectra (ESI+, V-mode (6,000 resolution), 1Hz, 650-2990 m/z) were processed by the AutoME program using the following conditions: 10 scans/ processing segment, input segment TIC threshold (10,000 counts), MaxEnt1 (0.75 Da width, 3,500-50,000 output mass range, 1 Da resolution, 15 iterations), Data reduction (Harmonic removal with 1 Da tolerance, 100 count minimum component intensity, Top30 peaks recorded/processing segment).

REFERENCES

- 1) Liu, H., Berger, S.J., Chakraborty, A.B., et al. (2002) J. Chromatogr B, 782: 267-289.
- 2) Ferrige, A.G., Seddon, M.J., et al. (1992) Rapid Commun Mass Spectrom 6: 707. and Ferrige, A.G., Seddon, M.J., et al. (1992) Rapid Commun Mass Spectrom 6: 765.



Component Level ID of Ribosomal Proteins



TOP: LC/MS Protein Map of E. coli Cytosol
LOWER RIGHT: Enlarged image of central map region
LOWER LEFT: Protein Mass Fingerprint (Deconvoluted Mass vs. Intensity) over the entire LC/MS analysis

CONCLUSIONS

- Intact protein LC/MS data sets are amenable to automated data processing and analysis approaches.
- Proper selection of data processing parameters permits intact protein mass data to be measured with high precision and accuracy.
- Processing parameters can be selected such that the underlying chromatographic profiles of individual components are retained or completely ignored.
- Processed data can be generated to compare protein profiles, component intensities, or whole analysis "fingerprints" between samples.
- Intact protein LC/MS analysis can be conducted on a "proteomic" samples with high complexity and a large dynamic range expression.