

OVERVIEW

- The feasibility of monitoring therapeutic proteins in serum samples over time is demonstrated
- A Linear response for alpha galactosidase A, between 10 and 500 fmol injected on a 75µM analytical column is presented
- The data from patient samples illustrates the potential for monitoring global protein concentration changes across multiple conditions

INTRODUCTION

Fabry disease is an X-linked recessive glycolipid storage disorder, caused by deficient activity of the lysosomal enzyme alpha-galactosidase A (alpha gal A). In affected patients progressive accumulation of the glycolipid substrate for this enzyme, globotriaosylceramide, occurs within vulnerable cells and tissues. Patients may show ocular deposits, febrile episodes, burning pain in the extremities and skin lesions (angiokeratomas). Administration of recombinant alpha-gal A has been shown to alleviate symptoms of the disease and is used as a long-term enzyme replacement therapy in patients with a confirmed diagnosis of Fabry's disease. Monitoring circulating levels of alpha gal A can determine the uptake and metabolism of the replacement enzyme and could potentially allow the therapeutic level of the drug to be determined in patient samples, allowing the dose of the drug to be adjusted to maintain safe, cost effective levels of the treatment. Here we detail a preliminary investigation into the potential of LC-MS for analyzing digested samples of human serum, containing spiked levels of recombinant alpha gal A. We have determined the limit of detection of the enzyme, present in the complex serum background, and have established the linear dynamic range over which measurements can be made. Finally, we have investigated the potential for relative quantification of alpha gal A from patient serum samples.

METHODS

Sample preparation.

The human serum samples and alpha gal A were diluted, and solubilized by incubation at 80°C for 15 mins in 0.1% RapiGest™ SF (Waters, Milford, MA) before reduction with 5 mM dithiothreitol and alkylation with 10 mM iodoacetamide. The proteins were then digested with 1% (w/w) sequence grade trypsin overnight (16 hr). RapiGest SF was cleaved by the addition of HCl, followed by centrifugation, and the supernatant subsequently collected. Samples were diluted with 0.1% formic acid to an appropriate concentration prior to analysis.

Experimental.

Experiments were conducted using a reversed phase gradient (5 to 40% acetonitrile in 120 minutes) on a Waters® CapLC™ System (Waters, Milford, MA) coupled to a Q-ToF™ mass spectrometer (Waters Micromass Manchester UK). Standards and samples were run in duplicate and triplicate, respectively. The Q-ToF mass spectrometer was programmed to acquire data as describes previously (1), where the energy applied to the collision cell is alternated between a low (10 eV) and elevated (23-33 eV) energy. Protein identifications and quantitative information were extracted by the use of specialized algorithms, and searching a Human species-specific database.

RESULTS

Initial feasibility experiments were conducted to determine the limits of detection for alpha gal A in the presence of a complex matrix, i.e. digested human serum.




Figure 1. Identification of alpha gal A from a species-specific database, spiked into human serum at the 100fmol level. The search uses the elevated energy fragment ion data.

Relative quantitation at the protein level

A tryptic digest of alpha gal A was spiked at the 10 to 500 fmol level in 0.2 µg/µl of serum. The result of these experiments which identify the proteins based on the high-energy fragment ion data, shown in Figure 1, prior to performing any quantification, is presented in Figure 2.

Accession	Gene Symbol	Description	100
-----------	-------------	-------------	---