

## INTRODUCTION

- One major application of the tandem LC/MS analysis is quantification. An optimum quantification protocol requires a fast LC separation with sufficient MS analysis in order to obtain good linearity and low limit of detection.
- Major challenges in achieving an ultra-sensitive quantification protocol include matrix interference, limitations of LC/MS system such as column efficiency, system volume, system back pressure etc., as well as the data collection speed and the stability of the MS.
- Improvement in LC/MS/MS quantification relies on either to improve LC separation (column efficiency, peak capacity, system volume etc) via increasing the separation efficiency, improving column chemistry as well as reducing column dimension and/or particle size, or to improve the data acquisition rates in the mass spectrometer.

### Overview

- LC-MS/MS for 5 drug mixture
  - Comparison of UPLC/MS/MS vs. HPLC/MS/MS
- Simultaneous Polarity switching
  - Analysis by both ESI+ and ESI- in a single injection
- Quantification performed in crashed Rat Plasma
  - Complete quantification limits defined for both UPLC/MS/MS and HPLC/MS/MS
- Complete method evaluation for UPLC/MS/MS
  - System robustness and carryover in crashed rat plasma


### Experimental Conditions

#### UPLC Conditions

- Waters ACQUITY Ultra Performance LC™ System (UPLC™)
- Waters ACQUITY Ultra Performance LC™ Column:
  - 2.1 x 50 mm, 1.7  $\mu$ m, 40°C
  - The Bridged Ethane-Silicon C<sub>18</sub> particles (BEHC<sub>18</sub>)
- Injection Volume: 5  $\mu$ L
- Flow Rate: 0.6 mL/min
  - Back Pressure at about 11,000 psi
  - No splitting for MS detection
- Mobile Phase:
  - 10 mM NH<sub>4</sub>OAc at pH 5.0 in 20% MeOH (A)
  - 10 mM NH<sub>4</sub>OAc at pH 5.0 in MeOH/Acetonitrile at 80/20 (B)
- Gradient:
  - 45% B to 50% in 0.4 minute
  - Hold at 95% B for another 0.1 minute
  - Total run time was 2 minute


#### MS Conditions

- Desolvation Temperature: 400°C
- Desolvation Gas Flow: 800 L/Hr
- Capillary Voltage: 0.5 kV
- Source Temperature: 130°C



## QUANTIFICATION RESULTS


### UPLC/MS/MS HPLC MRM Chromatograms



### Simultaneous Polarity Switching

- For UPLC/MS/MS, the total cycle time was 110 ms
  - ESI+, dwell time was 5 ms each for all ions
  - ESI-, dwell time was 20 ms for ibuprofen, and 30 ms for naproxen
  - Inter-channel delay was 5 ms
  - Inter-scan delay was 20 ms
- Both LC/MS/MS detections were performed with simultaneous polarity switching. As a result, analytical results of both ESI+ and ESI- were obtained with a single injection

### UPLC Calibration Curves 3 Injections per Concentration




- Standards were spiked into crashed rat plasma
  - Rat plasma was crashed with acetone at 2/1 ratio (2 AcN/1 rat plasma)
  - After vortex, centrifuged at 12000 rpm for 15 minutes in 4°C
  - Supernatant was diluted with water in 1:1 ratio
- The injection sequence was the following
  - Rat plasma blank (3 injections)
  - Standard rat plasma solution low to high concentration (0.1, 0.2, 0.5, 1.0, 2.0, 5.0, 10.0, 20.0, 50.0, 100, 200, 500, 1000 ng/mL)
  - Rat plasma blank (3 injections)
  - Rat plasma blank (3 injections)
  - Mid QC (6 injections)
  - Rat plasma (3 injections)
  - High QC (6 injections)
  - Rat plasma


### UPLC Quantification Limits

	Alprazolam	Diphenhydramine	Prednisolone	Ibuprofen	Naproxen
Linear Range ng/mL	0.1 - 1000	0.1 - 200	2 - 1000	0.1 - 1000	5 - 1000
R2	0.998	0.998	0.996	0.995	0.999
Slope	0.109	5249.22	10.68	7.04	2.21
LOD	0.0136	0.025	1.642	0.0847	0.541
LOQ	0.0828	0.150	9.852	0.508	3.246
Low QC RSD%	3.42	8.67	-----	7.92	7.02
Med QC 20 ng/mL RSD%	1.39	5.63	10.21	6.30	6.95
High QC 200 ng/mL RSD%	1.20	5.44	5.92	3.74	2.99
Bias%	-0.5	-4.71	-----	0.315	2.735
Med QC 20 ng/mL Bias%	-9.5	-21.5	-----	-----	-----
High QC 200 ng/mL Bias%	-5.35	-15.9	-8.30	-14.7	-3.15

### HPLC/MS/MS HPLC MRM Chromatograms




### HPLC Calibration Curves 3 Injections per Concentration



## THE COMPARATIVE ADVANTAGES OF UPLC™/MS/MS

### Higher Speed

- Faster separation was obtained by the ACQUITY UPLC™
  - High flow rate since the UPLC™ system can withstand up to 15,000 psi
  - Smaller system volume (150  $\mu$ L)
  - LC column with smaller particle size (1.7  $\mu$ m) plus innovated chemistry (the bridged ethane hybrid phase)
- The narrow UPLC™ peak demands shorter cycle time in MS/MS detection
  - Proper interpretation for quantification purpose requires minimum of 15 data point across the cycle time
  - The cycle time of an MS detection is the summation of the dwell time for all ions plus the inter-scan channel delays as well as the inter-scan delays
  - The shorter the cycle time, the more data points can be taken in terms of stability, and sensitivity
- Overlapping here on the left is the effect of the MS dwell time. The UPLC™ peak is 8x wide, if 100ms dwell time were used for the detection, only 7 data points were obtained across the peak, insufficient for quantification. However, if 5 ms dwell time was used, 60 data points were collected across the peak, for more sufficient quantification.
- The signal intensities were unaffected by the change of the dwell time.




### Increased Sensitivity

- Showing on the right is the comparison of the detected signals from UPLC/MS/MS and HPLC/MS/MS
  - With the same column loading (10 ng/mL at 5  $\mu$ L), a significant signal increase was observed in UPLC/MS/MS vs. the HPLC/MS/MS
- The chromatograms with blue labels were the LC results for diphenhydramine
  - ESI+ detection
  - Compared with HPLC, the signal increased 3.76 times was observed in the UPLC™
- The chromatograms with red labels were the LC results for ibuprofen
  - ESI- detection
  - Compared with HPLC, the signal increased 9.97 times was observed in the UPLC™
- Possible contributing factors for the signal increase in UPLC™ include
  - Narrower LC peak width due to the reduced system volume of the column (1.7  $\mu$ m for UPLC vs. 3.5  $\mu$ m for HPLC)
  - Shorter retention resulting in less sample dilution and dispersion during the concentrated in UPLC compared with HPLC
  - Concentrated in UPLC, the MS detector were more concentrated in UPLC compared with HPLC
- The scale of signal increase was compound dependent

### System Robustness

- The system robustness experiment was performed over the period of 24 hours with the following injection sequence
  - Rat plasma blank
  - Standards rat plasma solutions with concentrations low to high (1.0, 10.0, 100, and 1000 ng/mL)
  - Rat plasma blank
  - QC sequence (200 ng/mL) was injected after the above sequence was run 3 times.
  - Rat plasma blank
  - Standards rat plasma solution with concentrations low to high (1.0, 10.0, 100, and 1000 ng/mL)
  - Rat plasma blank
  - The sequence was also run 3 times after the 300 QC injection was completed
- The calibration curve shown on the figure consists 5 points at each concentration level
  - 3 points before the 300 QC injections
  - 3 points after the 300 QC injections
  - The signal was distinguishable before and after the QC injection indicating the stability of the system
- The two chromatograms shown in the figure were the very first QC injection on the top and the very last QC injection on the bottom
  - Consistent peak shape and peak area indicating the stability of the system



### Negligible Carryover

- The sequence of injections were as the following
  - Plasma blank
  - Standard rat plasma solution injections with concentration from low to high (the concentration sequence was 0.1, 0.2, 0.5, 1.0, 2.0, 5.0, 10, 20, 50, 100, 200, 500, and 1000 ng/mL)
  - Plasma blank
  - Showing here in the figure are from two consecutive injections
    - The standard solution at the highest concentration (1000 ng/mL)
    - The rat plasma blank injected immediately after the injection of the highest standard solution
  - No detectable carryover from rat plasma was observed for all analytes of interest

## CONCLUSIONS

- The system robustness for UPLC/MS/MS was demonstrated
  - For 300 injections of QC sample in crashed rat plasma over 24 hours, the retention time RSD% was 0.21%, the detected peak area RSD% was 4.61% for an externally calibrated compound
- The carryover for UPLC/MS/MS was studied
  - For the 5 compounds tested in this study, the carryover was negligible
  - The scale of sensitivity improvement was compound dependent