

Keith R. Compson, Emmanuelle Claude¹, Marten F. Snel¹, Richard Tyldesley¹, Thérèse McKenna¹, James Langridge¹, John Gebler²

¹Waters Corporation, Manchester, United Kingdom, ²Waters Corporation, Milford, MA, United States.

OVERVIEW

- This poster shows an evaluation of the use of novel MS technology for the identification of phosphopeptides
- Enhanced duty cycle (EDC) was used to increase the detection efficiency of low mass diagnostic fragment ions on a Q-ToF mass spectrometer.
- Detection of PO_3^- was evaluated as an indication of peptide phosphorylation.

INTRODUCTION

The analysis of phosphorylated peptides presents a unique challenge. It is well accepted now that phosphopeptides, often present at only low copy numbers per cell, ionize less efficiently than their non-phosphorylated counterparts, due to negative charge interference. This may result in poor detection of phosphopeptides relative to the 'background' of un-phosphorylated peptides. An effective method for improving the detection of phosphopeptides is the use of an affinity enrichment protocol, producing enriched samples e.g. using IMAC, TiO_2 etc. An example of data obtained from a sample purified using a TiO_2 protocol is shown in **Figure 1**. Following enrichment almost all un-phosphorylated species have been removed from the mixture.

The focus of the work presented here is the implementation of a classic "parent ion scanning" experiment using a MALDI quadrupole ToF geometry instrument. This approach has the potential to identify phosphopeptides from MALDI spots without enrichment, e.g. from LC separated digest samples. The technique employed makes use of enhanced duty cycle (EDC) technology to increase the detection efficiency of low mass diagnostic fragment ions. This approach has the potential to be expanded to cover other post translational modifications.

