

DEVELOPMENT OF A RAPID AND SENSITIVE METHOD FOR THE QUANTIFICATION OF BENZODIAZEPINES IN PLASMA AND LARVAE BY LC/MS/MS

Gert De Boeck¹, Nele Samyn¹, Karen Pien², Patrick Grootaert³ and Michelle Wood⁴.

¹ National Institute of Criminalistics and Criminology (NICC), Brussels, Belgium. ² Free University of Brussels, Belgium.

³ Royal Belgian Institute of Natural Sciences, Brussels, Belgium. ⁴ Waters Corporation, Manchester, UK.

INTRODUCTION

Benzodiazepines are the most widely prescribed psychoactive drugs in the world for the symptomatic treatment of anxiety and sleep disorders. However, misuse of these compounds has been reported and they are frequently encountered in postmortem blood analysis (suicide or accidental death).

Here we describe the development of a rapid and sensitive LC/MS/MS method for the quantification of 10 benzodiazepines. Limits of detection of 0.2 µg/L or better were achieved when just 25 µL plasma was used.

In addition, we present the application of this method to the analysis of benzodiazepines in *Calliphora vicina* larvae. Insects and their larvae are commonly used in the estimation of postmortem interval. Furthermore, they may serve as a reliable alternate source for toxicological analysis in the absence of suitable tissues and fluids that are normally taken for this purpose.

EXPERIMENTAL CONDITIONS

LC/MS/MS conditions

LC System: Waters Alliance® 2690
 Column: Conventional Phenyl Column
 (2.1 x 150 mm, 5 µm)
 Mobile phase : A =10:10:80 acetonitrile: methanol:
 20 mM ammonium acetate
 B = 95:5 acetonitrile:
 20 mM ammonium acetate
 Flow rate: 0.25 mL/min
 Injection volume: 10 µL

MS conditions:

Mass spectrometer: Quattro Ultima®

Ionisation Mode: ES positive ion

Capillary voltage: 3kV

MS/MS: MRM analysis (Table 1).

Collision gas Argon at 2.5×10^{-3} mbar

| Time (min) | A (%) | B (%) | Curve number |
|------------|-------|-------|--------------|
| 0 | 100 | 0 | 1 |
| 0.5 | 75 | 25 | 1 |
| 8 | 40 | 60 | 7 (concave) |
| 11 | 40 | 60 | 6 (linear) |
| 12 | 100 | 0 | 1 |
| 15 | 100 | 0 | 1 |

RESULTS AND DISCUSSION

Figure 1 shows the MS and MS/MS spectra for alprazolam. Table 1 summarizes the MRM transitions and conditions used for this and several other benzodiazepines (and their respective deuterated analogues). The latter were used as internal standards for quantification purposes.

A series of calibrators (1, 10, 40, 100, 200, 400 and 800 µg/L) were prepared by adding the benzodiazepines to drug-free plasma. Plasma samples were isolated from the matrix using a simple acetonitrile clean-up procedure (which also incorporates the addition of the internal standards).

Figure 2 shows the MRM chromatograms of the benzodiazepines obtained with a 10 µL injection of the 10 µg/L plasma calibrator. Quantification was performed by integration of the area under the specific MRM chromatograms. Figure 3 shows a typical standard curve for diazepam in plasma.

[APPLICATION NOTE]

Waters

THE SCIENCE OF WHAT'S POSSIBLE.TM

Responses were linear, in all cases, over the range investigated
(Coefficient of Determination > 0.99).

| Compound | Precursor ion (m/z) | Product ion (m/z) | Cone Voltage (V) | Collision energy (eV) |
|------------------|---------------------|-------------------|------------------|-----------------------|
| Alprazolam | 308.8 | 280.9 | 70 | 25 |
| Alprazolam-d5 | 313.8 | 285.8 | 100 | 25 |
| Clonazepam | 315.8 | 269.8 | 80 | 25 |
| Clonazepam-d4 | 319.9 | 273.8 | 100 | 25 |
| Diazepam | 284.9 | 154.0 | 60 | 25 |
| Diazepam-d5 | 289.8 | 153.7 | 80 | 25 |
| Flunitrazepam | 313.9 | 267.9 | 80 | 25 |
| Flunitrazepam-d7 | 320.8 | 274.8 | 80 | 25 |
| Lorazepam | 320.8 | 274.7 | 60 | 23 |
| Lorazepam-d4† | 326.8 | 280.8 | 80 | 23 |
| Nordiazepam | 270.9 | 139.8 | 80 | 25 |
| Nordiazepam-d5 | 275.9 | 139.8 | 80 | 25 |
| Oxazepam | 287.0 | 240.8 | 60 | 26 |
| Oxazepam-d5 | 291.7 | 245.8 | 80 | 26 |
| Prazepam | 324.9 | 270.9 | 80 | 25 |
| Prazepam-d5 | 330.0 | 276.0 | 80 | 25 |
| Temazepam | 300.9 | 255.0 | 60 | 25 |
| Temazepam-d5 | 305.8 | 259.8 | 60 | 25 |
| Triazolam | 342.9 | 307.7 | 60 | 25 |
| Triazolam-d4† | 349.0 | 313.9 | 60 | 25 |

Table 1. MRM transitions and conditions for the measurement of 10 benzodiazepines.

†Note that due to the isobaric nature between these benzodiazepines and their deuterated analogues alternative precursor ions were utilised.




Figure 1. MRM chromatograms for (top to bottom): lorazepam, temazepam, triazolam, prazepam, oxazepam, diazepam, alprazolam, flunitrazepam, nordiazepam and clonazepam. Responses were obtained with a 10 μL injection of the 10 $\mu\text{g}/\text{L}$ plasma calibrator.




Figure 2. MS and MS/MS spectra of alprazolam.

For all compounds, LOD's of 0.2 $\mu\text{g}/\text{L}$ (or better) and LOQ's of 1 $\mu\text{g}/\text{L}$ (or better) were achieved. The precision of the assay was assessed by performing replicate ($n=5$) extractions of plasma samples containing low, medium and high concentrations of the benzodiazepines (i.e. 2, 40 and 200 $\mu\text{g}/\text{L}$ respectively). Coefficients of variation (%CV's) were found to be highly satisfactory (<15%).

Compound name: Diazepam
 Coefficient of Determination: 0.999378
 Calibration curve: $1.22307^* x + 0.093412$
 Response type: Internal Std [Ref 10], Area* [IS Conc./IS Area]
 Curve type: Linear, Origin: Exclude, Weighing: 1/x, Axis trans: None




Figure 3. Typical response for plasma containing diazepam. Diazepam spiked plasma was firstly extracted using acetonitrile prior to analysis using LC/MRM. Benzodiazepines were quantified by reference to their deuterated internal standards.

The developed LC/MS/MS was subsequently applied to the analysis of *Calliphora vicina* larvae in a study to assess the feasibility of using insects and their larvae as alternate specimens in the absence of any suitable human specimens for toxicological analysis.

Larvae were reared on artificial foodstuff (beefheart) spiked with a range of concentrations of nordiazepam (0, 0.5, 1 and 2 µg/g). Post-feeding larvae were harvested (after 7 days) for analysis of drug content by LC/MS/MS. Figure 4 outlines the initial sample preparation method used for these specimens. All control larvae reared on spiked foodstuff were positive for nordiazepam and the metabolite oxazepam. All control samples were negative. Figure 5 shows the MRM chromatograms obtained following LC/MS/MS analysis of a control larva and a larva positive for nordiazepam. The method was sufficiently sensitive to measure benzodiazepines in single larvae whereas previous analytical techniques e.g. GC/MS, RIA, TLC have required pools i.e. typically 20 larvae.




Figure 4. Preparation of larvae for LC/MS/MS analysis.

CONCLUSION

We have developed a simple, rapid method that allows the simultaneous quantification of 10 benzodiazepines in plasma a single chromatographic run. LOD's were better than 0.2 µg/L when only 25 µL plasma was used. The method involves a simple protein precipitation step with acetonitrile followed by LC/MS/MS analysis.

The method was subsequently applied to the analysis of *Calliphora vicina* larvae in a study designed to assess the feasibility of using insects as alternate specimens in the absence of any suitable human tissues.

The sensitivity was such that it was possible to detect benzodiazepines in single larvae whereas previous methods have required pools.




Figure 5. MRM chromatograms obtained with the analysis of larvae that were reared on artificial foodstuff spiked with Nordiazepam at 0 and 1 µg/g (A and B respectively). Figure C shows the MRM chromatogram for the internal standard i.e. Nordiazepam-d5.

The MassTrak™ systems are in vitro diagnostic devices compliant with EU directive 98/79/EC. The systems are manufactured in accordance with US and European regulations. Waters Corporation is ISO 13485:2003 certified and operates in accordance with international regulations, United States Food and Drug Administration Quality System Regulations and current Good Manufacturing Practices.

Austria and European Export (Central South Eastern Europe, CIS and Middle East) 43 1 877 18 07, Australia 61 2 9933 1777, Belgium 32 2 726 1000, Brazil 55 11 5094 3788, Canada 1 800 252 4752 x2205, China 86 10 8586 8899, CIS/Russia 7 095 336 7000, Czech Republic 420 2 617 1 1384, Denmark 45 46 59 8080, Finland 358 9 5659 6288, France 33 1 30 48 72 00, Germany 49 6196 400600, Hong Kong 852 29 64 1800, Hungary 36 1 350 5086, India and India Subcontinent 91 80 2837 1900, Ireland 353 1 448 1500, Italy 39 02 27 421 1, Japan 81 3 3471 7191, Korea 82 2 820 2700, Mexico 52 55 5200 1860, The Netherlands 31 76 508 7200, Norway 47 6 384 60 50, Poland 48 22 833 4400, Puerto Rico 1 787 747 8445, Singapore 65 6273 1221, Spain 34 93 600 9300, Sweden 46 8 555 11 500, Switzerland 41 62 889 2030, Taiwan 886 2 2543 1898, United Kingdom 44 208 238 6100

All other countries: Waters Corporation U.S.A. 1 508 478 2000/1 800 252 4752



Waters
THE SCIENCE OF WHAT'S POSSIBLE.TM

Waters, XBridge and Alliance are trademarks of Waters Corporation.
All other trademarks are acknowledged.
©2007 Waters Corporation. Printed in the U.S.A.
March 2007 720001542EN MC-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com