

QUANTIFICATION OF MORPHINE, MORPHINE-3-GLUCURONIDE AND MORPHINE-6-GLUCURONIDE IN BIOLOGICAL SAMPLES BY LC/MS/MS

Michelle Wood¹ and Michael Morris¹.

¹Waters Corporation, Manchester, UK.

INTRODUCTION

Morphine is a potent analgesic isolated from the opium poppy papaver somniferum and traditionally used for the treatment of moderate to severe pain. Analgesia results from the action of morphine at the opioid receptors of the spinal cord and brain (Figure 1), where it attenuates both the speed of the impulse and the perception of pain.

In human subjects, morphine is extensively metabolised (primarily by conjugation with glucuronic acid) to form morphine-3-glucuronide (M3G) and morphine-6-glucuronide (M6G). Whilst, the principal metabolite i.e. M3G, has little or no analgesic effect, M6G has been shown to be highly effective and is believed likely to contribute significantly to the overall effectiveness of morphine¹. Hence, quantification of both the parent drug and metabolites is desirable for pharmacokinetic studies.

Previously we have described a LC/MS/MS method that allows the quantification of morphine and several other opiates in urine². Here we present a simple method that enables the quantification of morphine in plasma, whole blood and urine. Furthermore this procedure allows differentiation between two isobaric glucuronide metabolites.




Figure 1. Image of guinea-pig brain. The red areas represent the highest density of opioid receptors; yellow areas represent moderate density; whilst blue, purple and white represent low density.

METHODOLOGY

Sample preparation

Biological samples were prepared for LC/MS/MS analysis by means of a simple, solid-phase extraction (SPE) procedure. A Waters Oasis[®] HLB extraction Cartridge (1 cc/30 mg) was firstly conditioned with 1 mL volumes of each of the following: methanol, water and ammonium carbonate (10 mM, pH 8.8). Samples (100 µL, spiked with deuterated internal standards) were made up to a final volume of 1 mL with ammonium carbonate before applying to the pre-conditioned cartridge. The cartridge was then washed with 1 mL ammonium carbonate before elution of the sample using 100% methanol (0.5 mL). Eluents were dried using a Savant Speedvac Plus evaporator and then redissolved in 100 µL of mobile phase. Reconstituted samples were briefly vortex mixed before the analysis of 10 µL using LC in conjunction with multiple reaction monitoring (MRM).

LC/MS/MS

A Waters Quattro micro[™] triple quadrupole mass spectrometer fitted with ZSpray[™] ion interface was used for all analyses. Ionisation was achieved using electrospray in the positive ionisation mode (ES+). Details of the MRM conditions are given in Table 1.

Compound	Precursor ion (<i>m/z</i>)	Product ion (<i>m/z</i>)	Cone Voltage (V)	Collision energy (eV)
Morphine	286	165	45	38
Morphine-d3	289	165	45	40
Morphine-M3G-glucuronide	462	286	45	28
Morphine-M3G-d3-glucuronide	465	289	45	30
Morphine-M6G-glucuronide	462	286	45	28

Table 1: MRM transitions and conditions for the measurement of morphine and its metabolites. The deuterated analogues of morphine and morphine-3-glucuronide were also included for the purpose of internal standardisation.

APPLICATION NOTE

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

LC analyses were performed using a Waters 2795 separations module. Chromatography was achieved using a C₁₈ column (3.9 x 150 mm) eluted isocratically with 0.1% formic acid:acetonitrile (97:3) at a flow rate of 0.3 mL/min. Column temperature was maintained at 30 °C. All aspects of system operation and data acquisition were controlled using MassLynx™ 4.0 software with automated data processing using the QuanLynx™ program.

RESULTS

A series of calibrators (0.5-500 µg/L) were prepared in duplicate by adding standards to blank plasma, whole blood or urine. Samples were then extracted using the SPE method described above prior to LC/MRM analysis.

Following LC/MRM analysis, the areas under the specific MRM chromatograms were integrated.

Figure 2 shows the extracted MRM chromatogram of morphine, M3G and M6G obtained with a 10 µL injection of the 5 µg/L plasma calibrator. Opiates were quantified by reference to the integrated area of the deuterated internal standards. Responses were linear ($r = >0.999$) over the range investigated for all 3 compounds and in each matrix (Figure 3 shows a typical standard curve for M3G in urine).

SUMMARY

We present a sensitive method for the quantification of morphine and its glucuronide metabolites. The method involves a simple SPE purification prior to analysis using LC/MRM and is suitable for plasma, whole blood or urine samples.

REFERENCES

1. The Analgesic Effect of Morphine-6-Glucuronide. R Osborne, P Thomson, S Joel, D Trew, N Patel and M Slevin. Br J Clin. Pharmacol. 1992. 34 (2) 130-8.
2. Opiates: Use or Abuse? Quantification of Opiates in Human Urine. (Waters Application Brief WAB45). M Wood, K Rush*, M Morris and A Traynor*. Clinical Applications Development Group, UK Limited, Manchester, UK.

*Medscreen Ltd., 1A Harbour Quay, 100 Prestons Rd, London.




Figure 2. MRM chromatogram for morphine (MOR), M3G and M6G. The above responses were obtained with a 10 µL injection of the 5 µg/L plasma calibrator. Due to the isobaric nature of M3G and M6G chromatographic resolution is required to enable identification.

[APPLICATION NOTE]

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Compound name: Morphine-3-glucuronide
Correlation coefficient: $r = 0.999883$, $r^2 = 0.999766$
Calibration curve: $0.09404 * x + 0.100968$
Response type: Internal Std (Ref 4), Area * (IS Conc. / IS Area)
Curve type: Linear, Origin: Exclude, Weighting: 1/x Axis trans: None




Figure 3. Standard curve for M3G in urine. Responses (duplicates) were calculated in reference to the integrated area of the deuterated internal standards. The inserted figure shows the response for the range 0-10 µg/L.

The MassTrak™ systems are in vitro diagnostic devices compliant with EU directive 98/79/EC. The systems are manufactured in accordance with US and European regulations. Waters Corporation is ISO 13485:2003 certified and operates in accordance with international regulations, United States Food and Drug Administration Quality System Regulations and current Good Manufacturing Practices.

Austria and European Export (Central South Eastern Europe, CIS and Middle East) 43 1 877 18 07, Australia 61 2 9933 1777, Belgium 32 2 726 1000, Brazil 55 11 5094 3788, Canada 1 800 252 4752 x2205, China 86 10 8586 8899, CIS/Russia 7 095 336 7000, Czech Republic 420 2 617 1 1384, Denmark 45 46 59 8080, Finland 358 9 5659 6288, France 33 1 30 48 72 00, Germany 49 6196 400600, Hong Kong 852 29 64 1800, Hungary 36 1 350 5086, India and India Subcontinent 91 80 2837 1900, Ireland 353 1 448 1500, Italy 39 02 27 421 1, Japan 81 3 3471 7191, Korea 82 2 820 2700, Mexico 52 55 5200 1860, The Netherlands 31 76 508 7200, Norway 47 6 384 60 50, Poland 48 22 833 4400, Puerto Rico 1 787 747 8445, Singapore 65 6273 1221, Spain 34 93 600 9300, Sweden 46 8 555 11 500, Switzerland 41 62 889 2030, Taiwan 886 2 2543 1898, United Kingdom 44 208 238 6100

All other countries: Waters Corporation U.S.A. 1 508 478 2000/1 800 252 4752



Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Waters, XBridge and Alliance are trademarks of Waters Corporation.
All other trademarks are acknowledged.
©2007 Waters Corporation. Printed in the U.S.A.
March 2007 720001546EN MC-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com