

Introduction

The use of elevated temperature in liquid chromatography has recently been an active area of research (1-4). As the temperature of the mobile phase is increased, several variables that play a key role in chromatographic performance are affected. Mobile phase viscosity is decreased at higher temperatures. This means that with a constant mobile phase, flow rate, column dimension and particle size, operating pressure is reduced at higher temperature. Analyte diffusivity is increased at higher temperature. This fact results in a shift in the optimum flow rate to higher values. Horvath (5) has discussed the fact that under conditions where kinetics are fast, efficiency is independent of temperature, and this fact has been verified by Lee et al (3) and Sandra et al (4). Therefore, at least theoretically, we should not expect higher temperature to provide increased plate count or resolution for a given column length and particle size. Instead, due to the flattening of the van Deemter curve and the reduced mobile phase viscosity, higher temperature will allow faster separations to be obtained with equivalent resolution. However, the use of longer columns to provide increased plate count and resolution is enabled by higher temperature (4). While these longer columns can be used at ambient temperatures, the flow rates achievable under these conditions would lead to very long separation times. It is important to point out that for these longer columns the optimum flow rate is in most cases not achievable even at higher temperature, since with increasing temperature the optimum flow rate also increases.

Conclusion

- High temperature LC and UPLC allow faster separations, selectivity benefits, and sharper peaks.
- Efficiency is not improved at higher temperatures. It only allows longer columns which only then would you achieve greater N
- Increasing column temperature requires higher flow rates in order to run within the optimal conditions.
- Analyte degradation is possible and must be examined in every case
- Narrow peaks require fast scanning UV and MS detectors, but not necessarily above 40Hz.
- Alternative solvents can be used at higher temperatures for unique selectivity

References

1. X.Wang, W.Barber and P.Carr, *J.Chromatogr*, 1107 (2006) 139-151
2. Y.Xiang, Y.Lui and M.Lee, *J.Chromatogr*, 1104 (2006) 198-202
3. Y. Xiang, B. Yan, B. Yue, C. McNeff, P. Carr and M. Lee, *J.Chromatogr*, 983 (2003) 83-89
4. F. Lestremau, A. Cooper, R. Szucs, F. David and P. Sandra, *J.Chromatogr*, 1109 (2006) 191-196
5. F. D. Antia and C. Horvath, *J.Chromatogr*, 435 (1988) 1-15
6. Temperature and Pressure in Modern Chromatography Part 1: Theoretical Considerations Uwe D. Neue

Typical Chromatographic Phenomena Related to Using Elevated Temperatures : Predictable and Unpredictable

Elevated Temperature = Speed Benefits ≠ Increased Efficiency





Figure 5 - 6: Figure 5 demonstrates the isocratic separation of several alkyl phenones on a 2.1x150 mm ACQUITY UPLC BEH C18, 1.7 mm column. The mobile phase temperature was varied from 60 to 75 to 90°C, while the flow rate was varied to provide a roughly constant back pressure of approximately 7-8K and therefore constant efficiency (6). The resulting efficiency stayed constant at about 35,000 plates. At 60°C, the flow rate was 0.4 mL/min, resulting in a run time of about 4 minutes. At 75°C, the flow rate was 0.5 mL/min, resulting in a run time of about 2.6 minutes. At 90°C, the flow rate was 0.6 mL/min, resulting in a run time of about 1.8 minutes. This example demonstrates that for this sample mix, running at 90°C vs. 60°C results in a two-fold reduction in analysis time, with no loss in resolution. The data shown in the figure 6 shows that as the column temperature is increased, the gradient time can be reduced and the peak capacity retained, P_c 30°C = 283, 60°C = 284 and 90°C = 275.

Gradient



Effect of Temperature on Selectivity




Figure 7: As the column temperature is increased the rate of mass transfer and solubility of the analytes changes differentially. This can result in changes in the order of elution and the resolution between the peaks. This is of particular interest when transferring methodology from analytical to preparative scale. We can see that the temperature can be used as a tool to manipulate the peak separation profile such as changing the organic modifier to affect selectivity. Thus temperature can be viewed as a tool rather than a necessity.

Thermally Labile Compounds

The drug products are desired to have shelf life longevity with inexpensive storage conditions. Experiments show degradation of Lansoprazole on-column at temperatures excess of 90°C with significant breakdown above 100°C.




Figure 8: The combined extracted spectrum of the observed "bridged" phenomena tailing after the Lansoprazole analytic in the "blue" chromatogram is a result of on-column degradation yielding an m/z of 320 amu, a loss of ~48 amu.

Effect of Narrow Peaks on MS Data Quality

The combination of higher flow rates and longer columns gives rise to very narrow peaks approximately on average 1 second width peaks which places a significant strain on the detection system to collect data at a sufficiently fast rate to define the peak without affecting the sensitivity or accuracy of the mass spectrometry measurement. The spectra obtained from a rat urine analysis (figure 9) yield the endogenous metabolites xanthurenic acid (figure 10), hippuric acid (figure 11), kyurenic acid (2.3ppm), and pantothenic acid (1.4ppm) present in urine. We can see that analytes were collected with a mass accuracy of 3ppm or better.




Figure 9: Metabolite ID in Rat Urine. ACQUITY UPLC BEH C18 2.1 x 150mm; 1.7um. 90°C at 900uL/min. Peak capacity = 720. Peak widths ranged from just under 1 sec to 3sec.

Enabling Viscous Solvents with Higher Temperature

Solvents such as isopropyl alcohol and dimethyl-sulphoxide (DMSO) have been used rarely by chromatographers due to viscosity of these solvents. Raising the column temperature to 90°C significantly reduces the viscosity. The data shown below (figure 12) illustrates how the back pressure of a typical 2.1 x 100mm C18 ACQUITY UPLC column changes over the period of a 0-100% organic-aqueous gradient, for acetonitrile, methanol and isopropanol. We can see from this data that the isopropanol solvent reaches a back pressure maximum of 11,000psi. The benefits of using mobile phase modifiers such as isopropanol (IPA) are three fold; 1) faster analysis, 2) sharper peaks and 3) different analyte selectivity




Figure 12: Organic modifier pressure profile at 95°C




Figure 13: IPA offers speed. Compared to methanol and acetonitrile IPA gives a significant reduction in analysis time, this is because it has a greater elutropic strength.

Figure 14: The analysis of a more complex sample, in this case ginger root extract illustrates the benefits of using IPA as an organic modifier. We can see that these are sharper and better resolved in the IPA separation





Figure 15: The selectivity advantages of IPA. The analysis of the corticosteroid budesonide. The order of elution of the peaks has changed with IPA, the two large impurity peaks that eluted after the budesonide peak elute in the methanol and acetonitrile elute before the budesonide peak with the IPA gradient