

Polymer solutions were prepared using TCB solvent with 250mg/L San-INTRODUCTION tonox R antioxidant. The TCB was filtered using PTFE membrane (47mm diameter, 0.45µm pore size, Waters). Solutions were made by directly weighing samples and solvent into 60ml small neck glass bottles using an Gel Permeation Chromatography (GPC) of polyolefins uses 1,2,4electrobalance (Sartorius Analytic, Model A120S, Sartorius, Goettingen, trichlorobenzene (TCB) as the mobile phase at elevated temperatures Germany). Concentrations were calculated from densities and masses of (~140°C). This solvent is toxic and is on the Priority Chemicals List in the solvent and samples. Aluminum crimp caps with PETF/Aluminum/Silicone EPA National Waste Minimization Program. Replacing traditional GPC septa (Waters) were used to seal the bottles. A general-purpose laboratory columns with columns that enable lab personnel to run experiments in convection oven (VWR model 1330FM) set at 140°C was used to dissolve the less time and to generate less toxic waste (mobile phase) would benefit polymer solutions. 2 hours heating with gentle swirling (about 1 minute every any lab analyzing polyolefins. 30 minutes) was enough to dissolve polyethylene samples. After dissolution, polymer solutions were dispensed into 10 ml vials for GPC experiments.

The application of Waters® HSPgelTM series columns can reduce both run time and solvent consumption. With a single HSPgel (6.0x150 mm) column a sample run may finish within 7 minutes at a flow rate of 0.6 ml/min. The run time is ~5 X less and the solvent use per run is ~9 X less compared with a traditional analysis on a bank of three columns (7.8x300 mm). Using a cumulative matching calibration technique, the GPC results from the HSPgel columns are reproducible and accurate.¹

A key feature to evaluate when considering column change is that the resulting data are equivalent within acceptable limits. The molecular weight results obtained from a single HSPgel HT MB-M (5µm,6.0x150 mm) column were compared with the molecular weights from a traditional Styragel[®] HT (10 µm, 7.8x300 mm) three column bank. The GPC system used laser light scattering and viscometer detectors. The data shown here illustrate the usefulness of HSPgel HT series columns for polyolefin samples

EXPERIMENTAL

Waters Alliance [®] GPC/V2000 with
Precision Detectors PDI2040
690 nm
15 degree and 90 degree.
Alliance GPC2000 (version 2.0)
140°C (sample carousel, injector, and col-
compartments)
HSPgel HT MB-M (5µm, 6.0x150mm);
Styragel [®] HT column bank (2x HT6E and 1x
HT2,10 μm, 7.8x300mm)
1,2,4-trichlorobenzene (filtered, degassed)
0.6 ml/min for HSPgel HT column;
1.0 ml/min for Styragel HT columns
1Hz
52.5 µl for HSPgel HT MB-M column;
319.5 µl for Styragel HT column bank

<u>Materials</u>:

1,2,4-trichlorobenzene (TCB) is from Fisher Scientific (Pittsburgh, PA). Polystyrene narrow standards are from Toyo Soda (Japan). Broad polystyrene sample Dow1683 is from Dow Chemical. Polyethylene standards 1475a, 2886, and 2887 are from National Institute of Standard and Technology (NIST, Gaithersburg, MD). Santonox R antioxidant (CAS number 96-69-5) (Monsanto) was used.

Preparation of standards and samples

<u>Data processing</u>	
Software:	Empower™ (Build 1154)
Options:	GPC, GPCV, and LS
Light scattering calibran	its:
	PS narrow standards (Mw: 9,100 and 107,000)
dn/dc:	PS, 0.053 ml/g; PE, 0.109 ml/g
Refractive index:	1.520 (TCB)
Density:	1.31 g/ml (TCB at 140°C)
-	-



Figure 1. Chromatograms of a PS broad sample on an HSPgel HT MB-M column.

Sam

Dow16

NIST14

NIST28

NIST28

a: Results of conventional columns are the average of 2 injections for each sample except NIST 1475a, which is the average of 5 injections. b: The labeled intrinsic viscosity were measured at 130°C in TCB. c:140°C in TCB d: Results of HSPgel columns are the average of 3 injections for each sample.

2

results².

HIGH TEMPERATURE GPC ANALYSIS WITH HSPGELTM HT COLUMN

Table 1. Comparison of calculated molecular weight values with the vendor values for narrow polystyrene standards.

Vendor M _w	M _w	RD%
1,010,000	1,010,000	0
706,000	719,000	1.8
422,000	417,000	-1.2
186,000	198,000	6.5
107,000	107,000	0
65,900	64,800	-1.7



	Label Value			3 Conventional Columns ^a			1 HSPgel HT MB-M column ^d		
nple	Mn (1x10 ³)	Mw (1x10³)	[ŋ] ^ь (dL/g)	Mn (1x10³)	Mw (1x10³)	[ŋ] ^c (dL/g)	Mn (1x10³)	Mw (1x10 ³)	[ŋ]º (dL/g)
683	100	250	N/A	96	257	0.757	150	232	0.752
475a	18.3 ±0.4	53.1 ±0.6	1.010 ±0.009	17.6 ±1.6	54.0 ±1.0	1.023 ±0.008	18.4	51.2	0.979
886	N/A	87 ±6	1.58 ±0.02	77.6	81.1	1.54	77.6	80	1.517
887	N/A	196 ±14	2.77±0.03	166	190	2.71	162	179	2.765



MW

Figure 2. Mark-Houwink plots for Dow 1683 (PS) obtained from both HSPgel column and conventional Styragel HT columns and comparison with literature



Figure 3. Mark-Houwink plots for NIST 1475a (PE) from HSPgel column and conventional Styragel HT columns and comparison with literature values².

1) Jinchuan Yang, Rick Nielson, Peter Alden, Waters Chemistry Operation WW Meeting 2002, WA20727 2) Sun, T., Chance, R. R., Graessley, W. W., Lohse, D. J.; Macromolecules 2004, 37, 4304-4312



Jinchuan Yang and Alice J. Di Gioia Waters Corporation, 34 Maple Street, Milford, MA 01757

CONCLUSIONS

• HSPgel HT series columns are useful for high temperature GPC analyses with RI, viscometry and LS detectors.

• Molecular weights and intrinsic viscosities are within acceptable limits of the values obtained with a traditional 3-column bank

• Shorter run times result in less consumption of toxic solvents

• Run time here < 4.5 X compared with a traditional 3- column bank

• Solvent waste here < 7.5 X compared with a traditional 3- column bank

REFERENCE

	-			
	alliance	₽ GPCV 2000	ET.	
				*PD2040
Waters				Providence Determiner