

USE OF UPLC COUPLED WITH TOF MASS SPECTROMETRY FOR METABONOMIC INVESTIGATIONS OF PREDNISONE DOSING

Andrew Baker, Ignatius Kass, John Shockcor
Waters Corporation 34 Maple St Milford, MA

ABSTRACT

Fast UPLC separations coupled with high resolution TOF mass spectrometry are used to profile urine samples obtained after dosing with prednisone. Amino acid concentrations are monitored during therapy. PCA is used to identify components leading group separation.

INTRODUCTION

Glucocorticoid therapy is widely used for cancer chemotherapy, immunosuppression, and anti-inflammatory activity. Prednisone is one of the most common drugs used for the therapy, however there is a general activation of catabolic processes, including glycolysis and protein turnover. Several studies have attempted to describe the general mechanism for the protein turnover, finding an increase in urea production and excretion.¹⁻³

In addition to the classical *in vivo* drug metabolism experiments, metabonomic studies play an important role in assessing drug efficacy, drug safety, and mechanism of activity. Along with the quantification of the drug and its metabolites, profiling of other physiological fluid constituents can illuminate other metabolic processes that are affected by drug therapy.^{4,5} Metabolic profiling of physiological fluids such as plasma or urine is complicated by the number of endogenous components (sample complexity) and the number of analyses needed to gather meaningful statistical data.

Ultra performance liquid chromatography (UPLC) is characterized by the high chromatographic efficiency and high speed analysis accorded by the use of sub 2 μm particles.⁶ In addition to high efficiency and high speed chromatographic analysis, time of flight mass spectrometry is utilized to obtain high sensitivity full scan data that can be used for identification purposes. Accurate mass measurements generated from this data are used to deduce elemental compositions. In addition, the novel MS^E techniques are used to generate MS/MS like data for the whole dataset, reducing the necessity for dedicated MS/MS or DDA techniques.⁷

Because many hundreds to thousands of chromatographic components can be isolated during a short UPLC run, data reduction techniques must be used to detect, align, and filter out confounding variables so that biological components that are changing with respect to drug therapy can be identified. Principal component analysis and other statistical data treatments were used to characterize the time-course behavior of endogenous compounds found in urine during the drug therapy.⁸

METHODS

Prednisone Dosing and Sample Collection

Morning urine was collected from an adult male during treatment with Prednisone (60 mg/day for two days, with subsequent doses at 50 mg/day for two days, then 40 mg/day for two days). Aliquots were frozen at -40 °C until analysis. Samples were thawed, centrifuged at 14,000 g for 10 minutes, and then diluted with 10 volumes of water prior to analysis.

Liquid Chromatography/Mass Spectrometry

Waters ACQUITY UPLC

Mobile Phases

0.1% Aqueous Formic Acid
Acetonitrile

Column ACQUITY BEH C8 2.1x100 mm 1.7 μm d_p
Flow Rate 400 $\mu\text{l}/\text{min}$
Temperature 50 °C
5 μl Injections

Gradient:

Ramp from 5 % to 50% B in 5 minutes
Ramp to 95% B in 0.1 minutes
Hold 1 minute

Waters Q-ToF Premier TOF MS Analysis

ESI Positive or Negative Mode
Capillary= 3.5 kV
Cone Voltage= 25 V
Collision Energy = 4 eV
LockSpray using Leucine Enkephalin as a reference mass
W Mode resolution (fwhm) = 15,000 @ m/z 554

MS^E Experiments

ESI Positive or Negative Mode
Capillary= 3.5 kV
Cone Voltage= 25 V
Collision Cell 4 eV Low Energy
15 eV High Energy

LockSpray using Leucine Enkephalin as a reference mass
W Mode resolution (fwhm) = 15,000 @ m/z 554

Data Analysis

Quantitative data was processed using Quanlynx. Data was exported to Microsoft Excel for further analysis. Positive and negative ion LC/MS data sets were processed using Markerlynx to detect and quantify chromatographic components, time-align the chromatographic components, and normalize data from each sample. The data was processed using Principal Component Analysis to identify putative markers. The dataset was exported to SIMCA-P (Umetrics AB, Umea, Sweden) for further evaluation.

RESULTS




Figure 1. Markerlynx Browser Report for Negative Ion UPLC/TOF MS for time course study of prednisone dosing. Top Panel Sample Report; Middle Panel Table of identified components; Bottom Panes Chromatogram with identified components and extracted ion current plot for ion at 3.32 RT and m/z 481.2443, Marker Trend for ion at 3.32 RT and m/z 481.2443, PCA Scores Plot showing group separation, and PCA Loadings Plot of components leading to group separation. Scores Plot labels correspond to Day 0 (predose) through Day 6.

PLS-DA Scores and Loadings Plots




Figure 2. Partial Least Squares Discriminate Analysis plots for positive ion dataset for Day 0 vs. Days 1-6. Left Panel Scores Plot. Center Panel Loadings Plot. Right Plot Extracted Markers showing elevated response during the first three days of drug administration.

O-PLS Scores and Loadings S-Plots




Figure 3. Orthogonal Partial Least Squares plots for negative ion dataset for Day 0 vs. Days 1 and 2. Left Panel Scores Plot. Right Panel Loadings S-Plot. Markers with high or low pcorr values are responsible for the class separation between the pre-dose and post dose samples.

Identification of Markers




Figure 4. LC/MS Chromatograms for marker at 1.1 minutes and m/z 205.098. Top Trace XIC for Pre Dose Middle Trace XIC for Dose Day 2, Bottom Trace MSE Trace.

	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6
Carnitine	798.00	447.48	334.38	177.25	81.36	78.77
Tryptophan	112.58	148.85	116.56	137.62	121.28	112.01
Valine	86.21	n.d.	110.60	114.83	64.64	71.20
Leucine	n.d.	n.d.	47.11	n.d.	48.66	n.d.
Histidine	335.66	271.26	226.00	239.61	202.09	n.d.
Phenylalanine	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.
Glutamine	n.d.	119.93	74.75	113.18	105.60	89.68
Arginine	64.08	133.04	n.d.	164.11	62.82	54.34
Tyrosine	n.d.	99.42	n.d.	n.d.	80.06	n.d.
Glutamate	n.d.	n.d.	n.d.	n.d.	n.d.	n.d.
Methionine	128.34	233.54	56.29	n.d.	145.28	n.d.
Taurine	n.d.	55.35	n.d.	22.60	39.84	18.53
Threonine	109.86	143.71	81.61	124.51	141.03	113.56

Table 1. Statistically significant changes (97.5% confidence limits) in amounts of amino acid excretion. Data is expressed as percent change from the D0 level. n.d.= no statistically significant change between dosed level and pre-dose level




Figure 5. Background subtracted spectra for component at 1.1 minutes. Top Trace Spectrum for low energy acquisition (4 eV), Bottom Trace Spectrum from MSE high energy acquisition (18 eV).




Figure 6. Elemental composition report for low energy spectrum. Data-base searching of the elemental composition suggests tryptophan as a possible match.




Figure 7. Elemental composition report for MSE high energy spectrum at 1.1 minutes.

CONCLUSIONS

- UPLC has high resolution and peak capacity necessary for metabonomic profiling of complex biological separations.
- Use of supervised and unsupervised methods for determination of components responsible for group separation.
- Use of Q-ToF Premier MSE methodology gives comprehensive MS and MS/MS-like datasets to aid in small molecule biomarker characterization.
- Elevation of amino acids and carnitine in urine supports protein breakdown models for prednisone-induced muscle wasting.

References

- Lofberg et al, Effects of High Doses of Glucocorticoids on Free Amino Acids, Ribosomes, and Protein Turnover in Human Muscle. *Eur. J. Clin. Invest.* **2002**, *32*, 345-353.
- Short et al, Effect of Short Term Prednisone Use on Blood Flow, Muscle Protein Metabolism and Function. *J. Clin. Endocrinology and Metabolism* **2004**, *189*, 6198-6207.
- Wolthers et al, Effects of Budesonide and Prednisolone on Hepatic Kinetics for Urea Synthesis. *J. Hepatology* **2000**, *33*, 549-544.
- Nicholson et al, Metabonomics: a Platform for Studying Drug Toxicity and Gene Function. *Nature Reviews* **2000**, *1*, 153-161.
- Griffin, The Potential of Metabonomics in Drug Safety and Toxicology. *Drug Discovery Today* **2004**, *1*, 285-293.
- Mazzee et al, Advancing LC Performance with Smaller Particles and Higher Pressure. *Anal. Chem.* **2005**, *77*, 460A-467A.
- Silva et al, Absolute Quantification of Proteins by LC/MS^E: A Virtue of Parallel MS Acquisition. *Molecular and Cellular Proteomics* **2005**, *5*, 144-156.
- Holmes and Antti, Chemometric Contributions to the Evolution of Metabonomics: Mathematical Solutions to Characterizing and Interpreting Complex Biological NMR Spectra. *Analyst* **2002**, *127*, 1549-1557.