Waters

OVERVIEW

- Multiple reaction monitoring (MRM) by tandem quadrupole mass spectrometry is a powerful method for quantitatively screening protein biomarkers in biological samples. However, prediction of optimized MRM peptide transitions is challenging due to sample complexity and dynamic range.
- Empirical determinations of optimum transitions for MRM experiments are facilitated by alternate scanning (LC-MS^E) data acquisition which records a comprehensive inventory of peptide precursor and product ion accurate mass, retention time, and intensity information for a complex protein digest (Figure 1).
- The LC-MS^E discovery method was used to generate a targeted MRM assay for monitoring a panel of twenty serum proteins (Table 1).

INTRODUCTION

There is consensus that high and medium abundance serum proteins hold promise as potential clinical biomarkers. Broadband monitoring of these proteins is needed to improve clinical diagnostic capabilities across a diverse spectrum of human diseases. Mass spectrometry (MS) provides a means for such high throughput monitoring.

In this study, a comprehensive MS analysis of whole and depleted serum samples has been carried out via a label-free, alternate scanning LC-MS methodology (LC-MS^E).¹⁻⁴ The inventories of precursor and fragment ion information produced in this discovery method are a qualitative and quantitative record of the proteomic condition of the sample. The LC-MS^E data allows the empirical selection of optimum MRM transitions for any targeted peptide (Figures 2-3).

Optimum MRM transitions for peptides representing 20 serum proteins have been extracted for this work resulting in four monitored transitions per protein. In some cases, multiple charge states for a peptide were monitored. Limited cone voltage and collision energy profiling was conducted to crudely estimate preferred operating conditions for a given transition.

GOALS

- Demonstrate the utility of LC-MS^L for the development of optimized MRM quantitation of peptides associated with targeted proteins.
- Illustrate the sensitivity improvements associated with MRM for targeted peptide analysis.
- Compile retention time, peptide precursor and fragment ion information for the most sensitive and selective MRM transitions for serum proteins to be posted in the future at the Waters website.

REFERENCES

- 1. Bateman, RH. et al., *J Am Soc Mass Spectrom* **2002**, 13, 792-803
- 2. Purvine, S. et al., Proteomics 2003, 3, 847-997
- 3. J. Silva et al., Anal. Chem. 2005, 77, 2187-2200
- 4. J. Silva et al., *Mol. Cell. Proteomics* **2006**, 5, 144-156
- 5. J. Silva et al., Mol. Cell. Proteomics 2006, 5, 589-607
- 6. M. Hughes et al., *J. Proteome Res.* **2006**, 5, 54-63

METHODS

Depletion

Depletion of serum was carried out using a ProteoPrep® 20 Plasma Immunodepletion Kit (Sigma Product No. PROT20, Sigma-Aldrich, St. Louis. MO)

Samples

Whole and depleted serum protein samples were taken up in 50 mM NH₄HCO₃, 0.1% RapiGest[™] SF (Waters Corp., Milford, MA), pH 8.5 to a final concentration of ~ 1 μ g/ μ L. Reduction and alkylation was carried out with DTT and IAA, respectively. The proteins were digested with 1:25 (w/ w) proteomics grade, dimethylated trypsin (Sigma Product No. T6567, Sigma-Aldrich, St. Louis, MO) overnight (16 hr). Tryptic activity was degraded by the addition of 2 µL concentrated HCI, followed by centrifugation and the supernatant collected. Samples were spiked with known quantities of rabbit phosphorylase B and yeast enolase digest and diluted with 50 mM NH₄HCO₃ to a final concentration of 200 ng/µL of serum peptides for 600 ng loads.

LCMS^E Data Collection

The peptides were separated and analyzed by LC-MS^L using a nanoACQUITY[™]-Q-Tof Premier[™] (Waters Corp., Milford, MA), and then processed to produce a COMPREHENSIVE INVENTORY of identified peptides with their corresponding precursor and product intensities, their accurate mass and retention time information (AMRT), that was then used to select optimum MRM transitions.

MRM Data Collection

MRM analyses were performed on a Quattro Premier[™] XE tandem quadrupole mass spectrometer (Waters Corp., Milford, MA), monitoring four transitions per protein. They were selected from the AMRT inventories of each protein based upon product ion intensity and replication in at least two out of three experiments.



Figure 1. Alternate Scanning LC-MS (LC-MS^E). LC-MS^E data acquisition produces an inventory of all precursors and fragment ions detected throughout the entire experiment, with their corresponding accurate mass, retention time, and intensity.

TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS

INTELLIGENT DESIGN OF MRM ASSAYS FROM LC-MS^E INVENTORIES

Michael J. Nold, Gordon Fujimoto, Asish B. Chakraborty, Scott J. Berger, Scott J. Geromanos, Jeffrey C. Silva Waters Corporation, Milford, MA





Product Inventory







Figure 2. Extraction of MRM Transitions from AMRT Inventories. A. Total AMRTs from 3 replicates of whole and depleted serum tryptic digests B. Filtered AMRT precursor inventory - Found in at least 2 out of 3 experiments C. Vitamin D binding protein product inventory

RESULTS & DISCUSSION

Number Protein.Name

- beta 2 glycoprotein I precursor complement component 3 pred complement component 4 bind complement component 4B pre complement component 5
- complement component 6 pred complement component 7 prec
- hemopexin
- heparin cofactor II
- histidine rich glycoprotein pre-
- I factor complement
- inter alpha globulin inhibitor H
- inter alpha globulin inhibitor H
- inter alpha globulin inhibitor H **PREDICTED** similar to Apolipo
- serine or cysteine proteinase in
- serine or cysteine proteinase in
- serine or cysteine proteinase in
- vitamin D binding protein precu vitronectin precursor

number represents the fragment m/z.

Specificity & Signal to Noise Benefit of MRM

	LC-MS ^E , Discovery		MRM, Targeted		LC-MS ^E , Discovery			MRM, Targeted			
mber 12	Injection 1 Injection 2 Injection 3	$\begin{array}{cccccccccccccccccccccccccccccccccccc$	T-1 T-2 T-3 T-4 $\begin{bmatrix} 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 $	Sample Number Whole 12 13 14 Depleted 12	P-1 P-2	P-3 • •	P-4	T-1	T-2	T-3 • •	T-4
	S/N _{AVERAGE}	119 116 116 14	482 663 208 32	13 14	• •	•	•	•	•	•	•
13	Injection 1 Injection 2 Injection 3			Figure 4. Detection Array for inter alpha globulin inhibitors. Numbers 12-14 refer to the inter alpha globulins listed in Table 1. LC-MS ^E precurse ion detections or MRM transition detections made in at least 2 out of 3 of the respective experiments are colored green, otherwise they are colored red. The sensitivity an specificity of the MRM method allowed detection of all 12 transitions from the deplete serum sample that were not detected by LC-MS ^E .						rsor tive and eted	
	S/N_{AVERAGE}	208 55 118 122	873 343 452 156		COI	NCLUS	IONS	5			
14	Injection 1 Injection 2 Injection 3			 LC-MS^E data is a comprehensive, empirical guide to the best MRM transitions for quantitative monitoring of targeted proteins. MRM provides additional sensitivity to take advantage of the specificity revealed by LC-MS^E. 							
	S/N _{AVERAGE}	54 279 12 154	236 1040 65 265	 In controlled experiments, a 50x or greater improvement in detection limit was demonstrated in MRM acquisition mode relative 							

Figure 3. Signal to Noise for inter alpha globulin inhibitors from Whole Serum Numbers 12-14 refer to the inter alpha globulins listed in Table 1. P-1 through P-4 are selected ion chromatograms, with 50 mDa tolerance, of the precursors from the low energy LC-MS^E data that were selected for each inter alpha globulin inhibitor, respectively. T-1 through T-4 represent the corresponding MRM traces from the targeted assay. These data were collected with dedicated collision energies, and illustrate the sensitivity benefit (S/N) of MRM even before collision energy optimization.

	T-1	T-2	Т-3	T-4
	795.0 → 1842.8	638.7 → 927.4	638.7 → 987.6	511.8 → 652.3
ursor	595.8 → 567.4	895.4 → 602.3	595.8 → 666.4	595.8 → 624.3
ling protein	569.6 → 1550.8	564.8 → 591.3	569.6 → 1083.6	625.3 → 644.4
eproprotein	828.4 → 869.5	828.4 → 956.5	771.4 → 1129.6	771.4 → 574.3
	626.9 → 925.5	728.3 → 843.4	629.3 → 1739.8	626.9 → 1038.6
ursor	623.8 → 747.4	1005.5 → 811.5	621.0 → 538.3	1000.8 → 701.4
ursor	409.2 → 799.4	550.3 → 887.4	605.3 → 677.4	753.0 → 581.3
	610.8 → 959.5	571.7 → 586.3	833.4 → 663.3	610.8 → 775.4
	514.8 → 814.4	465.7 → 712.4	540.3 → 893.5	410.3 → 706.4
cursor	912.9 → 300.2	841.9 → 1058.4	841.9 → 1171.5	$\textbf{1008.5} \rightarrow \textbf{1106.6}$
	$\textbf{728.9} \rightarrow \textbf{300.2}$	579.8 → 924.4	728.9 → 856.4	719.9 → 841.5
1	579.3 → 902.5	669.3 → 775.4	669.3 → 874.4	437.3 → 631.4
2	791.9 → 1341.7	669.4 → 686.4	514.3 → 797.5	423.7 → 570.3
4	467.3 → 720.4	619.3 → 894.5	464.8 → 702.3	572.8 → 663.4
protein A I precursor	626.8 → 1025.5	626.8 → 1025.5	516.3 → 831.4	693.9 → 940.5
nhibitor clade A	888.5 → 718.4	508.3 → 829.5	888.5 → 605.3	$\textbf{624.3} \rightarrow \textbf{1134.5}$
nhibitor clade C	695.4 → 1161.6	625.6 → 687.3	454.7 → 666.3	454.7 → 795.4
nhibitor clade F	528.3 → 855.5	652.7 → 1771.9	489.3 → 656.4	$\textbf{780.4} \rightarrow \textbf{1134.5}$
ursor	903.1 → 1351.6	400.3 → 700.4	400.3 → 587.3	1259.5 → 926.4
	657.8 → 686.3	$\textbf{657.8} \rightarrow \textbf{629.3}$	556.3 → 691.3	$\textbf{835.4} \rightarrow \textbf{591.3}$

Table 1. Summary of 20 representative serum proteins and transitions derived from an LC-MS^E inventory of whole and depleted serum digests. This is a subset of LC-MS^E identified proteins from whole and depleted human serum which were then monitored by MRM. The MRM experiments were configured to monitor at least 4 transitions per protein (designated T-1 through T-4), with many peptides having more than 1 transition monitored. The first number of each transition represents the precursor m/z, and the second

Sensitivity of MRM

- to LC-MS^E acquisition mode.

FUTURE WORK

- Increase the panel of proteins within a single experiment.
- Optimize the LC-MS^E inventories through detailed evaluation of precursor and fragment ion energetics. 720001923EN