Developing Strategies for the Isolation and Purification of Synthetic Peptides

Monika Lamouri¹, Jo-Ann Jablonski², Ziling Lu², Joomi Ahn², Fang Xia², Thomas E. Wheat², and Jeffrey R. Mazzeo²

Waters Corporation¹ Rue Jacques Monod; 78280 Guyancourt; FRANCE; ² 34 Maple Street; Milford,

INTRODUCTION

The purification of synthetic peptides is important for many kinds of research and development. The required peptides present a problem in isolation because they can have a wide range of chemical and physical properties. Such differences often lead to specialized isolation methods that are adjusted for each sample. We have investigated the variables that can be manipulated to control a peptide separation so that an efficient isolation protocol can be developed. This includes selection of a column based on the properties of the peptides. The effects of column chemistry, pore size and particle size have been considered. Modifying method parameters, such as gradient slope, will also be described. These considerations can be used to suggest a general protocol for high throughput purification. Some samples still require optimized separations, either because of the requirements of subsequent experiments or because of the complexity of the synthetic mixture. These separations are optimized on smaller, pilot-scale columns. We will demonstrate an efficient process for scaling the optimized separation to a larger sample load. The combination of an optimized general protocol and systematic scaling procedures can improve efficiency in the isolation of synthetic peptides.

RESULTS FIGURE 1: COLUMN COMPARISON

FIGURE 1A: BEH130; 3.5µm; TFA



FIGURE 3: PORE SIZE EFFECTS FIGURE 3A: LARGE PEPTIDE



MATERIALS AND METHODS LC Conditions

LC System: Waters 2796 Separation Module Waters 2487 Dual Wavelength UV Detection: Absorbance Detector Wavelength 214 nm Waters ZQ[™] Mass Spectrometer MS System: Electrospray Ionization (+) Mobile Phase: With "TFA" modifier: A = 0.1% Trifluoroacetic Acid in Water B = 0.1% Trifluoroacetic Acid in Acetonitrile With "FA" modifier A = 0.1% Formic Acid in Water B = 0.1% Formic Acid in Acetonitrile

Gradient for Column Evaluation

Time	%A	%B	Flow	Curve
0	95	5	1.5	*
50	50	50	1.5	6
52	10	90	1.5	6
54	10	90	1.5	6
55	95	5	1.5	6
72	95	5	1.5	6

Columns for Evaluation

Peptide Separation Technology

250 5.00 7.50 10.00 12.50 15.00 17.50 20.00 22.50 25.00 27.50 30.00 32.50 35.00 37.50 40.00

FIGURE 3B: HYDROPHOBIC PEPTIDE



Figure 3: Effect of Pore Size. Large pore size packings are usually selected for peptide separations to minimize steric effects. With the largest peptide in the test set, narrow, symmetrical peaks are observed on both the small and large pore materials. It is also commonly reported that hydrophobic peptides behave better with larger pores, but no significant effect is observed with these peptides. While the 130Å is more retentive than the 300Å, no clear size or hydrophobicity limit is found in this test set. Steric effects certainly exist but probably only become important for significantly larger peptides.

FIGURE 4: Effect of Modifier



Figure 4: Effect of Modifier. The basic test peptide was separated on the XBridge[™]130 packing material in the presence of 0.1%TFA or 0.1%formic acid. As expected, the elimination of ion pairing reduces retention for all the peptides in the sample. Some peptides are affected more than others resulting in selectivity changes that should be useful in developing purification protocols.

XBridgeTM BEH130 C₁₈ 4.6mm x 150mm 3.5 μ m XBridgeTM BEH130 C₁₈ 4.6mm x 150mm 5 μ m XBridgeTM BEH130 C₁₈ 4.6mm x 150mm 5 μ m Prep XBridgeTM BEH130 C₁₈ 4.6mm x 150mm 10 μ m XBridgeTM BEH300 C₁₈ 4.6mm x 150mm 3.5 μ m XBridgeTM BEH300 C₁₈ 4.6mm x 150mm 5 μ m XBridgeTM BEH300 C₁₈ 4.6mm x 150mm 5 μ m

Test Peptides

Neutral and "Average ISQAVHAAHAEINEAGR 17 residues; 1772.89Da pl 6.4; B&B 5580; HPLC 26.5 Basic AIKEVKQTIVKHPRY 15 residues; 1809.06Da pl 10.3; B&B -860; HPLC 21.7 Hydrophobic WLTGPQLADLYHSLMK 16 residues; 1871.96Da pl 7.3; B&B -5200; HPLC 113.3

Large

Released Biopharmaceutical Dosage Form Approximately 40 residues; Over 4000 Da pl 4.7; B&B 2540; HPLC 89.6



Figure 1: Column Comparison. The four test peptides representing a wide range of chemical and physical properties gave symmetrical peaks on the large and small pore size columns in both mobile phase modifiers. Selectivity differences are observed both between the pore sizes and, as expected, with the different modifiers. This small set of columns can, therefore, provide the basis for developing methods for isolation.

FIGURE 2: Effect of Particle Size



Figure 2: Effect of Particle Size. The basic test peptide was separated on the XBridgeTM130 packing material in the presence of 0.1%TFA. The four different particle sizes were packed into identical columns and separated with the same method. The selectivity is identical in all four experiments so methods can be scaled. Resolution is nearly the same on the $3.5\mu m$ and the two $5\mu m$ particles. Larger resolution losses are observed with the $10\mu m$.

CONCLUSIONS

- The BEH Technology[™] packing material used as the basis for Peptide Separation Technology columns gives good peak shape for peptides over a wide range of properties.
- When TFA is replaced with formic acid as modifier, retention is reduced, as expected, but symmetrical peak shape is retained.
- Retention and selectivity are identical across a range of particle sizes, as required for scalability.
- Peaks are broader with larger particle size packings.
- There is no clear correlation between pore size and peptide molecular weight for this set of test peptides covering up to about 4kD.
- There is no clear correlation between pore size and peptide hydrophobicity for this set of test peptides covering up to about 40% acetonitrile.
- Changes in mobile phase modifier affect retention and selectivity.
- Secondary interactions between peptides and the stationary phase are minimized with Peptide Separation Technology columns.

©2006 Waters Corporation

TO DOWNLOAD A COPY OF THIS POSTER VISIT WWW.WATERS.COM/POSTERS