

TOP DOWN SEQUENCE ANALYSIS OF INTACT PROTEIN SPECIES USING ION MOBILITY COUPLED WITH TIME OF FLIGHT MASS SPECTROMETRY

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

Therese McKenna, James I Langridge and Christopher J Hughes
Waters Corporation MS Technologies Centre, Wythenshawe, Manchester, UK

INTRODUCTION

In bottom-up proteomics, intact proteins are subjected to a proteolytic enzyme, for example trypsin, to produce a mixture of peptides. These species are then usually separated by nanoscale reversed phase chromatography and introduced into a mass spectrometer where fragment ion information is generated using collision induced dissociation. Each peptide intact mass and associated fragment ion spectrum is then compared to the predicted fragment ions of tryptic peptides from proteins in a sequence database. Protein identification is accomplished by locating the protein in the database that best matches the observed tryptic peptides and / or their fragment ions. A potential drawback with bottom up proteomics is that information about the intact mass of the protein is lost during the digestion procedure and hence the determination of some Post Translational Modifications (PTMs) may be difficult if the peptide on which a PTM resides does not ionise well.

Top down proteomics is the identification and characterization of intact proteins from tandem mass spectrometry experiments and can enable the identification of PTMs that may not be identified in the bottom up approach. This is because the intact protein mass can be measured directly prior to any fragment ion information being generated. However, even with a small to medium sized protein the range of charge states produced during the fragmentation of a single charge state leads to a spectrum that can be extremely complex and challenging to interpret.

In this presentation, we have used ion mobility separation prior to Time of Flight mass analysis to enable the separation of different charged state fragments and by post processing of the data, produced a simplified top down fragment spectra containing ions of mainly one charge state. These ion mobility separations and region selections also reveal low abundance species that are masked by larger species in a spectrum where no mobility separation is used.

METHODS

Instrumentation

The Mass Spectrometer used in these studies was a Synapt HDMS (Waters Corporation), **Figure 1**. Briefly, ions produced by an ESI probe are sampled by a Z-spray source. They pass through a quadrupole which may be set to transmit a substantial mass range or to select a particular m/z. The Tri-Wave comprises three T-Wave devices [1]. The first device (trap T-Wave) accumulates ions and releases them in a short pulse

(100 μ s) every 20 ms into the next device (IMS T-Wave) in which the mobility separation is performed, the final device (transfer T-Wave) is used to transport the separated ions into the oa-ToF for subsequent analysis. Ions may be fragmented on entrance to the accumulation T-Wave and/or in the transfer T-Wave. The pressure in the accumulation and transfer T-Wave regions was $\sim 10^{-2}$ mbar of Argon and the pressure in the IMS T-Wave was 0.5 mbar of N₂. The T-Wave pulse velocity and voltage were optimised to provide ion mobility separation.




Figure 1 Schematic diagram of the Waters Synapt HDMS instrument.

Sample

Bovine Beta Lactoglobulin contains two disulphide bonds and was analysed in both the non-reduced and reduced (DTT treated) form. The samples were prepared in a solvent composition of 50/50 acetonitrile / water + 0.1% formic acid and the solution was introduced to the mass spectrometer at a flow rate of 400nL/min.

Experiment

Intact protein ToF and Mobility mass spectra were acquired in order to determine the mass of the protein and the m/z of the different charge states. The quadrupole was then set to transmit a single charge state and the trap collision energy raised to produce a fragment ion spectrum, **Figure 2**. Different charge states were subsequently isolated and further fragment ion spectra generated. From m/z vs drift time plots, regions were selected and exported to produce new datafiles of simplified spectra of mainly one charge state. For z>1, these spectra were deconvoluted using Maximum Entropy algorithms. Fragment ion stretches were used to identify the amino acid sequence of the protein.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS




Figure 2 Synapt Trap T-Wave Fragmentation. Fragment ions produced before the Ion Mobility device are separated based on charge state, size and shape. They are subsequently passed to the oa-ToF analyser for mass analysis. Typical collision energies used in this study were 35-55V.

RESULTS

