AN INTEGRATED INFORMATICS SOFTWARE APPROACH DURING METHOD DEVELOPMENT USING EMPOWER 2, NUGENESIS SDMS, AND SDMS VISION PUBLISHER

Chris L. Stumpf, Michael D. Jones, Carla Sehlmeyer, and Phil Kilby Waters Corporation, Milford, MA, U.S.A.

INTRODUCTION

Traditional chromatography method development evaluates numerous chromatographic variables that may include selection of the column, solvent, flow rate, etc.¹ Conducting a method development project often takes weeks-to-months and generates numerous chromatogram reports demonstrating the evaluation of each chromatographic condition. Organizing and maintaining the chromatographic reports for archiving purposes as well as documenting the results and writing the final method report presents an overwhelming administrative challenge to the scientist and consumes 20-40% of valuable research time that could have been applied to new projects.²

A Scientific Data Management System (SDMS) provides the ability to automatically organize and store scientific data from diverse chromatography systems as well as organize research documents. The benefits of this approach are three-fold: (1) it frees the chromatographer to focus on science while the informatics system takes care of the mundane data management, (2) electronic search tools make locating data easy, and (3) the data and results become accessible to multidisciplinary teams due to the simple user interface and global accessibility. This application note demonstrates an integrated approach to capturing and storing separation related reports during chromatography method development by using Empower™ 2, NuGenesis® SDMS, and SDMS Vision Publisher™. NuGenesis SDMS acts as a data repository for analytical raw data as well as printed reports. Whereas SDMS Vision Publisher provides a reporting and authoring module incorporated as standard into NuGenesis SDMS. In addition to capturing the Empower 2 separation reports within NuGenesis SDMS, the relevant findings are composed into a summary report by using SDMS Vision Publisher. Figure 1 illustrates the workflow followed in this discussion and the chromatography method was developed for simvastatin.^{3.4}

THE SCIENCE

WHAT'S POSSIBLE."



Figure 1. Method development workflow followed while employing NuGenesis SDMS to manage the chromatography reports and results.

[APPLICATION NOTE]

DISCUSSION

Searching

Performing a literature search provides a logical starting point when beginning a method development project. Searching an organization's private knowledge repository for all prior research on a focus compound may return information on previously developed analytical assays, compound stability information, etc. Figure 2 illustrates an electronic search for prior research conducted on simvastatin (a cholesterol lowering compound) by using the search tool provided by SDMS Vision Publisher.

• Maters eLob Notebook - [Expe	riments]			
File Experiment Chemistry Service	zes Administration Options Window Help			- 9.×
Messages Experiments	New Experiment			ч
Take Over	B File: N/A ShowDescription @ties @Deleta	Cearfiter 🔍 Filter 🔹		
8 🖱 🎽 🔥 - 🕊	Espeinent	- 04		
Laboratory / Project / Template Group	D Tile Status	Ceded ton To		
P40	Seventain	· • • • •	7/ 2/2007 -	
	Laboratory Project	Approved from To	_	
		7/ 2/2027 •	7/ 2/2007	
	Couch [Advant] Metala			
	Sagle [Advanced] Mrtadas G. Frid @Char @ OpenExpriment(s) + pl Ault Iral + @Hite/Ris Copenantia found 15	и И		
	Sagle (Adveced) Meades C, Find (Course) (Count-operated) - pl Ault Inst - pl Ault Inst Experiment Sound 15 Experiment ID Te	e Subtin	Desciption	Ceator Name
	Segle [Advance] Mittadds] [1] Ford [2] Char (# Conn-Lowinetch) - gl Audit Ind - grinds File (Expension) Touristics (S) [Connect Touristics (S)	e Sublike	Desciption	Cesto Name Mr. Michael Jones
	Stelle [Advecte] Metadas (-) End (Care: (E) Care Contraction - pi Auto Trat) gines inter Economication (E) (E) (Careline C) (E) (E) (Care) (Care) (E) (E) (E) (E) (E) (E) (E) (E) (E) (E	5.80% 2007_PHD_00033	Desciption	Ceator Name Mr. Michael Jones Mr. Paul Rainville
	State Antroite Nr. M. (2000) (2000) Vertical (2000) (2000) Control (2000) (1000) Control (2000) (1000) El 2001 (Prior (2000)) Simetafic Addres	w Subite 2007 (Perty jacobis Neg 2007	Description	Cesto Name Mr. Nichael Jones Mr. Paul Rainille Mr. Nichael Jones
	Sigit (Advect) (Mode)	т 5.666 2007./нос.,000035 Мар. 2007 1960	Desciption APPLICATER - Foced depadation of	Cesto Nane M: Michael Jones M: Paul Rain-Be M: Michael Jones M: Peter Alden
	State (Antones) Mendes) State (Conservation of Conservation o	* 5486 200 /Hor _00035 May 200 PAG_381 PAG_381	Decolption APPLCATION - Focuel departation of APPLCATION - Focuel departation of	Cesto Nune M. Mchael Jones M. Paul Raindle M. Mchael Jones M. Pater Alden M. Pater Alden
	Sigit (Advect) (Medic) -V.M. (Kow (@ conceptioned)) #/Auth Sol - #/Auth	и 54466 2007 /ног, 90028 Мар. 2007 Раб., 381 Раб., 381 Раб., 381	Description APPLCATON - Found departement APPLCATON - Found departement APPLCATON - Found departement	Cesto Name Mr. Michael Jones Mr. Paul Raindle Mr. Paul Raindle Mr. Pater Alden Mr. Pater Alden Mr. Pater Alden
	State Februerell Mendels Verder III: Concer III: Concer III: Concer Concentration III: Concer III: Concer III: Concer Concentration III: Concer III: Concer III: Concer III: Concer Concentration III: Concer III: Concer <td>* 5466 207_96_0033 He_307 Pe_38 Pe_38 Pe_38</td> <td>Deciption APECATON - Food department APECATON -</td> <td>Casto Nane M: Michael Jones M: Paul Rainvile M: Nichael Jones M: Peter Alden M: Peter Alden M: Peter Alden</td>	* 5466 207_96_0033 He_307 Pe_38 Pe_38 Pe_38	Deciption APECATON - Food department APECATON -	Casto Nane M: Michael Jones M: Paul Rainvile M: Nichael Jones M: Peter Alden M: Peter Alden M: Peter Alden
	Steph Februariti Mindels Number Control C	5.446 2007/Hz (20028 Mig 2007 Pic (2017 Pic (2017 Pic (2017 Pic (2017 Pic (2017) Pic (20	Description APELCATEN - Found departition of APELCATEN - Found departition of APELCATEN - Found departition of APELCATEN - Found departition of APELCATEN - Found departition of	Cento Nane M: Michael Jones M: Paul Ruinde M: Nichael Jones M: Peter Alden M: Peter Alden M: Peter Alden M: Peter Alden

Figure 2. A search performed on the keyword "simvastatin" within Vision Publisher returns previously performed intra-company research reports.

The search results produced several documents relating to previous separation studies and indicated that simvastatin degrades at pH > 6.5. Therefore, exploring pH conditions above 6.5 was an unreasonable approach. Hence, searching a globally accessible repository provided a good starting point for the proposed study.

Report Management Scenario

Chromatographers evaluate numerous chromatographic condition variables during a typical method development study and Chromatography Data Systems (CDS) manage the resulting data by either utilizing a relational database (e.g., Empower) or by storing files on a file server. However, the data and results are generally only accessible to separation specialists trained to use that particular CDS. In some cases chromatographers need to share unpolished results before writing the formal method report. For example, a NMR spectrometrist may need access to raw chromatograms while helping elucidate the structure of co-eluting chromatographic peaks. In such cases, publishing raw separation reports into a CDS independent data repository aids in sharing results and promotes collaboration.

Figure 3 illustrates raw chromatogram reports flowing from a CDS and into a globally accessible data repository. In this case, printing Empower 2 chromatogram reports into NuGenesis SDMS after evaluating each chromatographic variable provides an easy mechanism to organize the study results. NuGenesis SDMS automatically organizes printed Empower 2 reports based on scientifically relevant text found within each report, e.g., the column type, buffer, pH, etc. Hence, locating and highlighting important findings from the chromatogram reports while writing the final method report in SDMS Vision Publisher becomes simplistic due to the streamlined data storage.

[APPLICATION NOTE]

Empower Reports Demonstrating each chromatographic variable evaluated during method development.



Figure 3. Empower reports are printed into NuGenesis SDMS. Then a specific subset of results is highlighted within the SDMS Vision Publisher final method report.

Creating the final method report

Figure 4 shows the final method report including chromatograms extracted from the raw Empower 2 reports stored within NuGenesis SDMS. The structured organization of the reports within NuGenesis SDMS made finding relevant information to highlight in the final report straightforward (Figure 3).

The reports were organized into columns and rows based on scientifically relevant text (e.g., project, buffer, solvent, column selection, analyst, etc.). Sorting these descriptive columns provided an easy mechanism to evaluate one condition at a time. A particularly useful example involved a qualitative mathematical function (injection score) evaluated within Empower 2 and used to classify the quality of each chromatogram separation.

The injection score column ranked the separations from bestto-worst. See the NuGenesis SDMS Data Repository in Figure 3. Hence, evaluating the overall quality of the separations became streamlined and locating representative separations for highlighting in the final report was simplified.

The final report consisted of an automatically generated index, with direct links back to the complete Empower 2 reports that have been automatically captured and cataloged within NuGenesis SDMS, an electronic audit trail that tracked edits, and an electronic signoff. Hence, documentation and data management efficiency were enhanced by about 20%. Considering published studies indicate that scientists spend 20-40% of their research time on average maintaining research data and documentation,^{2,5} the outlined process in this application note highlights practical tools and procedures that have a direct impact on workflow and productivity.



Figure 4. The SDMS Vision Publisher document describing the simvastatin method development. The chromatograms from representative separations are shown along with the links back to the original completed chromatographic reports.

CONCLUSIONS

Conducting an intra-company electronic search on prior simvastatin research by using the SDMS Vision Publisher embedded tools provided a reasonable starting point for the chromatography method development study thus saving valuable time and effort.

NuGenesis SDMS automatically stores, organizes, and ranks separation quality of the Empower 2 reports upon evaluating each chromatographic condition. This reduced the data management burden on the chromatographer and improved productivity by about 20%.

By publishing the reports into NuGenesis SDMS, they became globally available within the organization thus promoting multidisciplinary collaboration.

Composing a final method report within SDMS Vision Publisher highlighted the optimal separation conditions, including links back to the original raw chromatogram reports and the electronic record meets the requirements for electronic records worldwide.

Finally, with the Waters Informatics integrated solution, the often tedious and inefficient Method Development workflow can be significantly streamlined, resulting in a reduction in associated costs, time, and effort.

[APPLICATION NOTE]

References

- 1. V G Dongre, et al. J Pharm Biomed Anal. 46: 236-242, 2008.
- Electronic Laboratory Notebooks-A Foundation for Scientific Knowledge Management, Edition II. Atrium Research & Consulting: 2005.
- 3. European Pharmacopoeia. Council of Europe, Strasburg Cedex. 4th Edition, 2002.
- 4. United States Pharmacopoeia. USP30-NF25, Page 3179, 2002. Pharmacopoeia Forum: Vol. 32(1) 141.
- 5. Ping Du, et al J Assoc Lab Autom. 12(3): 157-165, 2007.

UKAS

001



Waters and NuGenesis are registered trademarks of Waters Corporation. Empower, SDMS Vision Publisher, and The Science of What's Possible are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2008 Waters Corporation. Printed in the U.S.A. April 2008 720002489EN AG-AP



Waters Corporation 34 Maple Street Milford, MA 01757 U.S.A. T: 1 508 478 2000 F: 1 508 872 1990 www.waters.com