

WHISKY: TASTE, SMELL, AND MASS SPECTROMETRY

Waters

THE SCIENCE OF WHAT'S POSSIBLE™



Authors Emma Marsden-Edwards¹, Amy Bartlett¹, Peter Hancock¹, John Shockcor²
Affiliations, ¹Waters Corp, Manchester, UK ²Waters Corp, Milford, USA

INTRODUCTION

Whisky is Scotland's biggest export commodity

There are several different whisky regions each has characteristic flavours...

Whiskys from Islay are renowned for medicinal peatiness and maritime flavours

Speyside whiskys are renowned for elegance, exhibiting flowery, heathery-honey notes and a restrained, fragrant peatiness

Using a chemometric approach based upon principal component analysis (PCA) and the MarkerLynx application manager to mine exact mass UPLC- and GC-TOF/MS data we show that it is possible to distinguish between different whiskys.

This methodology could be for food quality, authenticity and traceability studies.

METHODS

GC/ToF MS

- CTC Analytics CombiPal autosampler with SPME
 - Fibre for flavours, volatiles & semi-volatiles
- Agilent 6890N GC
- GCT Premier
 - Electron Impact (EI+) m/z 50 → 500
 - lock mass = m/z 218.9856 Perfluorotributylamine,
 - Dynamic Range Enhancement (DRE) On
 - Analysis time = 30 min



Fig 1a: GCT Premier Fig 1b: UPLC/Q-ToF Premier

UPLC/Q-ToF Premier MS

- Waters Acuity UPLC column
 - Acuity BEH C18 1.7 μ m 2.1 × 100 mm Column
- Q-ToF Premier
 - Electrospray (EI-) m/z 50 → 1000
 - lock mass = m/z 554.2615 Leucine Enk.
 - Dynamic Range Enhancement (DRE) On
 - Analysis time = 10 min




Fig 2: Small molecule profiling work flow

Figure 2 shows a typical work flow used for small molecule profiling. This work flow was followed during this study.

CONCLUSIONS

- Using GC/MS and UPLC/MS different whiskys could easily be distinguished.**
- The whiskys used in this study do not separate on geographical location alone.**
- Whisky flavour is determined by many factors e.g. type of barrel it is matured in, type of still used etc**
- The application of UPLC/MS and GC/MS-chemometric analysis can be applied for:**
 - Monitoring of product quality
 - Routine monitoring of raw material quality
 - Detection of product adulteration, contamination & substitution
 - Replacement or augmentation of sensory evaluation with analytical instrument systems

RESULTS AND DISCUSSION




Fig 3: (a) Typical chromatograms obtained using GC/TOF/MS
(b) Typical chromatograms obtained using UPLC/TOF/MS




Fig 4: The MarkerLynx Browser containing, sample list, marker table, trend line, scores plot and loadings plot.

STATISTICAL ANALYSIS




Fig 5: The score plot obtained on analysis of the UPLC/TOF MS data




Fig 6: The loadings plot obtained on analysis of the UPLC/TOF MS data




Fig 7: The score plot obtained on analysis of the UPLC/TOF MS data




Fig 8: Screen shot showing the UPLC/TOF MS results table and the option to database search

- Within the results table (fig. 8) the elemental composition of the EMRTs can be calculated. Both accurate mass and elemental composition can be used in a database search to enable identification of the compounds.
- Local databases and on-line searching can be performed. ChemSpider (fig. 9) is a free online chemistry search engine that accesses over 80 different databases.
- Identification of EMRTs identified from the GC-MS analysis can be performed by library searching (fig. 10).




Fig 9: Online database search engine: ChemSpider




Fig 10: Screen shot showing some of the results obtained from searching the NIST library and the GC/TOF MS results table,