

Pretreatment

1. Homogenize 10 g of sample in 100 mL of 0.12 M hydrochloric acid.
2. Take 1 mL aliquot and treat with 400 µL of 50 mM 2-nitrobenzaldehyde in dimethylsulfoxide.
3. Hydrolyze/derivatize the sample for 16 hours at 37 °C.
4. Adjust the sample to pH 7.4 with potassium hydrogen phosphate.
5. Centrifuge sample for 5 minutes at 8000 rpm.

SPE Procedure

Oasis® HLB 3 cc/60 mg

<p>CONDITION/EQUILIBRATE:</p> <p>A. 1 mL methanol</p> <p>B. 1 mL water</p>
<p>LOAD:</p> <p>Approximately 100 mL of sample</p>
<p>WASH:</p> <p>A. 2 mL water</p> <p>B. 2 mL 30% methanol in water</p>
<p>Dry for 20 minutes</p>
<p>ELUTE:</p> <p>3 mL methyl-t-butyl/methanol/formic acid (89:9:2, v/v/v)</p>
<p>Evaporate and reconstitute in 200 µL mobile phase</p>

LC Conditions

Instrument: Waters Alliance® HPLC 2695 System
 Column: XTerra® MS C₁₈, 2.1 x 100 mm, 3.5 µm
 Flow rate: 0.2 mL/min
 Mobile phase: Isocratic 70% 20 mM ammonium formate pH 4, 30% acetonitrile
 Injection volume: 20 µL
 Column temperature: 30 °C


MS Conditions

Instrument: Waters Quattro micro™ API
 Ionization mode: Positive electrospray (ESI+) Multiple reaction monitoring


Analyte	MRM Transition
AOZ	236 → 134
AMOZ	335 → 291
SC	209 → 192
AH	249 → 178

MRM method parameters.

Results



Spiked chicken muscle (1 ng/g) metabolites as 2-nitrobenzaldehyde derivatives.



Ordering Information

Description	Part Number
Oasis HLB, 3 cc/60 mg, 30 µm, 100/box	WAT094226
XTerra MS C ₁₈ , 2.1 x 100 mm, 3.5 µm	186000404
Qsert™ Vials, LCGC Certified Combination Packs	186001126C