
Pretreatment

- 1. Weigh 50 g ground sample with 5 g sodium chloride and place in blender jar.
- 2. Add 100 mL 80:20 methanol: water (v/v) to jar.
- 3. Blend at high speed for 1 minute.
- 4. Filter extract with fluted filter paper. Collect filtrate in a clean vessel.
- 5. Pipet or pour 65 mL filtered extract into a clean vessel.
- 6. Dilute extract with 60 mL of phosphate buffer saline. Mix well.
- 7. Filter extract through glass microfiber filter into a clean vessel.

SPE Procedure


AflaTest[®] affinity column



LC Conditions

Instrument:	Waters Alliance® HPLC System		
Column:	XBridge [™] C ₁₈ , 4.6 x 250 mm, 5 µm		
Flow rate:	1 mL/min		
Mobile phase:	acetonitrile/water/methanol (17:54:29, v/v/v)		
Injection:	100 µL		
Detector:	2475 Multi Wavelength Fluorescence		
Detection:	excitation wavelength: 333 nm emmission wavelength: 460 nm		





Chromatogram of aflatoxins in red pepper extract. Recovery: 76% at 20 ppb (7B1:1B2:3G1:1G2 aflatoxin mix).

Action Levels for Aflatoxins

United States (FDA) action levels (B1, B2, G2, G2, M1)

Food Stuff	Level	Regulation
All products - except milk - designated for humans	20 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
Milk	0.5 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
Corn for immature animals and dairy cattle	20 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
Corn for breeding beef cattle, swine and mature poultry	100 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
Corn for finishing swine	200 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
Corn for finishing beef cattle	300 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
Cottonseed meal (as feed ingredient)	300 ng/g	Policy Guides 7120.26, 7106.10, 7126.33
All feedstuff other than corn	200 ng/g	Policy Guides 7120.26, 7106.10, 7126.33

Aflatoxin regulatory action limits.

Ordering Information

Description	Part Number
VICAM aflaOchra HPLC [™] Columns, 25/box	G1017
VICAM Glass Cuvette	34000
XBridge C ₁₈ , 4.6 x 250 mm, 5 µm	186003117
LCMS Certified Combination Packs	600000751CV

©2008 Waters Corporation. Waters, Alliance, AflaTest, aflaOchra HPLC, and XBridge are trademarks of Waters Corporation. All other trademarks are acknowledged.