

THE DEVELOPMENT OF AN MRM ASSAY FOR QUANTITATION OF LOW ABUNDANCE PROTEINS

A.Davies¹, T.Mckenna¹, C.Hughes¹, J.P.C.Vissers¹, J. Connolly¹, S.Geromanos², C.Donneau², J.Langridge¹
¹. Waters Corporation (MS Technologies Centre), Manchester, United Kingdom. ² Waters Corporation, Milford, MA, United States.

INTRODUCTION

The development of a biomarker assay typically involves the detection, identification and quantitation of multiple proteotypic peptides. Given that, in any one biomarker study, there may be numerous samples containing vast numbers of proteins, a holistic approach to identification and subsequent quantitation is advantageous.

Cytochromes P450 (CYP450) are a large group of heam containing proteins that reside within the cell's endoplasmic reticulum. CYP450s are extensively studied due to their major role in oxidative metabolism and their potential to be induced or inhibited by a wide range of drugs, foods and disease states. The ability to quantify and monitor levels of specific CYP450 enzymes is a useful therapeutic tool. However, high sequence homology between CYP isoforms can make this a challenging task.

This study describes the complete workflow for development of a biomarker assay, from protein identification through to peptide quantification. Here the workflow is applied to identification and quantitation of key CYP450 enzymes in control and chemically induced rat microsomes (liver cells), the latter of which should display up or down regulation of the CYP450 enzymes (see figure 1).

Using the IDENTITY^E high definition proteomics system, CYP450 proteins, in a control rat microsome tryptic digest, were confidently identified from their tryptic peptides. LC-MS^E data from this analysis provided precursor and fragment ion information which was then systematically filtered through a series of criteria to optimise parameters for quantitation. Multiple reaction monitoring (MRM) methods were automatically built from the sorted IDENTITY^E data. Finally, the MRM methods were implemented on a triple quadrupole mass spectrometer and used for quantitation of the peptides.

Chemical Inducer	Main Induced CYP450
CONTROL	
AROCLOR_1254	1A2
NAPTHOFLAVONE	1A2
DEXAMETHAHSONE	3A1
3-METHYLCHLOANTHERE	2A1
PHENOBARBITAL	2B2 and 2C6 and 3A2

Figure 1. The Key CYP450s induced in the Chemically Induced Rat Liver Microsomes used in this study.

PEPTIDE QUANTITATION METHOD DEVELOPMENT WORKFLOW

