

DATA PROCESSING IMPROVEMENTS IN MS PEAK DETECTION FOR TRACE QUANTITATION FROM ACCURATE MASS LC/MS PEPTIDE MAPS

BETH L. GILLECE-CASTRO¹; Marc V. Gorenstein¹; Daniel Golick¹; Keith Richardson²; Barry Dyson²; Scott Berger¹; Jeff Mazzeo¹; Thomas E. Wheat¹; Diane Diehl¹
Waters Corporation; ¹Milford MA, USA and ²Manchester UK

OVERVIEW

- Software tools were developed for application to protein characterization.
- Two novel algorithms were developed to detect and deconvolute peptide signals in LC/MS maps.
- A matching algorithm relates the observed peptides to the protein structure, including modifications.

METHODS

Peak Detection
The first data-analysis algorithm detects the ions obtained in an LC/MS separation. A convolution-based technique measures three key properties of each ion: retention time, mass-to-charge ratio, and intensity. The second algorithm simplifies spectra by selecting only those ions whose retention times fall within restricted ranges. As an example, consider a peptide that elutes at retention time *tr*. All its ions must also elute at *tr*. Variations from *tr* are due only to measurement error. By selecting ions that have the same retention time (to within measurement error), the algorithm simplifies spectra. Such simplified spectra can more clearly reveal the unique, multi-ion signature of peptides.

RESULTS

Figure 2 — Analysis of LC/MS Peptide Map




Figure 2A: Contour Map showing the time-m/z-intensity data matrix for a tryptic digest of bovine hemoglobin between 27 and 33min. Interpretation of this data set and its relation to




Figure 2B: Processed LC/MS Peptide Map from Figure 2A. Each peptide is represented as a “stick” located at the apex retention time. The intensity of the stick is the sum of all of the m/z ions for isotopes and charge states for each component peptide.




Figure 2C: Coverage Map from Figure 2B. Each peptide is matched to the known sequence using exact mass measurement and a hierarchical search.

Figure 3 — Comparing LC/MS Peptide Maps

Two samples of a tryptic digest of hemoglobin were compared. One sample was spiked with nine synthetic peptides. Spike levels were at 2% for the control and 0.2% for the analyte relative to the hemoglobin.




Figure 3A and 3B: Processed LC/MS Peptide Maps of Spiked Digests. The two digests have been spiked with different amounts of peptide standard to create a simulation of two batches of protein with differences in structure. The stacked and mirror image plots are time-aligned to provide alternative views for manual inspection.




Figure 3C: Comparative Coverage Maps LC/MS Peptide Maps of Spiked Digests. The visual comparisons in Figures 3A and 3B are more easily interpreted when each map is subjected to the sequence matching algorithm. Color coding is used to identify peptides that are common to both digests as well as those that are only observed in the control digest or the analyte digest.

Figure 4 — Fragment and Dimer Identification

Source conditions can be altered to enhance fragmentation of peptides. In this example a low level of fragmentation is present to confirm the sequence as shown. Electrospray source conditions may also produce dimer ions when concentrations of peptides are high.




Figure 4: Two fragments of T8, y5 and y6, confirm the identification suggested by the accurate mass measurement. Low to moderate levels of dimer formation in the ion source are observed for peptides T8 and T9. Both fragments and dimers are automatically interpreted by the software.

Figure 5 — Trace Contaminant Detection




Figure 5A: Comparative LC/MS Peptide Maps of Spiked Digests. The specific differences between the digests as shown in the coverage maps can be related to the chromatogram by extracting a difference plot. Spiked peptides are highlighted in purple in the difference plot. The automatically labeled peaks are higher in the control than the analyte as expected since the control was spiked at 2% and the analyte at 0.2%.




Figure 5B: Quantitative Comparisons of LC/MS Peptide Maps of Spiked Digests. The differences between the digests as shown in Figure 5A are clearly identified by the exact mass measurement. The relative intensities for the two different spike levels are in reasonable agreement with the expected values.

Figure 6 — Modification and Disulfide Identification

Forced degradation, such as oxidation, produce highly modified peptides. Multiple oxidations combined with deamidation were observed in this tryptic digest of rabbit phosphorylase b.

Disulfide bonds are identified from non-reduced digest samples and can be compared to a reduced sample.




Figure 6A: The oxidized sample, “analyte”, shows high intensity peptides with multiple oxidations on both Met and Trp residues. Lines which are green show that the control sample was partially oxidized. Lines which are gold show that some new oxidation states are observed after treatment with peroxide.




Figure 6B: Disulfide-bonded peptides are shown in this example of a Lys-C digest. Lines in the table for the control non-reduced sample are blue because the analyte reduced sample had no matching peaks.

CONCLUSIONS

- Retention time alignment and intensity normalization allow accurate comparisons between runs and batches.
- The protein coverage of peptide maps can be readily assessed and compared using software tools.
- Less than 0.5% of trace contaminants and modifications can be detected, and the amounts compared between runs and batches.
- MS source conditions can be set to take advantage of in-source fragmentation for sequence confirmation.
- Disulfide bond linkages can be confirmed.
- Software tools provide efficient matching of chromatographic peaks with known structural features of the protein as well as modifications to that structure.

METHODS

Samples
MassPREP™ Phosphorylase b Digestion Standard
MassPREP™ Hemoglobin Digestion Standard
MassPREP™ Peptide Standards

Add 250 µL water to each vial. Vortex. Final concentration of digest is 4pmol/µL
Oxidized digest
Add 250µL 0.01% hydrogen peroxide to 1 vial Phosphorylase b MassPREP™ Digestion Standard. Vortex. Incubate at room temperature for 2 hours.

Instruments
Waters® ACQUITY UPLC™System including ACQUITY UPLC™BEH300 C₁₈ 1.7µm 2.1x100mm Column
Waters Micromass LCT Premier Mass Spectrometer or Waters Micromass Q-tof Premier Mass Spectrometer

Software
Masslynx 4.1
Biopharmalynx 1.1, Alpha Build 6

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

