

COMPLEMENTARY MASS SPECTROMETRY OR JUST QUESTIONABLE IDENTIFICATIONS?

Scott J. Geromanos
Waters Corporation, Milford, MA 01757

OVERVIEW

- Is it reasonable to assume that the number of identifiable peptides to a protein is directly proportional to its molecular weight and concentration?
- Should replicate analyzes of the same sample or similar samples result in the same proteins being identified?
- Should the same sample or similar sample analyzed in a different laboratory identify the same proteins?
- Is the concentration of certain proteins consistent across different samples?
- Should the fragmentation pattern of the same peptide identified from different samples be internally consistent?
- Or is absence of replication proof that everything that does not replicate is by definition Complementary?

INTRODUCTION



It is generally accepted that the reason for poor reproducibility both within and between platforms on similar if not identical samples is that each analysis is complementary. Basically the belief is that data-directed/dependant analyzes are by nature serendipitous and as such it should not be considered uncommon that different proteins are identified from replicate injections on the same instrument or between instruments or across platforms. However, the protein complement, concentration and ionization efficiency of the enzymatically produced peptides of these samples are for all intents and purposes not only constant but similar in like samples. Careful inspection of the DDA switch lists reveals that ~67% of the ions switched on are the same across replicate analyzes, instruments and platforms so why the lack of reproducibility?

METHODS

All samples were enzymatically digested with trypsin utilizing the protocols described by Dorschel et.al ASMS 2004.



The resulting polypeptide pools were analyzed in triplicate on a QToF Premier and/or Synapt HDMS operating in a data-independent mode (MS^E) of acquisition. Both mass spectrometers were interfaced with a NanoAcuity UPLC system. Column diameters, flow rates and gradient conditions varied with each sample.

For some of the samples Data Directed/Dependent Analyses were acquired in triplicate on both the QToF Premier and Synapt HDMS as well as an LTQ, LTQ Orbitrap, QStar Elite and Agilent 6510.



Though the duty cycle for the trapping instruments was much higher than the orthogonal-time-of-flight instruments, ~15K switches per hour of separation compared to 3–4K (instrument dependant). In all cases ~two-thirds of the switched on ions were the same in all replicate injections.

Peptide Ion Intensity Distribution



The above bar chart illustrates the percent of the total ion intensity detected in a given ion intensity bin for a properly loaded column. As is clearly illustrated ~80% of the detected ions are 2 orders of magnitude lower in intensity than the most intense ions with ~68% being 2.5 orders lower.

The above scatter plot illustrates the femtomolar concentration of common (blue) and uniquely identified (red) proteins between the Synapt MS data and the DDA Orbitrap data. The data is consistent with the peptide results. DDA analyzes have problems identifying the lower intensity peptides efficiently as such the do not identify the lower abundant proteins very well either.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS

	femtmoles of Protein on Column - Peptides Matched											
Accession	Protein MW	Amniotic Fluid	Breast Cancer Depleted	Glioma	Heart Cells	HL60 Cells	NK Cells	Non Depleted	Platelets	Stem Cells	Urine	White BloodCells
1433B_HUMAN	47137	9-8	12-7	7-4	38-5	32-10	19-10	22-13	28-14			
1433E_HUMAN	55209	12-6	21-5	9-5	72-9	152-29	25-9	31-12	25-8			
1433F_HUMAN	22110	14-9	12-8	7-3	34-3	11-10	21-11	19-7				
1433G_HUMAN	28776	10-6	15-7	9-5	42-9	18-6	22-12	26-7				
1433T_HUMAN	90178	13-5	16-4	7-3	39-6	21-11	21-7	26-8	19-8			
1433Z_HUMAN	41605	9-13	18-7	7-5	55-10	36-7	60-20	38-10	110-21			
A26CA_HUMAN	160797	22-9	42-13	23-3	21-10	215-14	249-9	386-12	288-17	1259-23		
A26CB_HUMAN	41774	21-15	40-13	223-5	20-9	206-7	239-7	370-18	276-11	1206-17		
A2MG_HUMAN	53520	209-100	215-112	1-4		128-36	4-18					
ACTA_HUMAN	121370	31-9	41-11	316-6	30-17	209-11	247-5	498-11	351-15	2070-18		
ACTB_HUMAN	28087	39-20	57-13	414-23	22-11	300-23	444-22	823-43	536-38	10-7	2747-48	
ACTBL_HUMAN	28371	19-16	35-6	147-3	16-4	141-10	318-4	298-10	269-11	7-12	1753-13	
ACTC_HUMAN	27764	32-12	48-10	330-6	31-15	218-10	257-3	520-10	366-17	2160-19		
ACTG_HUMAN	29173	41-20	60-13	431-18	23-13	312-21	462-25	857-36	558-34	11-7	2861-51	
ACTH_HUMAN	27745	32-9	42-8	330-8	31-13	218-9	257-4	520-10	366-13	2160-16		
ACTK_HUMAN	82959	45-4	189-2	158-4	330-7	314-5	228-7	980-8				
ACTN1_HUMAN	36557	8-25	17-6	3-7	16-19	55-85	15-23	87-73				
ACTN2_HUMAN	32818		13-4	3-19	9-16	39-21	8-22	33-24				
ACTN3_HUMAN	28390		13-5	2-7	6-14	28-26	7-27	32-25				
ACTN4_HUMAN	50432		16-19	2-8	19-37	48-31	11-22	69-37				
ACTS_HUMAN	121363	31-11	46-7	316-7	30-13	209-8	247-4	498-11	351-14	2070-19		
ALBU_HUMAN	42003	376-53	350-49	141-13	36-9	995-75	625-62	601-60	69-7			
ALDOA_HUMAN	67244	16-22	17-4	32-24	131-23	123-34	17-11	18-15	109-23			
ANXA2_HUMAN	28972	9-12	24-16	47-17	5-15	31-20		33-19	17-10			
ANXAS_HUMAN	22540	8-7	6-4	11-8	22-20		30-15	17-9				
AP0A1_HUMAN	22711	203-43	26-26	264-41	11-10	491-30	19-15					
ATPA_HUMAN	14923	8-18	11-5	27-21	19-14	6-14	21-8	8-15				
ATPB_HUMAN	44843	7-19	15-8	32-24	15-15	7-16	21-11	40-28				
CALM_HUMAN	49907		12-3	79-8	16-9	15-6	27-7	38-9				
CO3_HUMAN	50135	232-147	8-35	112-92	328-162	112-88	6-25					
COF1_HUMAN	18736	11-7	11-4	4-4	45-3	59-14	55-19	28-5	63-14			
COF2_HUMAN	32850		4-2	21-5	9-5	19-6	24-5	28-9				
EF1A1_HUMAN	49895	17-14	25-2	5-4	103-10	60-17		188-21				
ENO1_HUMAN	84528	13-10	25-4	154-8	212-8	14-10	40-15	40-20				
ENO4_HUMAN	103294	15-23	36-12	3-3	174-22	207-37	21-22	60-30	68-36			
ENO5_HUMAN	226401	8-8	5-6	64-6	43-8	12-6	29-10	21-8				
ENOG_HUMAN	36507	9-12	7-3	65-4	44-13	10-7	31-16	17-20				
ENPL_HUMAN	46855	6-11	17-7	2-5	53-12	31-19	6-25	22-21	32-21			
FIBA_HUMAN	49916	28-24	6-10	123-56	2-4		32-36	49-21				
G3P_HUMAN	41784	19-10	53-12	31-14	203-11	463-19	51-21	80-32	198-27			
GELS_HUMAN	92877	15-16	4-16	26-36	3-4	6-4	26-37	8-12	14-7	53-47		
GRP78_HUMAN	18371	7-14	27-14	5-3	88-20	89-20	15-17	42-29	31-22			
GSP1_HUMAN	17881		8-6	3-4	15-7	36-13	6-11	24-9	41-12			
HEMO_HUMAN	123668	193-26	13-10	149-29	5-9	78-8	20-5					
HPT_HUMAN	50151	410-41	17-18	7-12		20-7	2-10	13-13				
H517L_HUMAN	51769	10-12	26-3	6-3	92-8	119-10	18-11	50-8	12-24			
H590A_HUMAN	38472	9-14	18-4	4-7	146-13	148-32	9-18	78-20	27-18			
H590B_HUMAN	41816	12-15	22-9	3-8	158-18	161-43	11-14	92-28	27-17			
H5P71_HUMAN	50140	8-19	17-2	6-8	85-9	74-6	13-18	36-22	13-19			
H5P72_HUMAN	28961	9-7	24-4	6-6	111-15	110-15	15-9	53-18	26-23			
H5P7												