







UTILIZING EXPERIMENTALLY-GENERATED PROTEIN ION MAPS FROM DATA-INDEPENDENT LC-MS ACQUISITIONS FOR IDENTIFYING LOW ABUNDANT PROTEINS IN COMPLEX MIXTURES

Roy Martin¹; J. Will Thompson²; Arthur Moseley³; Scott Geromanos¹
¹Waters Corporation, Beverly, MA; ²Duke University, Igsp, Morrisville, NC; ³Duke University Medical Center, Durham, NC



The above Venn diagrams illustrate the number of peptides and proteins identified in the depleted (green), non-depleted (yellow) and extracted datasets (rose). Querying the Ion Maps from the lower abundant proteins in the depleted serum analysis resulted in a 64% and 57% increase in the number of proteins and peptides identified in the non-depleted analysis.



The scatter plot above shows the comparison of the matched ions from proteins in the non-depleted serum to the extracted ions in the depleted serum sample. The extracted proteins did not pass identification criteria during discovery.

Sample	Ion Type	mhp	RT	Int
Depleted	Extracted	1161.538	32.97	42305
Depleted	Extracted	-235.116	32.99	3121
Depleted	Extracted	-349.19	32.96	1839
Depleted	Extracted	-583.287	32.95	1288
Depleted	Extracted	-1161.53	32.93	1209
Depleted	Extracted	-769.381	32.96	936
Depleted	Matched	-1014.49	32.92	514
Depleted	Matched	-218.148	32.96	476
Depleted	Matched	-420.232	32.97	257
Non-Depleted	Matched	1161.547	33.46	27498
Non-Depleted	Matched	-235.119	33.62	2434
Non-Depleted	Matched	-349.194	33.59	1269
Non-Depleted	Matched	-583.294	33.58	863
Non-Depleted	Matched	-1161.55	33.56	870
Non-Depleted	Matched	-769.39	33.59	730

The above Table illustrates the extracted ions matched to the peptide sequence FSGTWYAMAK from RETBP.HUMAN. The ratio of matched precursor and product ion intensities confirm a correct match.

DISCUSSION

- Identity^E generates high fidelity protein identification based on a reproducible pattern of accurate mass and retention time correlated peptide fragment ions.
- We have shown that these patterns are also reproducible by intensity; the relative intensity of peptides are consistent, and the relative intensity of their fragments are consistent.
- These patterns, consisting of a RT, precursor, and 5 fragment ions, can be used as a "spectral library" to query other data sets.

CONCLUSION

- Once a protein is securely identified extracting the best ionizing peptides and preferred product ions provide a means for identifying that protein at concentration levels below that necessary for database search algorithms.
- Extracted protein identifications can be validated by comparing the ion intensity ratios of "matched" precursor and product ions.
- Using fragmentation and retention-time modeling preferred precursor and product ions generated from an *in-silico* digested can be queried against a spectra database.
- The VerifyE tool can also be used to generate methods for subsequent analyses by MRM techniques.