Sean M. McCarthy, Vera Ivleva, Martin Gilar, Ying-Qing Yu, John C. Gebler Waters Corporation, 34 Maple Street, Milford, MA 01757

## **OVERVIEW**

- There is an increased desire for MS compatible mobile phases for single stranded and duplex oligonucleotide separations.
- A systematic comparison of the resolving power and MS compatibility of mobile phase modifiers for oligonucleotide separations is presented.
- Choice of mobile phase modifier for ion-paring (IP) reversed phase liquid chromatography requires consideration of many factors including IP ability, ease of removal following separation, cost of purchase and disposal, and toxicity.
- We studied various mobile phase compositions for their ability to separate homo and heteromeric oligonucleotides.
- Successful candidates should yield predictable elution order based on oligo length and adequate compatibility with ESI-MS. detection.

# HOMOMERIC OLIGONUCLEOTIDE SEPARATIONS **UV 260nm DETECTION**



Figure 1: Separations were accomplished with a Waters ACQUITY UPLC® System using a Waters Oligonucleotide Separations Technology column (ACQUITY UPLC<sup>®</sup> OST C18, 1.7µm, 2.1x50) maintained at 60 °C. On column loading was 20 pmol/oligo of Waters MassPREP<sup>™</sup> OST standard.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT WWW.WATERS.COM/POSTERS



## SINGLE QUADRUPOLE MS DETECTION



Figure 2: (A) Peak capacities for mobile phases at 100 mM for data shown in Figure 1. (B) Signal to Noise determined from three separate analysis of 25-mer oligo dT via single quadrupole MS.

• Elution conditions (ACN gradient) were adjusted so 15-mer eluted at  $\approx$  5 minutes and 35-mer eluted at  $\approx$  10 min. Separations were run in triplicate with exceptional reproducibility.

• Generally, the more hydrophobic IP agents (longer alkyl chain) required higher ACN content for oligo elution.

 Mobile phases were prepared by the addition of equimolar ratios of acetic acid and appropriate base. The pH was adjusted to  $\approx$  7.0 by the addition of either acetic acid or amine as needed.





• Figure 3: (A) Separation of a digested mixture of 24-mer, 19-mer and 14-mers (5'- CCC CTT GGT TAA CCA AGG TTC CAA-3') Gradients and conditions are the same as those in Figure 1.

## • HAA and TPAA give predictable retention of heteromeric oligonucleotides, but HAA provides better resolution.

• Increased IP efficiency makes the retention pattern of heteromeric oligos more regular.

# **RNAI CHARACTERIZATION WITH LC-MS**

- Use LC-MS compatible mobile phases, such as 15mM TEA, 400 mM **HFIP** ion-pairing buffer.
- Sum mass spectra under the peaks of interest, deconvolute mass failed sequences based on  $\Delta$  mass (Figure 4).



Figure 4: UPLC-MS analysis of crude 21nt RNAi, UUC UGU AAU CUC UUG UCU ATT (5'-3'), ACQUITY UPLC® OST C18, 1.7µm, 2.1x50 mm column, 60 °C, UV 260 nm. m.p. A: 15mM TEA, 400 mM HFIP, pH 7.9, B: 50% m.p. A, 50% MeOH. Gradient 20-40 % B in 10 minutes, 0.2 mL/min.

# **RNA DUPLEX ANALYSIS**

- LC analysis of siRNA can be difficult due to on-column melting and single stranded impurities.
- OST column technology coupled with HPLC and UPLC, using duplex compatible mobile phases, solves these problems by giving predictable retention and non-denaturing separation conditions.
- Use of IP RP HPLC and UPLC analysis of duplexes can be interfaced with MS to provide mass data not available with other methods such as PAGE and AX-HPLC



## Figure 5. Separation of mis-matched duplexes via UPLC with HAA mobile phase. ACQUITY UPLC® OST C18, 1.7µm, 2.1x50 mm column, 60 °C, UV 260 nm.

# Waters THE SCIENCE OF WHAT'S POSSIBLE.™

with MaxEnt1, calculate  $\Delta$  mass between adjacent peaks, assign

incomplete resolution of full length from mis-matched duplexes and



Figure 6: Raw and MaxEnt 1 deconvoluted spectra collected with single quadrupole detector from truncated siRNA separation shown in Figure 5. Data shown is for full length upper strand hybridized with N-1 lower strand.

## **RNA DUPLEX PURIFICATION**

- Purification of siRNA is necessary for therapeutics to ensure specificity and unwanted gene silencing
- Purification accomplished by on-column annealing of siRNA • Inject first strand under initial gradient conditions immediately followed by injection of second complementary strand
- After injection of second strand, gradient started to elute products
- Purification is scalable and siRNA is formed guantitatively
- siRNA is easily collected by proper hearth cutting of peak
- Analysis of collected peak confirms that collected fraction contains desired full length duplex composed of equimolar ratio of each complementary strand



Figure 7. Semi-preparative purification of siRNA via on column annealing of complementary single straned RNA and verification of duplex purity at 20° C (B) and 60° C (C).

## CONCLUSIONS

- HAA and TPAA mobile phases provide exceptional resolution and peak capacities for homo and hetero-oligonucleotide analysis
- IP RP LC allows for MS analysis and sequencing of RNA
- HAA is siRNA and MS compatible and resolves partial duplexes
- IP RP LC combined with non-denaturing column conditions allows for scalable on column annealing and purification of siRNA

