

RAPID, SPECIFIC ANALYSIS OF MELAMINE CONTAMINATION IN INFANT FORMULA AND LIQUID MILK BY UPLC/MS/MS

Jeremy Shia¹, Claude Mallet¹, Michael Young¹, Jianzhong Li², Ying Meng², and Cai Qi²

¹Waters Corporation, Milford, MA U.S.A.; ²Waters Technologies Ltd., Shanghai, China

INTRODUCTION

An increased incidence of kidney stones and renal failure in infants has been reported in China, believed to be associated with the ingestion of infant formula contaminated with melamine. It appears that nitrogen rich melamine was added to raw milk to boost its apparent protein content, which is assessed through determination of the nitrogen content by the Kjeldahl method.

Melamine has many industrial uses that include the production of laminates, adhesives, and melamine resins. Some of these products may become food contact materials, yielding detectable residues in food. Additionally, there are reports that melamine is a metabolite of the pesticide cyromazine. Taking these widespread sources into account, an oral uptake estimate for melamine is 0.007 mg/kg body weight/day.¹

A tolerable daily intake (TDI) has been established by the U.S. Food and Drug Administration (FDA) at 0.63 mg/kg body weight², while the TDI quoted by the European Food Safety Authority (EFSA) is in broad agreement at 0.50 mg/kg body weight.³

Maximum permitted concentrations for melamine in food for adults are approximately 2.5 mg/kg (EU, U.S., and Hong Kong). However, Hong Kong set a tolerance at 1 mg/kg in infant foods. The U.S. FDA stated that no tolerance can be set in infant formula because there are too many uncertainties to rule out public health concerns at a specific level⁴ and Taiwanese authorities have stated that melamine should not be detected in any food using the most sensitive instrumentation.

Hence, there is a need for methodology capable of detecting and quantifying melamine in milk products from parts per million (ppm) levels down to very low parts per billion (ppb) levels.

This application note details a rapid method that combines a specific solid phase extraction (SPE) method for melamine in liquid milk and infant formula with UPLC®/MS/MS detection that can enable detection and quantitation of melamine across a wide range of concentrations.

EXPERIMENTAL

CLICK ON PART NUMBERS FOR MORE INFORMATION

A stock melamine standard (TCI America) was prepared by dissolving 10 mg melamine in 10 mL of 2% ammonium hydroxide. Working standards were prepared by further dilution in mobile phase.

A stable isotope labeled internal standard for melamine (¹³C₃¹⁵N₃ melamine) was obtained from Cambridge Isotope Laboratories, Inc. and a working internal standard was prepared by dilution in mobile phase.

Whole-fat liquid milk and infant formula were obtained from a local supermarket for use in recovery experiments.

Extraction protocol for liquid milk and milk formula:

- Weigh 0.5 g of powder infant formula (10 g liquid milk).
- Add 10 mL 0.2M HClO₄.
- Dissolve by ultrasonication for 10 minutes and centrifuge at > 10,000 rpm.
- Condition Oasis® MCX cartridges (3 cc 60 mg, P/N [186000254](#)) with:
 - 3 mL NH₄OH (5% v/v) in methanol.
 - 3 mL conc. HCl (12 M) (1% v/v) in methanol.
- Equilibrate Oasis MCX cartridges with 3 mL methanol followed by 3 mL water.
- load 3 mL of sample supernatant.
- Wash cartridge with 3 mL water followed by 3 mL methanol.
- Elute using 3 mL NH₄OH (5% v/v) in methanol.
- Evaporate to dryness under nitrogen and reconstitute with 1 mL mobile phase.
- Add 10 µL of ¹³C₃¹⁵N₃ melamine working internal standard.

[APPLICATION NOTE]

UPLC conditions

LC system: Waters® ACQUITY UPLC® System
 Column: ACQUITY UPLC BEH HILIC 2.1 x 100 mm, 1.7 µm
 Part Number: [186003461](#)
 Mobile Phase A: 10 mM ammonium acetate
 Mobile Phase B: 10 mM ammonium acetate in 95:5 acetonitrile/H₂O
 Weak needle wash: 600 µL 10 mM ammonia acetate in 95:5 acetonitrile/H₂O
 Strong needle wash: 200 µL 2% ammonium hydroxide
 Needle type: PEEK

Gradient table

Time (min)	Flow rate (mL/min)	%A	%B	Curve
1. Initial	0.6	0.0	100.0	
2. 0.80	0.6	0.0	100.0	6
3. 2.30	0.6	22.0	78.0	6
4. 2.80	0.6	22.0	78.0	6
5. 2.90	0.6	0.0	100.0	6
6. 4.00	0.6	0.0	100.0	6

MS conditions

MS system: Waters ACQUITY TQ Detector
 Software: Waters MassLynx™ Software v.4.1
 Polarity: ESI+
 Capillary voltage: 3.00 kV
 Cone: 40.00 V
 Source temp.: 150 °C
 Desolvation temp.: 400 °C
 Cone gas flow: 50 L/hr
 Desolvation gas flow: 900 L/hr
 Collision gas flow: 0.25 mL/min

Two Multiple Reaction Monitoring (MRM) transitions were monitored to meet relevant criteria for identification and confirmation of melamine. The MRM transitions, dwell times, cone voltages, and collision energies are listed in Table 1.

MRM transitions	Dwell time (sec)	Cone voltage (V)	Collision energy (eV)
127 > 85	0.04	40	17
127 > 68	0.04	40	25

Table 1. ESI positive MRM conditions for melamine.

RESULTS AND DISCUSSION

To evaluate extraction efficiency of the method, the infant formula was fortified at ten concentrations in duplicate across the range 0.02 ppm to 20 ppm (0.01 ppm to 10 ppm for liquid milk). Although internal standard is usually added prior to extraction to correct for extraction efficiency, during this experiment, it was added after extraction, only correcting for efficiency of ionization. During this small study, average absolute recovery was approximately 80% for infant formula and approximately 90% for whole liquid milk, as shown in Figure 1.




Figure 1. Recoveries of milk and infant formula (50 ppb).

[APPLICATION NOTE]

By including stable isotope-labeled internal standard in the extracts, the linearity of the method could effectively be extended to cover the range of concentration stated above for liquid milk, permitting the use of this procedure for detection and quantitation of melamine in foods destined for both adults and infants (Figure 2).

Employing UPLC in combination with MS/MS detection on the TQ Detector permitted concentrations equivalent to 1 ppb or less to be detected and confirmed, by MRM transitions, as shown in Figure 3.




Figure 2. Calibration curve ranging from 1 ppb to 10 ppm.




Figure 3. Melamine standard 1 ppb – two transition ions.




Figure 4. Infant formula fortified with melamine at 20 ppb.

[APPLICATION NOTE]

Using MRMs for melamine and by the inclusion of a stable isotope-labeled internal standard, infant formula and liquid milk could be screened, confirmed, and the concentration accurately quantified in a single analysis, as shown in Figure 5.




Figure 5. Confirmation of fortified melamine in A) standard and B) liquid milk at 50 ppb.

CONCLUSION

There is a need for a rapid, sensitive method for the analysis of melamine in infant formula and other milk products. This work describes a method that employs a specific solid phase extraction of melamine, using Oasis MCX, from these products following acidification. It is worth noting that the method described here used perchloric acid, but preliminary studies indicate that trichloroacetic acid may equally be substituted.

Following extraction, the ACQUITY UPLC System combined with the TQ Detector provides rapid chromatography (retention time < 2 minutes), with tandem MS detection. Application of UPLC offers significant advantages for speed of analysis resulting in high throughput. Combining UPLC with the TQ detector also offers excellent sensitivity. The addition of a stable isotope internal standard affords good linearity across a wide range of concentrations, resulting in confident quantification of melamine in infant milk and formula.

With the desire of food producers worldwide to demonstrate due diligence regarding the safety and quality of milk products, the method described offers significant business advantages in meeting the challenge of timely, uninterrupted supply of product whilst simultaneously and unequivocally ensuring the safety of consumers.

References

1. World Health Organization: Melamine and Cyanuric Acid: Toxicity, Preliminary Risk Assessment and Guidance on Levels in Food. September 25, 2008. http://www.who.int/foodsafety/fs_management/Melamine.pdf
2. U.S. FDA: Interim Melamine and Analogues Safety/Risk Assessment. May 25, 2007. US FDA/CFSAN - Interim Melamine and Analogues Safety/Risk Assessment
3. EFSA: EFSA'S provisional statement on a request from the European Commission related to melamine and structurally related compounds such as cyanuric acid in protein-rich ingredients used for feed and food. June 7, 2007. http://www.efsa.eu.int/cs/BlobServer/Statement/efsa_statement_melamine_en_rev1.pdf?ssbinary=true
4. FDA Issues Interim Safety and Risk Assessment of Melamine and Melamine-related Compounds in Food October 3, 2008. [FDA Issues Interim Safety and Risk Assessment of Melamine and Melamine-related Compounds in Food](#)

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



Waters, UPLC, Oasis, ACQUITY UPLC, and ACQUITY are registered trademarks of Waters Corporation, MassLynx and The Science of What's Possible are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2008 Waters Corporation. Produced in the U.S.A.
October 2008 720002823EN LB-SC-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com